

MICROSCOPES, CAMERAS & MAGNIFIERS

Compound Microscopes

Upright



JC-211

Jenco

- BK-211, Binocular, Brightfield** \$1,156.00 ↓
- Widefield 10X/20 eyepieces
 - Binocular, 30° inclined, Siedentopf design viewing head with single diopter; 55-75 mm interpupillary adjustment
 - DIN, plan achromat 4X, 10X, 40XR, and 100XR objectives on 5-place, inward facing nosepiece
 - Ball-bearing mechanical stage is 20.9 cm(W) x 14 cm(D)
 - Coaxial coarse and fine focus controls; coarse has tension and up-stop adjustments
 - Centerable 1.25 N.A. Abbe condenser with rack and pinion focusing, and iris diaphragm
 - 20 watt halogen Koehler-type illuminator
 - 27.8 cm(D) x 20 cm(W) x 40.9 cm(H)
 - 5-year mechanical & optical/2-year electrical warranty

MT5300H



Meiji

- MT5300H, Trinocular, Brightfield** \$2,977.00
- Widefield 10X/22 eyepieces
 - Trinocular, 30° inclined, Siedentopf design viewing head has single diopter; integrated male C-mount; 53-75 mm interpupillary adjustment
 - Infinity-corrected, U Plan achromat 4X, 10X, 40XR, and 100XR objectives on 5-place, inward facing nosepiece
 - Ball-bearing, ceramic coated mechanical stage is 19.1 cm(W) x 12.6 cm(D)
 - Coaxial coarse and fine focus controls; coarse has tension and up-stop adjustments
 - Centerable 1.25 N.A. Abbe condenser has rack and pinion focusing, and iris diaphragm
 - 30 watt halogen Koehler-type illuminator
 - 39 cm(D) x 21 cm(W) x 41 cm(H)
 - Lifetime mechanical & optical/90-day electrical warranty



168-ASC



169-ASC

National

- 168-ASC, Binocular, Brightfield** \$738.00
- Widefield 10X/18 eyepieces
 - Binocular, 30° inclined, Jentsch design viewing head with dual diopters; 54-76 mm interpupillary adjustment
 - DIN, super high contrast achromat 4X, 10X, 40XR, and 100XR objectives on 4-place, inward facing nosepiece
 - Rackless x-axis mechanical stage is 14 cm(W) x 13.5 cm(D); up-stop adjustment is provided
 - Coaxial coarse and fine focus controls; coarse has tension adjustment
 - 1.25 N.A. Abbe condenser with rack and pinion focusing, iris diaphragm, and swing-out filter holder
 - 3 watt variable LED illuminator, corded
 - 24 cm(D) x 18.5 cm(W) x 40 cm(H)
 - Lifetime mechanical & optical/3-year electrical warranty
- 169-ASC, Trinocular, Brightfield** \$942.00
- Trinocular version of above
 - 0.5X C-mount adapter included

M15B-P



Swift

- M15B-P, Binocular, Brightfield** \$989.00
- Widefield 10X/20 eyepieces
 - Binocular, 30° inclined, Siedentopf design viewing head has single diopter; 54-76 mm interpupillary adjustment
 - Infinity-corrected, plan achromat 4X, 10X, 40XR, and 100XR objectives on 4-place, inward facing nosepiece
 - Ball-bearing, mechanical stage is 14.2 cm(W) x 14 cm(D); up-stop adjustment is provided
 - Coaxial coarse and fine focus controls; coarse has tension adjustment
 - 1.25 N.A. Abbe condenser with rack and pinion focusing, iris diaphragm, and swing-out filter holder
 - 20 watt halogen illuminator
 - Carrying handle
 - 39 cm(D) x 18.5 cm(W) x 37 cm(H)
 - Lifetime mechanical and optical/3-year electrical warranty
- M15T-P, Trinocular, Brightfield** \$1,237.00
- Trinocular version of above

Continued