


LAB, SURGICAL & DENTAL INSTRUMENTS

Legend

- Economy:** stainless steel; forged and finished in Pakistan
Standard: surgical-grade (OR) German steel; forged in Germany and finished in Pakistan and/or Germany
German: surgical-grade (OR) German steel; made entirely (forged and finished) in Germany
GT: George Tiemann & Co, Germany
RS-: Roboz Surgical
-TC: Tungsten Carbide inlays; identified by gold handles

Forceps (Tweezers) & Tongs



General

Forceps

- Fine Point, serrated handles, stainless steel
 Straight, serrated tips
 10 cm (4"), Economy, guide pin. **26-275** \$7.55
 10 cm (4"), Standard, guide pin. **40-720** \$14.40
 10 cm (4"), German, guide pin. **420-802** \$32.80
 ★ 11.5 cm (4½"), Economy, guide pin. **300-130** \$2.60
 11.5 cm (4½"), Economy, No Pin. **20-373** \$2.95
 14 cm (5½"), German. **57.598** \$24.20
 20.5 cm (8"), Economy, guide pin. **300-141** \$7.50
 25.4 cm (10"), Economy, guide pin. **300-143** \$8.50
 Straight, Smooth Tips
 20.5 cm (8"), economy, guide pin. **300-111** \$8.50
 Angled, serrated tips, 15.5 cm (6"), guide pin. **3600** \$4.00
 Half Curved (45°), serrated tips
 10 cm (4"), Economy, guide pin. **26-280** \$7.70
 10 cm (4"), Standard, guide pin. **40-722** \$14.40
 10 cm (4"), German, guide pin. **420-810** \$33.90
 ★ 11.5 cm (4½"), Economy, No Pin. **4190** \$3.30
 ★ 12 cm (4¾"), Economy, guide pin. **300-131** \$2.70
 Full Curved (90°), serrated tips
 10 cm (4"), Economy, guide pin. **24-665S** \$6.10
 Medium-Fine Point, serrated handles, stainless steel
 Curved, smooth tips, no pin, 15.2 cm (6"). **300-011** \$6.90

- Medium Point (Thumb Dressing), serrated tips & handles
 Straight, no guide pin, stainless steel
 9 cm (3½"), Economy. **17-315** \$2.25
 11.5 cm (4½"), Economy. **300-038 (17-035)** \$2.65
 11.5 cm (4½"), Standard. **540-035** \$10.50
 12 cm (4¾"), German. **420-455** \$15.50
 15.5 cm (6"), Economy. **300-041 (17-050)** \$3.15
 15.5 cm (6"), Standard. **540-050** \$12.60
 15.5 cm (6"), German. **420-562** \$19.75
 16 cm (6¼"), Economy, black matte. **TWZ-703.74** \$3.95
 20.5 cm (8"), Economy. **2240** \$4.05
 20.5 cm (8"), Standard. **540-055** \$21.00
 25.4 cm (10"), Economy. **2250** \$4.80
 25.4 cm (10"), Standard. **540-060** \$31.50
 30.5 cm (12"), Economy. **2260** \$5.85
 51 cm (20"), Economy. **19-380** \$17.90

Exclusive!

Extra wide tips ideal for removing specimens from storage containers.

- Medium Point, serrated tips & handles, stainless steel
 Curved, 11.5 cm (4½"), Economy, guide pin. **17-655 (3150)** \$3.30



smooth handles do not have guide pin

- Gooseneck, blunt serrated tips, 15 cm (6"). **300-001** \$4.90



- Student Grade, stainless steel, 11.5 cm (4½"). **18-915** \$1.30
 straight, serrated tips & handles



- Plastic, straight, serrated tips & handles, 13 cm. **97-3901** \$0.45
 Corrosion and heat resistant ABS plastic is nonconductive, and noncontaminating. Will not scratch optics.

Micro



Economy Set: styles 5, 3c, 6, 7, 532, SS, AA, 3, AF, 2, OOB, 2A

Forceps, Micro, smooth tips*, no guide pins

- Economy, stainless steel
 #2A Style, straight, 12 cm. **57.0752A** \$5.90
 #3 Style, straight, 12 cm. **19-361** \$5.20
 ★ #5 Style, straight, 11.5 cm. **19-364** \$5.20
 #7 Style, curved, 11.5 cm. **19-366** \$5.20
 #7B Style, curved, *serrated, nonmag, 11.5 cm. **57.838** \$10.90
 Economy Set, 12-pc, nonmagnetic. **TWZ-700.12** \$33.60
 styles: 5, 3c, 6, 7, 532, SS, AA, 3, AF, 2, OOB, 2A



Style #3 (very sharp)



Style #5 (needle sharp)



Style #7 (sharp)

Dumont Economy, Inox, magnetic, Swiss

- #3 Style, 12 cm (4¾"). TWZ-301.16 \$25.30
- #5 Style, 11 cm (4¼"). 57.520 \$26.80
- #7 Style, 11.5 cm (4½"). TWZ-301.28 \$29.40
- SS Style, nonmag, straight, 14 cm (5½"). TWZ-302.56..... \$22.65

Dumont, Inox, high-precision, magnetic, Swiss

- #3 Style, 12 cm (4¾"). RS-5042 \$35.75
- #5 Style, 11 cm (4¼")
 - Standard Tips (0.10W x 0.06 mm). RS-5045 \$29.50
 - Biologie Tips (0.05W x 0.01 mm). RS-5015 \$35.60
- #55 Style, light shanks, 11 cm (4¼")
 - Biologie Tips (0.05W x 0.01 mm). RS-5010 \$42.20
- #7 Style, 11.5 cm (4½"). TWZ-201.28 \$33.10



Dumont, Dumoxel, standard tips (0.10W x 0.06 mm), Swiss
#5/45 Style, tips angled 45°, 11 cm (4¼"). RS-5058 \$37.35



Rubis Swiss Tweezer®, perforated handle (1 side only), Swiss
#5G Style, straight, 11 cm (4¼"). 57.184 \$31.00
nonmagnetic stainless steel

Specialty



Forceps, stainless steel
Adson Dressing, serrated tips & handles, 12 cm. 17-555 \$6.00



Balance, curved plastic tips, 10 cm (4"). 300-006 \$4.80



Bishop Harmon, fine serrated tips, 9 cm (3½"). 4211 \$18.00



Cupped, epoxy-coated tips, 14.5 cm (5¾"). TWZ-104.84 \$7.25



Cutting, economy
Style 15A, 11.5 cm (4½"). TWZ-801.60 \$9.95



Diamond, grid-serrated tips, 15.5 cm (6"), Italy. TWZ-171.74 \$9.80
serrated handles, medium point



Filter, smooth tips & polypropylene grips, 11.5 cm. 300-025 \$4.35



Foerster, octagonal grips, 0.5 mm serrated tips, 9.5 cm (3¾")
Curved, economy. 4202 \$18.00



Graefe Eye Dressing, serrated tips & handles, 7 cm (2½")
Straight, standard. 4215 \$18.00
Curved, standard. 4216 \$18.00



Insect Pinning, rounded & serrated tips, 11 cm. 17-682 \$4.15



Paddle Tip, grid-serrated pads, economy
22 x 22 mm Pad, 16 cm (6¼"). 33-1045 \$14.50
30 x 35 mm(W) Pad, 15.2 cm (6"). 33-1045B \$21.00
30 x 40 mm(W) Pad, 15.2 cm (6"). 33-1045A \$21.00



Wide Handle (Paddle Handle), economy
serrated and perforated grips
Straight, Serrated Tips
14 cm (5½"), 1.5 mm tips. 15-600 \$24.75
Straight, Smooth Tips
14 cm (5½"), 0.6 mm tips. 15-625 \$24.75
14 cm (5½"), 1.5 mm tips. 15-585 \$24.75
17.5 cm (7"), 3 mm tips. 15-500 \$24.75
Angled, Serrated Tips
14 cm (5½"), 1.5 mm tips. 15-595 \$24.75
Angled, Smooth Tips
14 cm (5½"), 1.5 mm tips. 15-590 \$24.75

Spade Tips (Kuehne Pan)

Angled



Stainless Steel, 11 cm (4 1/4"), economy. **33-1041**..... \$4.15
(3270)

Teflon® Coated, 11 cm (4 1/4"), Bel-Art. **379420000**..... \$20.90
Straight, stainless steel, 10.5 cm, economy. **57.926** \$5.85

Square Point Tips, straight



Stainless Steel, 12 cm (4 3/4"), economy. **1244**..... \$2.25
Teflon® Coated, 11 cm (4 1/4"), Bel-Art. **379430000**..... \$22.90

Locking



300-066

Forceps, stainless steel

Self-Closing (cross-locking), squeeze handles to open tips

Fine Point, smooth tips

Economy, 13 cm (5 1/8"). **300-066**..... \$4.95

Peer-Vigor, 11 cm (4 1/4"), Swiss. **57.324**..... \$42.40

Medium Point, serrated tips, 17 cm (6 5/8"). **TWZ-948.00** \$3.95



TWZ-101.76

Slide-Locking, grid-serrated tips

Fine Point, Economy, 14 cm (5 1/2"). **57.0357A** \$9.00

Fine Point, Dumont, 16 cm (6 1/4"). **TWZ-101.72** \$45.75

Medium Point, Economy, 16 cm (6 1/4"). **203** \$9.00

Medium Point, Dumont, 16 cm (6 1/4"). **TWZ-101.76**..... \$43.60

Snap-Closing, angled, fine points, 15.5 cm (6"). **3610** \$5.80
serrated tips & handles, guide pin

Tissue



17-574

Forceps, stainless steel, economy


Adson-Brown, 9 x 9 side grasping teeth, 12 cm. **17-574**..... \$15.75

Adson Tissue, 1 x 2 teeth, 12 cm (4 3/4"). **17-560** \$6.00

Iris Tissue, 1 x 2 teeth, 10 cm (4"). **26-295** \$7.70



2275



17-145

Russian Tissue, circular serrated tips, 15.5 cm (6"). **2275**..... \$5.25

Tissue, 1 x 2 teeth, 14 cm (5 1/2"). **17-145** \$4.00

Tongs




320-101

Tongs, Crucible, serrated jaws, stainless steel

24.5 cm (9 3/4"). **320-101**..... \$11.60

51 cm (20"). **320-120**..... \$27.75


Hemostats & Clamps



Kelly, Straight

Kelly, Curved

Mosquito, Curved



Polypropylene Hemostat

Hemostats, serrated jaws

Hartman Mosquito, stainless steel

Straight, Economy, 9 cm (3 1/2"). **12-180S** \$6.05

Straight, Standard, 10 cm (4"). **525-025**..... \$25.20

Curved, Economy, 9 cm (3 1/2"). **14-030** \$5.80

Curved, Standard, 10 cm (4"). **20-360**..... \$21.70

Halstead Mosquito, 12.5 cm (5"), stainless steel

Straight, Economy. **14-065** \$6.80

Straight, Standard. **20-501** \$22.45

Curved, Economy. **14-070**..... \$6.80

Curved, Standard. **20-519** \$22.85

Kelly, stainless steel

Straight, Economy, 14.5 cm (5 3/4"). **14-225** \$6.20

Straight, Standard, 14.5 cm (5 3/4"). **525-225**..... \$22.40

Curved, Economy, 14 cm (5 1/2"). **14-230** \$6.20

Curved, Standard, 14.5 cm (5 3/4"). **20-568**..... \$22.40

Rochester-Pean, 20.5 cm (8"), stainless steel

Straight, Economy. **1920** \$6.50

Straight, Standard. **525-435** \$44.10

Curved, Economy. **1930** \$6.50

Curved, Standard. **525-440** \$44.10

Polypropylene, 12.5 cm (5"). **HEM-300.00** \$1.50



Specimen Forceps, 25 cm (9 3/4"), serrated jaws. **300-500** \$17.25