


LAB, SURGICAL & DENTAL INSTRUMENTS

Legend

- Economy:** stainless steel; forged and finished in Pakistan
Standard: surgical-grade (OR) German steel; forged in Germany and finished in Pakistan and/or Germany
German: surgical-grade (OR) German steel; made entirely (forged and finished) in Germany
GT: George Tiemann & Co, Germany
RS-: Roboz Surgical
-TC: Tungsten Carbide inlays; identified by gold handles

Forceps (Tweezers) & Tongs



General

Forceps


- Fine Point, serrated handles, stainless steel
 Straight, serrated tips
 10 cm (4"), Economy, guide pin. **26-275** \$7.55
 10 cm (4"), Standard, guide pin. **40-720** \$14.40
 10 cm (4"), German, guide pin. **420-802** \$32.80
 ★ 11.5 cm (4½"), Economy, guide pin. **300-130** \$2.60
 11.5 cm (4½"), Economy, No Pin. **20-373** \$2.95
 14 cm (5½"), German. **57.598** \$24.20
 20.5 cm (8"), Economy, guide pin. **300-141** \$7.50
 25.4 cm (10"), Economy, guide pin. **300-143** \$8.50
 Straight, Smooth Tips
 20.5 cm (8"), economy, guide pin. **300-111** \$8.50
 Angled, serrated tips, 15.5 cm (6"), guide pin. **3600** \$4.00
 Half Curved (45°), serrated tips
 10 cm (4"), Economy, guide pin. **26-280** \$7.70
 10 cm (4"), Standard, guide pin. **40-722** \$14.40
 10 cm (4"), German, guide pin. **420-810** \$33.90
 ★ 11.5 cm (4½"), Economy, No Pin. **4190** \$3.30
 ★ 12 cm (4¾"), Economy, guide pin. **300-131** \$2.70
 Full Curved (90°), serrated tips
 10 cm (4"), Economy, guide pin. **24-665S** \$6.10
 Medium-Fine Point, serrated handles, stainless steel
 Curved, smooth tips, no pin, 15.2 cm (6"). **300-011** \$6.90

- Medium Point (Thumb Dressing), serrated tips & handles
 Straight, no guide pin, stainless steel
 9 cm (3½"), Economy. **17-315** \$2.25
 11.5 cm (4½"), Economy. **300-038 (17-035)** \$2.65
 11.5 cm (4½"), Standard. **540-035** \$10.50
 12 cm (4¾"), German. **420-455** \$15.50
 15.5 cm (6"), Economy. **300-041 (17-050)** \$3.15
 15.5 cm (6"), Standard. **540-050** \$12.60
 15.5 cm (6"), German. **420-562** \$19.75
 16 cm (6¼"), Economy, black matte. **TWZ-703.74** \$3.95
 20.5 cm (8"), Economy. **2240** \$4.05
 20.5 cm (8"), Standard. **540-055** \$21.00
 25.4 cm (10"), Economy. **2250** \$4.80
 25.4 cm (10"), Standard. **540-060** \$31.50
 30.5 cm (12"), Economy. **2260** \$5.85
 51 cm (20"), Economy. **19-380** \$17.90

Exclusive!

Extra wide tips ideal for removing specimens from storage containers.

- Medium Point, serrated tips & handles, stainless steel
 Curved, 11.5 cm (4½"), Economy, guide pin. **17-655 (3150)** \$3.30




smooth handles do not have guide pin

- Gooseneck, blunt serrated tips, 15 cm (6"). **300-001** \$4.90



- Student Grade, stainless steel, 11.5 cm (4½"). **18-915** \$1.30
 straight, serrated tips & handles



- Plastic, straight, serrated tips & handles, 13 cm. **97-3901** \$0.45
 Corrosion and heat resistant ABS plastic is nonconductive, and noncontaminating. Will not scratch optics.

Micro



Economy Set: styles 5, 3c, 6, 7, 532, SS, AA, 3, AF, 2, OOB, 2A

Forceps, Micro, smooth tips*, no guide pins

- Economy, stainless steel
 #2A Style, straight, 12 cm. **57.0752A** \$5.90
 #3 Style, straight, 12 cm. **19-361** \$5.20
 ★ #5 Style, straight, 11.5 cm. **19-364** \$5.20
 #7 Style, curved, 11.5 cm. **19-366** \$5.20
 #7B Style, curved, *serrated, nonmag, 11.5 cm. **57.838** \$10.90
 Economy Set, 12-pc, nonmagnetic. **TWZ-700.12** \$33.60
 styles: 5, 3c, 6, 7, 532, SS, AA, 3, AF, 2, OOB, 2A



Style #3 (very sharp)



Style #5 (needle sharp)



Style #7 (sharp)

Dumont Economy, Inox, magnetic, Swiss

- #3 Style, 12 cm (4¾"). TWZ-301.16 \$25.30
- #5 Style, 11 cm (4¼"). 57.520 \$26.80
- #7 Style, 11.5 cm (4½"). TWZ-301.28 \$29.40
- SS Style, nonmag, straight, 14 cm (5½"). TWZ-302.56..... \$22.65

Dumont, Inox, high-precision, magnetic, Swiss

- #3 Style, 12 cm (4¾"). RS-5042 \$35.75
- #5 Style, 11 cm (4¼")
 - Standard Tips (0.10W x 0.06 mm). RS-5045 \$29.50
 - Biologie Tips (0.05W x 0.01 mm). RS-5015 \$35.60
- #55 Style, light shanks, 11 cm (4¼")
 - Biologie Tips (0.05W x 0.01 mm). RS-5010 \$42.20
- #7 Style, 11.5 cm (4½"). TWZ-201.28 \$33.10




Dumont, Dumoxel, standard tips (0.10W x 0.06 mm), Swiss
#5/45 Style, tips angled 45°, 11 cm (4¼"). RS-5058 \$37.35



Rubis Swiss Tweezer®, perforated handle (1 side only), Swiss
#5G Style, straight, 11 cm (4¼"). 57.184 \$31.00
nonmagnetic stainless steel

Specialty



Forceps, stainless steel
Adson Dressing, serrated tips & handles, 12 cm. 17-555 \$6.00



Balance, curved plastic tips, 10 cm (4"). 300-006 \$4.80



Bishop Harmon, fine serrated tips, 9 cm (3½"). 4211 \$18.00



Cupped, epoxy-coated tips, 14.5 cm (5¾"). TWZ-104.84 \$7.25



Cutting, economy
Style 15A, 11.5 cm (4½"). TWZ-801.60 \$9.95



Diamond, grid-serrated tips, 15.5 cm (6"), Italy. TWZ-171.74 \$9.80
serrated handles, medium point



Filter, smooth tips & polypropylene grips, 11.5 cm. 300-025 \$4.35



Foerster, octagonal grips, 0.5 mm serrated tips, 9.5 cm (3¾")
Curved, economy. 4202 \$18.00



Graefe Eye Dressing, serrated tips & handles, 7 cm (2½")
Straight, standard. 4215 \$18.00
Curved, standard. 4216 \$18.00



Insect Pinning, rounded & serrated tips, 11 cm. 17-682 \$4.15



Paddle Tip, grid-serrated pads, economy
22 x 22 mm Pad, 16 cm (6¼"). 33-1045 \$14.50
30 x 35 mm(W) Pad, 15.2 cm (6"). 33-1045B \$21.00
30 x 40 mm(W) Pad, 15.2 cm (6"). 33-1045A \$21.00



Preferred by histologists 15-595

Wide Handle (Paddle Handle), economy
serrated and perforated grips

- Straight, Serrated Tips
 - 14 cm (5½"), 1.5 mm tips. 15-600 \$24.75
- Straight, Smooth Tips
 - 14 cm (5½"), 0.6 mm tips. 15-625 \$24.75
 - 14 cm (5½"), 1.5 mm tips. 15-585 \$24.75
 - 17.5 cm (7"), 3 mm tips. 15-500 \$24.75
- Angled, Serrated Tips
 - 14 cm (5½"), 1.5 mm tips. 15-595 \$24.75
- Angled, Smooth Tips
 - 14 cm (5½"), 1.5 mm tips. 15-590 \$24.75

Spade Tips (Kuehne Pan)

Angled



Stainless Steel, 11 cm (4 1/4"), economy. **33-1041**..... \$4.15
(3270)

Teflon® Coated, 11 cm (4 1/4"), Bel-Art. **379420000**..... \$20.90

Straight, stainless steel, 10.5 cm, economy. **57.926** \$5.85

Square Point Tips, straight



Stainless Steel, 12 cm (4 3/4"), economy. **1244**..... \$2.25

Teflon® Coated, 11 cm (4 1/4"), Bel-Art. **379430000**..... \$22.90

Locking



Forceps, stainless steel

Self-Closing (cross-locking), squeeze handles to open tips

Fine Point, smooth tips

Economy, 13 cm (5 1/8"). **300-066**..... \$4.95

Peer-Vigor, 11 cm (4 1/4"), Swiss. **57.324**..... \$42.40

Medium Point, serrated tips, 17 cm (6 5/8"). **TWZ-948.00** \$3.95



Slide-Locking, grid-serrated tips

Fine Point, Economy, 14 cm (5 1/2"). **57.0357A** \$9.00

Fine Point, Dumont, 16 cm (6 1/4"). **TWZ-101.72** \$45.75

Medium Point, Economy, 16 cm (6 1/4"). **203** \$9.00

Medium Point, Dumont, 16 cm (6 1/4"). **TWZ-101.76**..... \$43.60

Snap-Closing, angled, fine points, 15.5 cm (6"). **3610** \$5.80
serrated tips & handles, guide pin

Tissue




Forceps, stainless steel, economy

Adson-Brown, 9 x 9 side grasping teeth, 12 cm. **17-574**..... \$15.75

Adson Tissue, 1 x 2 teeth, 12 cm (4 3/4"). **17-560** \$6.00


Iris Tissue, 1 x 2 teeth, 10 cm (4"). **26-295** \$7.70



Russian Tissue, circular serrated tips, 15.5 cm (6"). **2275**..... \$5.25

Tissue, 1 x 2 teeth, 14 cm (5 1/2"). **17-145** \$4.00

Tongs




Tongs, Crucible, serrated jaws, stainless steel

24.5 cm (9 3/4"). **320-101**..... \$11.60

51 cm (20"). **320-120**..... \$27.75

Hemostats & Clamps



Hemostats, serrated jaws

Hartman Mosquito, stainless steel

Straight, Economy, 9 cm (3 1/2"). **12-180S** \$6.05

Straight, Standard, 10 cm (4"). **525-025**..... \$25.20

Curved, Economy, 9 cm (3 1/2"). **14-030** \$5.80

Curved, Standard, 10 cm (4"). **20-360**..... \$21.70

Halstead Mosquito, 12.5 cm (5"), stainless steel

Straight, Economy. **14-065** \$6.80

Straight, Standard. **20-501** \$22.45

Curved, Economy. **14-070**..... \$6.80

Curved, Standard. **20-519** \$22.85

Kelly, stainless steel

Straight, Economy, 14.5 cm (5 3/4"). **14-225** \$6.20

Straight, Standard, 14.5 cm (5 3/4"). **525-225**..... \$22.40

Curved, Economy, 14 cm (5 1/2"). **14-230** \$6.20

Curved, Standard, 14.5 cm (5 3/4"). **20-568**..... \$22.40

Rochester-Pean, 20.5 cm (8"), stainless steel

Straight, Economy. **1920** \$6.50

Straight, Standard. **525-435** \$44.10

Curved, Economy. **1930**..... \$6.50

Curved, Standard. **525-440** \$44.10

Polypropylene, 12.5 cm (5"). **HEM-300.00** \$1.50



Specimen Forceps, 25 cm (9 3/4"), serrated jaws. **300-500** \$17.25

Accessories



Bibulous Paper, 10 x 15 cm, Bk/50 sheets. **BIB-46-50** bk cs/48 bk
 \$3.65 169.00

Forceps

Spade Tips

- Angled (Kuehne)
 Stainless Steel, 11 cm (4¼"), economy. **3270** \$4.15
 Teflon® Coated, 11 cm (4¼"), Bel-Art. **379420000** \$20.90
 Straight, stainless steel, economy. **57.926** \$4.90

Square Point Tips, straight

- Stainless Steel, 12 cm (4¾"), economy. **1244** \$2.25
 Teflon® Coated, 11 cm (4¼"), Bel-Art. **379430000** \$22.90



Minimize contamination during slide & coverslip handling. Will not scratch; corrosion and heat resistant.

Hemostat, Polypropylene, 12.5 cm. **HEM-300.00** \$1.50

Cytology Funnels

Single,
Metal Clips
Required



Double,
All Plastic



White filter cards will accept up to 0.5 mL of sample, while tan filters will accept volumes up to 0.4 mL.

Cytology Funnels, Simport CytoSep™

Shandon (Thermo) Cytospin® Cytocentrifuges

Single. *Metal slide clips required*

White Filter & Cap, Cs/500. **M964-10FW** \$382.00
 (59 91 040, 110 25 48)

Tan Filter & Cap, Cs/500. **M964-10FT** \$382.00

Single. *All plastic - metal slide clips not required*

White Filter & Cap, Cs/480. **M965-10FW** \$503.00
 (A787 100003)

Tan Filter & Cap, Cs/480. **M965-10FT** \$503.00

Double. *All plastic - metal slide clips not required*

White Filter & Cap, Cs/480. **M965-20FW** \$626.00
 (A787 100005)

Immersion Oils & Mounting Media



Cargille immersion oils A, B, and NVH can be mixed with each other for custom viscosities. USA

Immersion Oil

Type A, low viscosity (150 cSt). *Short W.D. optics*

- 7.5 mL in LDPE dispensing bottle. **18101** \$5.35
 20 mL in HDPE dispensing bottle. **18102** \$10.00
 120 mL in amber bottle. **18103** \$22.45
 480 mL in amber bottle. **18104** \$49.60

Type B, high viscosity (1250 cSt). *Long W.D. optics*

- 7.5 mL in LDPE dispensing bottle. **18105** \$5.35
 20 mL in HDPE dispensing bottle. **18106** \$10.00
 120 mL in amber bottle. **18107** \$22.45
 480 mL in amber bottle. **18108** \$49.60


Type NVH, very high viscosity (21,000 cSt)

- 480 mL in amber bottle. **18113** \$57.75

Type FF, fluorescence free (170 cSt)

- 7.5 mL in LDPE dispensing bottle. **18121** \$7.35
 480 mL in amber bottle. **18124** \$81.75

A long nozzle, textured sides, and slim profile make this chemically inert, high-density polyethylene (HDPE) container ideal for dispensing immersion oil.



Bottle, Dispensing, HDPE, 20 mL, each. **18145** \$6.90



Store immersion oil or mounting media in this borosilicate glass bottle with an applicator rod and ground glass cap.

Bottle, Balsam, glass, 45 mL (empty), USA. **616** \$17.30

K-750

Suitable for fluorescence microscopy.



K-700

Takes only 10 minutes to dry.

Mounting Media

Mount-Quick, Kiyota

- 30 mL in Dropper Bottle, aqueous-based. **K-750** \$22.40
 30 mL in Tube, solvent-based, 30 cps. **K-700** \$15.40

Very Fast flow (360 mL/min), 35 µm retention, creped, Grade 617	
9 cm, Pk/50. 6170-0900	\$6.05
12.5 cm, Pk/50. 6170-1250	\$8.15
15 cm, Pk/50. 6170-1500	\$9.50
20.5 cm, Pk/50. 6170-2050	\$14.45
33 cm, Pk/50. 6170-3300	\$30.95

Qualitative, Wet-Strengthened Grades (for vacuum filtration)

Medium flow (65 mL/min), 3 µm retention, Grade 909 (<i>Whatman 91; S&S 0860, 1574, 3003; Reeve-Angel 201</i>)	
32 cm, Pk/100. 9090-3200	\$61.40
Fast flow (245 mL/min), 27 µm retention, creped, Grade 962 (<i>Whatman 114; S&S 0856, 0905, 3002; Reeve-Angel 202</i>)	
11 cm, Pk/100. 9620-1100	\$9.15

Qualitative, Pleated (Fluted)

Slow flow (30 mL/min), 2 µm retention, Grade 542 (<i>Whatman 2V</i>)	
12.5 cm, Pk/100. 5420-1250	\$65.15
18.5 cm, Pk/100. 5420-1850	\$73.75
Medium flow (60 mL/min), 6 µm retention, Grade 513	
32 cm, Pk/100. 5130-3200	\$86.65

Quantitative, Ashless

(Designed for use in analytical and gravimetric analysis)

Slow flow (4 mL/min), 1.5 µm retention, Grade 94 (<i>Whatman 42; S&S 589/3</i>)	
9 cm, Pk/100. 0940-0900	\$26.35
15 cm, Pk/100. 0940-1500	\$45.35
Medium flow (30 mL/min), 2 µm retention, Grade 74 (<i>Whatman 40; S&S 589/2, 589/4</i>)	
9 cm, Pk/100. 0740-0900	\$26.35
15 cm, Pk/100. 0740-1500	\$45.35
Fast flow (200 mL/min), 10 µm retention, Grade 54 (<i>Whatman 41; S&S 589/1</i>)	
9 cm, Pk/100. 0540-0900	\$26.35
12.5 cm, Pk/100. 0540-1250	\$36.05
15 cm, Pk/100. 0540-1500	\$45.35
Fast flow (220 mL/min), 15 µm retention, Grade 55 (<i>Whatman 54, 541; S&S 1505</i>)	
9 cm, Pk/100. 0550-0900	\$29.75
15 cm, Pk/100. 0550-1500	\$47.80

Glass Microfiber


Medium flow (102 mL/min), 1.1 µm retention, Grade 161 (<i>Whatman 934-AH; S&S GF30,GF50; Reeve-Angel 934-AH</i>)	
2.4 cm, Pk/100. 1610-0240	\$21.95
4.7 cm, Pk/100. 1610-0470	\$31.10
5.5 cm, Pk/100. 1610-0550	\$35.90
7 cm, Pk/100. 1610-0700	\$44.10
9 cm, Pk/100. 1610-0900	\$57.20
11 cm, Pk/100. 1610-1100	\$76.15

Blotting/Chromatography/Electrophoresis Papers

0.19 mm thick, 125 mm/30 min absorption, Grade 601 (<i>Whatman 1, 1CHR; S&S 593, 597, 2043A</i>)	
46 x 57 cm, Pk/100. 6018-4657	\$107.00
0.34 mm thick, 140 mm/30 min absorption, Grade 238 (<i>Whatman 3MM, S&S 598,2316</i>)	
7 x 10 cm, Pk/100. 2388-0710	\$25.30
20 x 20 cm, Pk/100. 2388-2020	\$55.60
46 x 57 cm, Pk/100. 2388-4657	\$203.00
0.42 mm thick, 150 mm/30 min absorption, Grade 237 (<i>Whatman 3, 3CHR; S&S 22/2, 591, 2469, 5703</i>)	
46 x 57 cm, Pk/100. 2378-4657	\$233.00
0.83 mm thick, 220 mm/30 min absorption, Grade 222 (<i>Whatman 31ETCHR; S&S 470, 740E</i>)	
7 x 8 cm, Pk/100. 2228-0708	\$17.75
7 x 10 cm, Pk/100. 2228-0710	\$19.75
10 x 15 cm, Pk/100. 2228-1015	\$36.35
20 x 20 cm, Pk/100. 2228-2020	\$56.35
0.95 mm thick, 180 mm/30 min absorption, Grade 243 (<i>Whatman 17CHR; S&S 2668</i>)	
11 x 14 cm, Pk/100. 2438-1114	\$24.05
20 x 20 cm, Pk/100. 2438-2020	\$44.60
46 x 57 cm, Pk/100. 2438-4657	\$342.00
2.48 mm thick, 240 mm/30 min absorption, Grade 320 (<i>S&S 2727</i>)	
11 x 14 cm, Pk/50. 3208-1114	\$38.00
16 x 18 cm, Pk/50. 3208-1618	\$52.95
20 x 20 cm, Pk/50. 3208-2020	\$64.60



Forceps, Filter, 11.5 cm. **300-025** \$4.35



Flask Weights, vinyl-coated lead	
C-Shaped lead rings, blue, Heathrow	
Fits 125-500 mL Flasks, 42 mm ID, 227 g. HS8881A	\$23.75
Fits 250-1000 mL Flasks, 51 mm ID, 472 g. HS8881B	\$27.05
Fits 500-2000 mL Flasks, 54 mm ID, 599 g. HS8881C	\$29.80
O-Shaped lead ring, blue, Heathrow	
Fits 1000-4000 mL Flasks, 70 mm ID, 970 g. HS8882D	\$29.80
Lead Wire, 300 mm long, 114 g. 183240000	\$30.10