

**Spade Tips (Kuehne Pan)**

Angled



Stainless Steel, 11 cm (4 1/4"), economy. **33-1041**..... \$4.15  
(3270)

Teflon® Coated, 11 cm (4 1/4"), Bel-Art. **379420000**..... \$20.90  
Straight, stainless steel, 10.5 cm, economy. **57.926** ..... \$5.85

**Square Point Tips, straight**



Stainless Steel, 12 cm (4 3/4"), economy. **1244**..... \$2.25  
Teflon® Coated, 11 cm (4 1/4"), Bel-Art. **379430000**..... \$22.90

**Locking**



**Forceps, stainless steel**

Self-Closing (cross-locking), squeeze handles to open tips

Fine Point, smooth tips

Economy, 13 cm (5 1/8"). **300-066**..... \$4.95

Peer-Vigor, 11 cm (4 1/4"), Swiss. **57.324**..... \$42.40

Medium Point, serrated tips, 17 cm (6 5/8"). **TWZ-948.00** ..... \$3.95



**Slide-Locking, grid-serrated tips**

Fine Point, Economy, 14 cm (5 1/2"). **57.0357A** ..... \$9.00

Fine Point, Dumont, 16 cm (6 1/4"). **TWZ-101.72** ..... \$45.75

Medium Point, Economy, 16 cm (6 1/4"). **203** ..... \$9.00

Medium Point, Dumont, 16 cm (6 1/4"). **TWZ-101.76**..... \$43.60

Snap-Closing, angled, fine points, 15.5 cm (6"). **3610** ..... \$5.80  
serrated tips & handles, guide pin

**Tissue**




**Forceps, stainless steel, economy**

Adson-Brown, 9 x 9 side grasping teeth, 12 cm. **17-574**..... \$15.75

Adson Tissue, 1 x 2 teeth, 12 cm (4 3/4"). **17-560** ..... \$6.00


Iris Tissue, 1 x 2 teeth, 10 cm (4"). **26-295** ..... \$7.70



Russian Tissue, circular serrated tips, 15.5 cm (6"). **2275**..... \$5.25

Tissue, 1 x 2 teeth, 14 cm (5 1/2"). **17-145** ..... \$4.00

**Tongs**





**Tongs, Crucible, serrated jaws, stainless steel**

24.5 cm (9 3/4"). **320-101**..... \$11.60

51 cm (20"). **320-120**..... \$27.75

**Hemostats & Clamps**



**Hemostats, serrated jaws**

**Hartman Mosquito, stainless steel**

Straight, Economy, 9 cm (3 1/2"). **12-180S** ..... \$6.05

Straight, Standard, 10 cm (4"). **525-025**..... \$25.20

Curved, Economy, 9 cm (3 1/2"). **14-030** ..... \$5.80

Curved, Standard, 10 cm (4"). **20-360**..... \$21.70

**Halstead Mosquito, 12.5 cm (5"), stainless steel**

Straight, Economy. **14-065** ..... \$6.80

Straight, Standard. **20-501** ..... \$22.45

Curved, Economy. **14-070**..... \$6.80

Curved, Standard. **20-519** ..... \$22.85

**Kelly, stainless steel**

Straight, Economy, 14.5 cm (5 3/4"). **14-225** ..... \$6.20

Straight, Standard, 14.5 cm (5 3/4"). **525-225**..... \$22.40

Curved, Economy, 14 cm (5 1/2"). **14-230** ..... \$6.20

Curved, Standard, 14.5 cm (5 3/4"). **20-568**..... \$22.40

**Rochester-Pean, 20.5 cm (8"), stainless steel**

Straight, Economy. **1920** ..... \$6.50

Straight, Standard. **525-435** ..... \$44.10

Curved, Economy. **1930**..... \$6.50


Curved, Standard. **525-440** ..... \$44.10

Polypropylene, 12.5 cm (5"). **HEM-300.00** ..... \$1.50



Specimen Forceps, 25 cm (9 3/4"), serrated jaws. **300-500** ..... \$17.25

Continued




Sponge Forceps, Foerster, serrated jaws, stainless steel  
 Straight  
 17.5 cm (7"), economy. 14-660 ..... \$14.75  
 24 cm (9½"), economy. 14-670 ..... \$17.85  
 Curved, 17.5 cm (7"), economy. 14-665 ..... \$14.75



**Clamps**  
 Dieffenbach Serrefines, stainless steel  
 Straight  
 48 mm (17/8"), economy. 12-247 ..... \$8.90  
 60 mm (23/8"), economy. 300-120 ..... \$6.80  
 Curved  
 48 mm (17/8"), economy. 12-248 ..... \$9.15  
 60 mm (23/8"), economy. 300-121 ..... \$6.80




Schwartz Vessel Clips, 1.7 x 8 mm jaws, 2.5 cm long, Germany  
 Straight Jaws. RS-5452 ..... \$19.00  
 Slightly Angled Jaws. RS-5457 ..... \$19.00



8.6 cm shaft; 5 mm serrated jaw  
 Micro Alligator Forceps, stainless steel. 2769 ..... \$50.00

**Needle Holders & Wound Closure**



Needle Holders, stainless steel  
 Baumgartner, cross-serrated jaws, 14 cm (5½"). 16-285S ..... \$11.05  
 Castroviejo, cross-serrated jaws w/lock, 14 cm. 15-302 ..... \$39.10  
 Crile-Wood, cross-serrated jaws, 15 cm (6")  
 Economy, central groove. 16-293S ..... \$13.45  
 Standard. 32-485 ..... \$33.75  
 German. 320-366 ..... \$60.80  
 Derf, 12 cm (4¾"), economy. 16-238 ..... \$9.70  
 serrated jaws with central groove

Continued

Olsen-Hegar, cross-serrated jaws with scissors, 14 cm (5½")  
 Economy, central groove. 16-296 ..... \$14.30  
 Economy, No Groove. 16-296A ..... \$14.45  
 Standard, central groove. 32-500 (555-160) ..... \$35.25  
 German, central groove. 320-556 ..... \$71.50  
 Webster, smooth jaws, 12 cm (4¾")  
 Economy. 15-035 ..... \$13.30  
 Standard. 555-035 ..... \$42.00  
 Tungsten Carbide. 0817 ..... \$45.00  
 German. 320-101 ..... \$60.50

**Reflex Wound Closure System**

Reflex Clip Applier, USA  
 For 7 mm Clips. 204-1000 ..... \$149.00  
 For 9 mm Clips. 202-1000 ..... \$149.00


Reflex Wound Clips, nonsterile, stainless steel, USA  
 7 mm, 5 packs of 20 clips, Bx/100. 203-1000 ..... \$55.00  
 9 mm, 5 packs of 20 clips, Bx/100. 201-1000 ..... \$50.00

Reflex Clip Remover, fits 7 mm & 9 mm clips, USA. 205-1000 ... \$99.00

**Michel Wound Closure System**




Michel Clip Applier/Remover, double-ended. 3700 ..... \$13.75  
 11.5 cm (4½"), stainless steel




Michel Wound Clips, stainless steel  
 11 mm, Bx/1000. 380-519 ..... \$74.00  
 14 mm, Bx/1000. 380-543 ..... \$89.00  
 16 mm, Bx/1000. 380-568 ..... \$97.00  
 18 mm, Bx/1000. 380-584 ..... \$98.00  
 22 mm, Bx/1000. 380-626 ..... \$116.00

**Inoculating Loops/Needle/Picks & Cell Spreaders**



Inoculating Loops  
 Disposable  
 1 µL Loop, flexible polypropylene  
 Cs/1000, sterile bg/20. HS81121A ..... \$50.20  
 10 µL Loop, polystyrene, Ino-Loop™  
 Cs/1000, sterile bg/25. L200-2 ..... \$64.60  
 Nichrome, 3 mm twisted loop of 25-gauge wire. 3001 ..... \$1.90  
 21 cm aluminum handle  
 Inoculating Needle. 3000 ..... \$1.60  
 76 mm nichrome wire, 21 cm aluminum handle

**Accessories**



Bibulous Paper, 10 x 15 cm, Bk/50 sheets. **BIB-46-50** bk \$3.65 cs/48 bk 169.00

**Forceps**

**Spade Tips**



- Angled (Kuehne)
  - Stainless Steel, 11 cm (4¼"), economy. **3270**..... \$4.15
  - Teflon® Coated, 11 cm (4¼"), Bel-Art. **379420000**..... \$20.90
  - Straight, stainless steel, economy. **57.926**..... \$4.90

**Square Point Tips, straight**



- Stainless Steel, 12 cm (4¾"), economy. **1244**..... \$2.25
- Teflon® Coated, 11 cm (4¼"), Bel-Art. **379430000**..... \$22.90



Minimize contamination during slide & coverslip handling. Will not scratch; corrosion and heat resistant.

Hemostat, Polypropylene, 12.5 cm. **HEM-300.00** ..... \$1.50

**Cytology Funnels**

Single, Metal Clips Required



Double, All Plastic



White filter cards will accept up to 0.5 mL of sample, while tan filters will accept volumes up to 0.4 mL.

Cytology Funnels, Simport CytoSep™

**Shandon (Thermo) Cytospin® Cytocentrifuges**

- Single. *Metal slide clips required*
  - White Filter & Cap, Cs/500. **M964-10FW** ..... \$382.00  
(59 91 040, 110 25 48)
  - Tan Filter & Cap, Cs/500. **M964-10FT** ..... \$382.00
- Single. *All plastic - metal slide clips not required*
  - White Filter & Cap, Cs/480. **M965-10FW** ..... \$503.00  
(A787 100003)
  - Tan Filter & Cap, Cs/480. **M965-10FT** ..... \$503.00
- Double. *All plastic - metal slide clips not required*
  - White Filter & Cap, Cs/480. **M965-20FW** ..... \$626.00  
(A787 100005)

**Immersion Oils & Mounting Media**



Cargille immersion oils A, B, and NVH can be mixed with each other for custom viscosities. USA

**Immersion Oil**

**Type A**, low viscosity (150 cSt). *Short W.D. optics*

- 7.5 mL in LDPE dispensing bottle. **18101** ..... \$5.35
- 20 mL in HDPE dispensing bottle. **18102** ..... \$10.00
- 120 mL in amber bottle. **18103** ..... \$22.45
- 480 mL in amber bottle. **18104** ..... \$49.60



**Type B**, high viscosity (1250 cSt). *Long W.D. optics*

- 7.5 mL in LDPE dispensing bottle. **18105** ..... \$5.35
- 20 mL in HDPE dispensing bottle. **18106** ..... \$10.00
- 120 mL in amber bottle. **18107** ..... \$22.45
- 480 mL in amber bottle. **18108** ..... \$49.60


**Type NVH**, very high viscosity (21,000 cSt)

- 480 mL in amber bottle. **18113** ..... \$57.75

**Type FF**, fluorescence free (170 cSt)

- 7.5 mL in LDPE dispensing bottle. **18121** ..... \$7.35
- 480 mL in amber bottle. **18124** ..... \$81.75

A long nozzle, textured sides, and slim profile make this chemically inert, high-density polyethylene (HDPE) container ideal for dispensing immersion oil.



Bottle, Dispensing, HDPE, 20 mL, each. **18145** ..... \$6.90



Store immersion oil or mounting media in this borosilicate glass bottle with an applicator rod and ground glass cap.

Bottle, Balsam, glass, 45 mL (empty), USA. **616**..... \$17.30

K-750



Suitable for fluorescence microscopy.

K-700



Takes only 10 minutes to dry.

**Mounting Media**

Mount-Quick, Kiyota

- 30 mL in Dropper Bottle, aqueous-based. **K-750**..... \$22.40
- 30 mL in Tube, solvent-based, 30 cps. **K-700**..... \$15.40