


**Rongeurs & Bone Instruments**




Blumenthal Rongeur

- Rongeurs, stainless steel**  
 Blumenthal, 15 cm (5¾")  
 Economy. **410-571** ..... \$39.20  
 Standard. **650-571** ..... \$94.40  
 Blumenthal, Slightly Angled Jaws, 15 cm (5¾"). **9571** ..... \$37.50  
 Mini-Blumenthal, 11.5 cm (4½"). **410-560** ..... \$50.00  
 Friedman, 14 cm (5½")  
 Economy. **410-574** ..... \$39.20  
 Standard. **650-574** ..... \$94.40




Liston Bone Cutter


- Bone Cutters, stainless steel, economy**  
 Liston  
 14 cm (5½"), 20 mm blades. **320-086** ..... \$23.80  
 21.5 cm (8½"), 30 mm blades. **320-087 (1095)** ..... \$28.75  
 Scissor Style, 10 cm (4"), 18 mm blades. **320-025** ..... \$8.75  
 straight, strong tips




- Bone Nipper, 10 cm (4") long, 6 mm cutting edge. B286**..... \$12.00  
 angled offset jaws



- Bone Shears, Constable, 23 cm (9") long, 50 mm (2") cutting jaws**  
 Economy. **33-1001** ..... \$45.00  
 German. **965-620** ..... \$160.00



- Postmortem Hammer, hook handle, 25 cm (10") long**  
 Economy. **33-1032** ..... \$35.70  
 Standard. **74-043** ..... \$80.00



16 mm (5/8") cutting blade;  
 9 cm (3½") long

- Skull Breaker, Virchow, economy. 33-1011**..... \$42.40

**Scalpels**



Economy

**Scalpels**

- Disposable, plastic handle, sterile  
 Economy, gray handle
- |  | 1 bx   | cs/50 bx |
|--|--------|----------|
| ★ #10 carbon steel blade, Bx/10. <b>9410</b> ..... | \$9.60 | 416.00   |
| #11 carbon steel blade, Bx/10. <b>9411</b> .....   | \$9.60 | 416.00   |
| #22 carbon steel blade, Bx/10. <b>9422</b> .....   | \$9.60 | 416.00   |
- Lance, green handle, England  
 #10 carbon steel blade, Bx/10. **92410** ..... \$13.00  
 #11 carbon steel blade, Bx/10. **92411** ..... \$13.00  
 #21 carbon steel blade, Bx/10. **92421** ..... \$13.80  
 #22 carbon steel blade, Bx/10. **92422** ..... \$13.80



- Miniature, plastic thumb handle, nonsterile, England  
 #11P stainless steel blade, each. **2201111**..... \$0.95




Smooth Handle




Fluted Handle

- Reusable, one-piece stainless steel, 1½" blade  
 Smooth Handle, 15.5 cm (6"). **14-350** ..... \$3.15  
 (1445)  
 Fluted Handle, narrow blade, 15.5 cm (6"). **14-360**..... \$3.15  
 (1460)

Continued



- Scalpel Blades**, individually wrapped
- Carbon Steel, Economy, Bx/100 ..... \$18.95
    - #10 blade. **9310**                      #21 blade. **9321**
    - #11 blade. **9311**                      #22 blade. **9322**
    - #12 blade. **9312**                      #23 blade. **9323**
    - #15 blade. **9315**                      #24 blade. **9324**
    - #20 blade. **9320**
  - Carbon Steel, Bard-Parker, Rib-Back® design, nonsterile, USA
    - #60 blade, 25 strips of 6, Bx/150. **371340** ..... \$343.00
  - Carbon Steel, Germany, nonsterile
    - #60 blade, 25 strips of 4, Bx/100. **2200260** ..... \$46.00
    - #70 blade, 25 strips of 4, Bx/100. **2200270** ..... \$63.00
  - Carbon Steel, Swann-Morton, nonsterile
    - #60 blade, Bx/50. **22002SM60** ..... \$29.60
  - Stainless Steel, Economy, Bx/100 ..... \$21.15
    - #10 blade. **36-010S**                      #21 blade. **36-021S**
    - #11 blade. **36-011S**                      #22 blade. **36-022S**
    - #20 blade. **36-020S**
  -  Stainless Steel, Lance, Bx/100 ..... \$38.00
    - #11 blade. **92111**                      #15 blade. **92115**
  - Stainless Steel, Economy
    - #60 blade, Bx/100. **052-620** ..... \$45.00
    - #70 blade, Bx/100. **052-622** ..... \$45.00



Safely remove and store up to 100 used blades in this disposable box.

- Scalpel Blade Remover. 22048** ..... \$3.60
- Stainless Steel, Solid, #3 Style      Stainless Steel, Solid, #4 Style
  - Stainless Steel, Solid, #7 Style      Stainless Steel, Hollow
  - PVC Coated Handle, #3 Style      Plastic Handle, #5/6/8 Styles


- Scalpel Handles**, reusable
- Stainless Steel, solid
    - #3 Style, accepts blades #10-17
      - Economy. **36-405** ..... \$3.80
      - Standard. **10-305** ..... \$11.50
    - #4 Style, accepts blades #18-36
      - Economy. **36-410** ..... \$3.80
      - Standard. **565-410** ..... \$11.20
    - #7 Style, accepts blades #10-17
      - Economy. **36-415** ..... \$4.50
      - Standard. **565-415** ..... \$14.00

- Stainless Steel, Hollow
  - #8 Style, accepts blades #60 & #70. **HAN8** ..... \$8.00
- PVC Coated Stainless Steel. *Not autoclavable*
  - #3 Style, accepts blades #10-17. **KNF-225.00** ..... \$5.35
- Plastic, contoured to fit palm of hand
  - #5 Style, accepts blades #10-17. **KNF-221.00** ..... \$5.00
  - #6 Style, accepts blades #18-36. **1605** ..... \$5.00
  - #8 Style, accepts blades #60 & #70
    - Economy. **HAN8P** ..... \$6.00
    - Bard-Parker™, USA. **371080** ..... \$31.35


**Scissors**



- Scissors, Spring**, ambidextrous, stainless steel
- Barraquer (butterfly), 6 cm long, economy. **3737** ..... \$34.50
  - 10 mm angled-on-flat blades



- Castroviejo, straight, sharp tips, economy
  - 7 cm (2¾"), 5 mm blades. **36205** ..... \$45.00
  - Effective cutting edge is 3 mm; tip-to-tip is 3 mm.
  - 8 cm (3"), 7 mm blades. **5-254** ..... \$45.00
  - Effective cutting edge is 6 mm; tip-to-tip is 4 mm.
  - 8.5 cm (3¼"), 5 mm blades. **5-855** ..... \$45.00
  - Effective cutting edge is 4 mm; tip-to-tip is 3 mm.
  - 11.5 cm (4½"), 17 mm blades. **M3720** ..... \$26.00
  - Effective cutting edge is 12 mm.



- Uniband™      Lifetime guarantee against rust and defects.
- Featherlite™
- Uniband™, sharp tips, Laschal Surgical, USA
  - Straight, 12 mm Blades. **LA-4** ..... \$32.00
  - Straight, 20 mm Blades. **LA-1** ..... \$32.00
  - Straight, 20 mm Serrated Blades. **LA-2** ..... \$32.00
  - Curved, 12 mm Blades. **LA-4C** ..... \$32.00
  - Curved, 20 mm Blades. **LA-1C** ..... \$32.00
- Featherlite™, Castroviejo, sharp tips, Laschal Surgical, USA
  - Greater tactile sensitivity than Uniband™ scissors.
  - Straight, 12 mm Blades. **MPF-4** ..... \$80.80
  - Straight, 20 mm Blades. **MPF-1** ..... \$80.80
  - Curved, 12 mm blades. **MPF-4C** ..... \$80.80

- Scissors, Bonn**, 15 mm blades, 9 cm (3½"), stainless steel
- Straight, German. **744-074** ..... \$73.15
  - Curved, German. **744-076** ..... \$76.90

Continued