


- Scalpel Blades**, individually wrapped
- Carbon Steel, Economy, Bx/100 \$18.95
 - #10 blade. **9310** #21 blade. **9321**
 - #11 blade. **9311** #22 blade. **9322**
 - #12 blade. **9312** #23 blade. **9323**
 - #15 blade. **9315** #24 blade. **9324**
 - #20 blade. **9320**
 - Carbon Steel, Bard-Parker, Rib-Back® design, nonsterile, USA
 - #60 blade, 25 strips of 6, Bx/150. **371340** \$343.00
 - Carbon Steel, Germany, nonsterile
 - #60 blade, 25 strips of 4, Bx/100. **2200260** \$46.00
 - #70 blade, 25 strips of 4, Bx/100. **2200270** \$63.00
 - Carbon Steel, Swann-Morton, nonsterile
 - #60 blade, Bx/50. **22002SM60** \$29.60
 - Stainless Steel, Economy, Bx/100 \$21.15
 - #10 blade. **36-010S** #21 blade. **36-021S**
 - #11 blade. **36-011S** #22 blade. **36-022S**
 - #20 blade. **36-020S**
 -  Stainless Steel, Lance, Bx/100 \$38.00
 - #11 blade. **92111** #15 blade. **92115**
 - Stainless Steel, Economy
 - #60 blade, Bx/100. **052-620** \$45.00
 - #70 blade, Bx/100. **052-622** \$45.00



Safely remove and store up to 100 used blades in this disposable box.

- Scalpel Blade Remover. 22048** \$3.60
- Stainless Steel, Solid, #3 Style Stainless Steel, Solid, #4 Style
 - Stainless Steel, Solid, #7 Style Stainless Steel, Hollow
 - PVC Coated Handle, #3 Style Plastic Handle, #5/6/8 Styles


- Scalpel Handles**, reusable
- Stainless Steel, solid
 - #3 Style, accepts blades #10-17
 - Economy. **36-405** \$3.80
 - Standard. **10-305** \$11.50
 - #4 Style, accepts blades #18-36
 - Economy. **36-410** \$3.80
 - Standard. **565-410** \$11.20
 - #7 Style, accepts blades #10-17
 - Economy. **36-415** \$4.50
 - Standard. **565-415** \$14.00

- Stainless Steel, Hollow
 - #8 Style, accepts blades #60 & #70. **HAN8** \$8.00
- PVC Coated Stainless Steel. *Not autoclavable*
 - #3 Style, accepts blades #10-17. **KNF-225.00** \$5.35
- Plastic, contoured to fit palm of hand
 - #5 Style, accepts blades #10-17. **KNF-221.00** \$5.00
 - #6 Style, accepts blades #18-36. **1605** \$5.00
 - #8 Style, accepts blades #60 & #70
 - Economy. **HAN8P** \$6.00
 - Bard-Parker™, USA. **371080** \$31.35


Scissors



- Scissors, Spring**, ambidextrous, stainless steel
- Barraquer (butterfly), 6 cm long, economy. **3737** \$34.50
 - 10 mm angled-on-flat blades



- Castroviejo, straight, sharp tips, economy
 - 7 cm (2¾"), 5 mm blades. **36205** \$45.00
 - Effective cutting edge is 3 mm; tip-to-tip is 3 mm.
 - 8 cm (3"), 7 mm blades. **5-254** \$45.00
 - Effective cutting edge is 6 mm; tip-to-tip is 4 mm.
 - 8.5 cm (3¼"), 5 mm blades. **5-855** \$45.00
 - Effective cutting edge is 4 mm; tip-to-tip is 3 mm.
 - 11.5 cm (4½"), 17 mm blades. **M3720** \$26.00
 - Effective cutting edge is 12 mm.



- Uniband™, sharp tips, Laschal Surgical, USA
 - Lifetime guarantee against rust and defects.
 - Straight, 12 mm Blades. **LA-4** \$32.00
 - Straight, 20 mm Blades. **LA-1** \$32.00
 - Straight, 20 mm Serrated Blades. **LA-2** \$32.00
 - Curved, 12 mm Blades. **LA-4C** \$32.00
 - Curved, 20 mm Blades. **LA-1C** \$32.00
- Featherlite™, Castroviejo, sharp tips, Laschal Surgical, USA
 - Greater tactile sensitivity than Uniband™ scissors.
 - Straight, 12 mm Blades. **MPF-4** \$80.80
 - Straight, 20 mm Blades. **MPF-1** \$80.80
 - Curved, 12 mm blades. **MPF-4C** \$80.80


- Scissors, Bonn**, 15 mm blades, 9 cm (3½"), stainless steel
- Straight, German. **744-074** \$73.15
 - Curved, German. **744-076** \$76.90

Continued



Scissors, Delicate, stainless steel

Straight, sharp tips, 22 mm blades	
9 cm (3½"), Standard. 13-001	\$19.00
10.5 cm (4"), Standard. 13-003	\$21.30
10 cm (4"), German. 103-500	\$44.70
Curved, sharp tips, 22 mm blades	
10 cm (4"), Standard. 13-059	\$22.45
10 cm (4"), German. 103-503	\$46.90




Scissors, Dissecting, 11.5 cm (4½"), stainless steel

Straight, sharp tips	
25 mm Blades, Economy. 12-410	\$5.05
25 mm Blades, Standard. 13-002	\$19.00
25 mm Blades, German. 103-259	\$38.50
30 mm Blades, Economy. 310-040	\$4.50
35 mm Blades, Economy. 310-045	\$4.35
Curved, sharp tips	
25 mm Blades, Economy. 310-035	\$4.35
25 mm Blades, Standard. 13-010	\$19.40




Angled, 25 mm angled-on-side blades, economy. 310-066 \$6.00



Scissors, Deaver, 14 cm (5½"), stainless steel

Straight, S/S (sharp/sharp) tips, Standard. 10-632	\$21.80
Straight, S/B (sharp/blunt) tips, Standard. 10-630	\$21.30



Scissors, Mayo, stainless steel

Straight	
14 cm (5½"), Economy. 310-086	\$9.15
14.5 cm (5¾"), Standard. 15-575	\$22.55
14.5 cm (5¾"), Standard-TC. 16-851	\$51.50
17.5 cm (6¾"), Economy. 310-088	\$11.25
17.5 cm (6¾"), Standard. 15-579	\$26.95
17.5 cm (6¾"), Standard-TC. 16-854	\$58.50
17.5 cm (6¾"), German-TC. 126-755	\$123.00
23 cm (9"), Standard. 15-585	\$42.35


Continued

Scissors, Mayo (cont)

Curved	
14 cm (5½"), Economy. 310-087	\$9.65
14.5 cm (5¾"), Standard. 15-576	\$22.95

Scissors, Mayo-Stille, rounded blades, stainless steel

Straight	
14 cm (5½"), Standard. 15-577	\$25.85
17 cm (6¾"), Standard. 15-583	\$36.85



Scissors, Operating, stainless steel

Straight, S/S (sharp/sharp) tips	
11.5 cm (4½"), Standard. 10-456	\$17.85
12 cm (4¾"), German. 101-451	\$33.60
14 cm (5½"), Economy. 12-055	\$5.15
14 cm (5½"), LR (large rings), Economy. 11-107SSLR	\$5.20
14 cm (5½"), Standard. 10-552	\$21.30
14 cm (5½"), Standard-TC. 16-800	\$47.50
14 cm (5½"), German. 101-550	\$36.50
15.5 cm (6"), Economy. 12-060	\$5.15
15.5 cm (6"), Standard. 510-060	\$23.80
16.5 cm (6½"), Economy. 12-065	\$6.75
16.5 cm (6½"), Standard. 10-602	\$24.75
16.5 cm (6½"), German. 101-600	\$40.20

Straight, S/B (sharp/blunt) tips

11.5 cm (4½"), Economy. 12-015	\$4.90
12 cm (4¾"), Standard. 10-453	\$17.85
12 cm (4¾"), German. 100-453	\$33.60
14 cm (5½"), Economy. 12-025	\$5.15
14 cm (5½"), Standard. 10-550	\$21.30
14 cm (5½"), Standard-TC. 16-805	\$47.60
14 cm (5½"), German. 100-552	\$36.50
16.5 cm (6½"), Economy. 12-035	\$6.75
16.5 cm (6½"), Standard. 10-600	\$24.75
16.5 cm (6½"), German. 100-602	\$40.20

Straight, B/B (blunt/blunt) tips

11.5 cm (4½"), Economy. 12-075	\$4.90
14 cm (5½"), Economy. 12-085	\$5.15
14 cm (5½"), Standard. 10-558	\$21.30
14 cm (5½"), Standard-TC. 16-815	\$47.50
14 cm (5½"), German. 102-558	\$36.50
16.5 cm (6½"), Economy. 12-095	\$6.75
16.5 cm (6½"), Standard. 10-608	\$24.75
16.5 cm (6½"), German. 102-608	\$40.20


Curved, S/S (sharp/sharp) tips

14 cm (5½"), Economy. 12-155	\$5.60
14.5 cm (5¾"), Standard. 10-566	\$21.80
16.5 cm (6½"), Economy. 12-165	\$7.10
16.5 cm (6½"), Standard. 10-616	\$25.25

Continued

- Curved, S/B (sharp/blunt) tips
- 11.5 cm (4½"), Economy. 12-115.....\$4.90
 - 14 cm (5½"), Economy. 12-125.....\$5.60
 - 14 cm (5½"), Standard. 10-560.....\$21.80
 - 14 cm (5½"), German. 100-560.....\$37.90
 - 16.5 cm (6½"), Economy. 12-135.....\$7.10
 - 16.5 cm (6½"), Standard. 10-610.....\$25.25
 - 16.5 cm (6½"), German. 100-610.....\$41.30

- Curved, B/B (blunt/blunt) tips
- 11.5 cm (4½"), Economy. 12-175.....\$4.90
 - 14 cm (5½"), Economy. 12-185.....\$5.60
 - 14 cm (5½"), Standard. 10-568.....\$21.80
 - 16.5 cm (6½"), Economy. 12-195.....\$7.10
 - 16.5 cm (6½"), Standard. 10-618.....\$25.25




Scissors, Metzenbaum, stainless steel

- Straight
- 15 cm (5¾"), Economy. 12-225.....\$9.00
- Curved
- 14.5 cm (5¾"), Economy. 12-230.....\$9.35
 - 14.5 cm (5¾"), Standard. 15-552.....\$23.05
 - 14 cm (5½"), German. 105-569.....\$46.20




- Scissors, Stevens Tenotomy, 11.5 cm (4½"), stainless steel**
- Straight, blunt tips, Standard. 13-520.....\$24.75




Scissors, Strabismus, stainless steel

- Straight
- 10 cm (4"), Economy. 26-105.....\$6.25
 - 11.5 cm (4½"), Standard. 13-524.....\$23.65
- Curved
- 11 cm (4¼"), Economy. 26-110.....\$6.25
 - 11.5 cm (4½"), Standard. 13-526.....\$24.15



- Scissors, Student, 11.5 cm (4½"), steel wire, USA. 22.....\$0.90**



- Scissors, Bandage, stainless steel, economy**
- Knowles, straight, sharp probe tips, 14 cm (5½"). 1070.....\$6.25
 - Lister, angled, blunt probe tips, 14 cm (5½"). 1050.....\$5.00
- Can be used as enterotomy scissors.



Cut through everything from wire to bone with these durable scissors. Bottom blade is serrated and has a footplate for safety. Blades are stainless steel & handles are plastic. Autoclavable to 120° C.

- Utility Shears, 18 cm (7¼"), economy. 310-341.....\$5.25**

Scrapers

Remove paraffin, ice, etc., safely from your laboratory equipment with these plastic scrapers.



- Scrapers, plastic**
- 5 cm (2") Wide Blade. 311-245.....\$3.70
 - 10 cm (4") Wide Blade. 311-250.....\$4.70
 - 15 cm (6") Wide Blade. 311-255.....\$5.70

Spatulas & Scoops



- Spatulas, Hardwood Handle, stainless steel blade**
- 10 cm (4") Blade. 14-365B.....\$4.25
 - 15 cm (6") Blade. 14-365D.....\$4.80
 - 20 cm (8") Blade. 20-497.....\$6.05



- Spatulas, Heavy-Duty, stainless steel**
- Round End/Square End, both flat
 - 30.5 cm (12"). 320-425.....\$14.95
 - 46 cm (18"). 320-430.....\$18.90