


- Curved, S/B (sharp/blunt) tips
- 11.5 cm (4½"), Economy. 12-115.....\$4.90
 - 14 cm (5½"), Economy. 12-125.....\$5.60
 - 14 cm (5½"), Standard. 10-560.....\$21.80
 - 14 cm (5½"), German. 100-560.....\$37.90
 - 16.5 cm (6½"), Economy. 12-135.....\$7.10
 - 16.5 cm (6½"), Standard. 10-610.....\$25.25
 - 16.5 cm (6½"), German. 100-610.....\$41.30

- Curved, B/B (blunt/blunt) tips
- 11.5 cm (4½"), Economy. 12-175.....\$4.90
 - 14 cm (5½"), Economy. 12-185.....\$5.60
 - 14 cm (5½"), Standard. 10-568.....\$21.80
 - 16.5 cm (6½"), Economy. 12-195.....\$7.10
 - 16.5 cm (6½"), Standard. 10-618.....\$25.25




Scissors, Metzenbaum, stainless steel

- Straight
- 15 cm (5¾"), Economy. 12-225.....\$9.00
- Curved
- 14.5 cm (5¾"), Economy. 12-230.....\$9.35
 - 14.5 cm (5¾"), Standard. 15-552.....\$23.05
 - 14 cm (5½"), German. 105-569.....\$46.20




Scissors, Stevens Tenotomy, 11.5 cm (4½"), stainless steel

- Straight, blunt tips, Standard. 13-520.....\$24.75



Scissors, Strabismus, stainless steel

- Straight
- 10 cm (4"), Economy. 26-105.....\$6.25
 - 11.5 cm (4½"), Standard. 13-524.....\$23.65
- Curved
- 11 cm (4¼"), Economy. 26-110.....\$6.25
 - 11.5 cm (4½"), Standard. 13-526.....\$24.15



- Scissors, Student, 11.5 cm (4½"), steel wire, USA. 22.....\$0.90



Scissors, Bandage, stainless steel, economy

- Knowles, straight, sharp probe tips, 14 cm (5½"). 1070.....\$6.25
 - Lister, angled, blunt probe tips, 14 cm (5½"). 1050.....\$5.00
- Can be used as enterotomy scissors.



Cut through everything from wire to bone with these durable scissors. Bottom blade is serrated and has a footplate for safety. Blades are stainless steel & handles are plastic. Autoclavable to 120° C.

- Utility Shears, 18 cm (7¼"), economy. 310-341.....\$5.25

Scrapers

Remove paraffin, ice, etc., safely from your laboratory equipment with these plastic scrapers.



Scrapers, plastic

- 5 cm (2") Wide Blade. 311-245.....\$3.70
- 10 cm (4") Wide Blade. 311-250.....\$4.70
- 15 cm (6") Wide Blade. 311-255.....\$5.70

Spatulas & Scoops




Spatulas, Hardwood Handle, stainless steel blade

- 10 cm (4") Blade. 14-365B.....\$4.25
- 15 cm (6") Blade. 14-365D.....\$4.80
- 20 cm (8") Blade. 20-497.....\$6.05



Spatulas, Heavy-Duty, stainless steel

- Round End/Square End, both flat
- 30.5 cm (12"). 320-425.....\$14.95
 - 46 cm (18"). 320-430.....\$18.90



Spatulas, Micro, stainless steel

- Square End/Spoon End
 - 19 cm (7½"). 320-411 \$4.15
 - 23 cm (9"). 320-414 \$4.25
- Round End/Tapered End, 18 cm (7") long. 320-370 \$4.65
 - Round end: 8 x 51 mm; Tapered end: 6 x 51 mm
- V-Shaped End/Round End, 16 cm (6¼") long. HS15906 \$5.25



Aluminum



Polycarbonate

Scoops

- Aluminum, ribbed handle, round
 - 350 mL (12 oz). 105ASC12 \$10.45
 - 700 mL (24 oz). 105ASC24 \$13.45
 - 1100 mL (38 oz). 105ASC38 \$17.95
 - 1700 mL (58 oz). 105ASC58 \$23.95
- Polycarbonate (PC), transparent, rectangular
 - 175 mL (6 oz). 1052882 \$9.95
 - 945 mL (32 oz). 1059F75 \$17.95
 - 1890 mL (64 oz). 1059F76 \$19.95



Stainless Steel

- V-Shaped (Reagent Digger), 16.5 cm (6½"). 600LSCSS6 \$1.70
- Handle, wooden, 13 cm long, each. 600LSHND1 \$2.00

Teasing Needles, Probes, Picks & Pin Vises

97-1302




Teasing Needles

- Angular, wooden handle, 14.5 cm (5¾"), Pk/12. 97-1302 \$4.15
- Straight, plastic handle, 14.5 cm (5¾"), Pk/12. 97-1313 \$3.90



- Stainless Steel, straight, 14 cm. M28 \$4.00

Blunt tip is angled and tapered. Hexagonal handle is 6 mm wide.



Probes, Light Mall, 15.5 cm (6")

- Nickel Plated Steel. 15-800 \$2.90
(5800)
- Stainless Steel. 15-800SS \$3.80



Probes, Bowman Lacrimal, double-ended

- Economy, sterling silver, 12.7 cm (5") long
 - Size 0000-000 (~0.4 – 0.5 mm). 24-706 \$11.65
 - Size 00-0 (~0.6 – 0.7 mm). 24-708 \$11.65
 - Size 1-2 (~0.9 – 1.0 mm). 24-710 \$11.65
 - Size 3-4 (~1.1 – 1.2 mm). 24-712 \$11.65



Probe/Seekers, stainless steel

- Double-Ended, 15.5 cm. 320-216A \$4.50
 - One end tapered to a point; one end curved with a rounded tip.
- Single End, 12.5 cm (5"). 320-217A \$4.40
 - Curved, fine point tip with hexagonal handle.



Probe with Bulb Tips, 25 cm (10"), stainless, German. 541-807 \$17.50



Probe with Eye & Bulb Tip, flexible, 15.5 cm (6"). 23-496 \$2.25 (5850)

Picks (Explorers) Call regarding availability

PIN-220.00



PIN-212.00



Hold needles, probes, wires, drill bits, etc.

Pin Vise, Swivel-Head, 0-3.3 mm opening. PIN-220.00 \$7.95
Two double-ended collets are provided; one in screw chuck, one inside handle.

Pin Vise, Mini, 0-0.8 mm opening. PIN-212.00 \$4.95
blackened tool steel