

Inverted

Jenco

- CP-2A1S, Trinocular, Brightfield/Phase \$2,295.00
- Widefield 10X/20 eyepieces
 - Trinocular, 30° inclined, Siedentopf design viewing head with single diopter; 55-75 mm interpupillary adjustment
 - Three DIN, LWD plan achromat objectives: brightfield 4X; plan phase 10X and 25X
 - 5-place, reverse facing nosepiece
 - Fixed stage
 - Coaxial coarse and fine focus controls; coarse has tension and up-stop adjustments
 - 1.25 N.A. Abbe condenser with rack and pinion focusing, iris diaphragm and filter holder; working distance of 55 mm for phase contrast, or 125 mm if swung out for brightfield
 - Centerable, dedicated phase annuli mounted on a slider
 - 30 watt halogen Koehler-type illuminator
 - 48.25 cm (19")D x 21.6 cm (8½")W x 47 cm (18½")H
 - 5-year mechanical & optical/2-year electrical warranty



CP-2A1

Jenco

- CP-2A1, Trinocular, Brightfield/Phase..... \$2,975.00
- Widefield 10X/20 eyepieces
 - Trinocular, 30° inclined, Siedentopf design viewing head with single diopter; 55-75 mm interpupillary adjustment
 - Six DIN, LWD plan achromat objectives: brightfield 10X, 25X, and 40X; plan phase 10X, 25X, and 40X
 - 5-place, reverse facing nosepiece
 - Attachable mechanical stage includes three inserts
 - Coaxial coarse and fine focus controls; coarse has tension and up-stop adjustments
 - 1.25 N.A. Abbe condenser with rack and pinion focusing, iris diaphragm and filter holder; working distance of 55 mm for phase contrast, or 125 mm if swung out for brightfield
 - Centerable, dedicated phase annuli mounted on sliders
 - 30 watt halogen Koehler-type illuminator
 - 48.25 cm (19")D x 21.6 cm (8½")W x 47 cm (18½")H
 - 5-year mechanical & optical/2-year electrical warranty



BC-3A1

Jenco

- BC-3A1, Trinocular, Brightfield/Phase \$5,100.00



BC-366

Jenco

- BC-366, Trinocular, BF/Phase/Fluor \$8,960.00
- Widefield 10X/22 eyepieces
 - Binocular, 45° inclined, Siedentopf design viewing head with single diopter; 55-75 mm interpupillary adjustment
 - Side camera port
 - Six infinity, LWD plan achromat objectives: brightfield/fluorescence 10X, 25X, and 40X; phase 10X, 25X, and 40X
 - 6-place, angled nosepiece
 - Attachable mechanical stage
 - Coaxial coarse and fine focus controls; coarse has tension and up-stop adjustments
 - 1.25 N.A. Abbe condenser with rack and pinion focusing, iris diaphragm and filter holder; centerable, dedicated phase annuli mounted in turret
 - 30 watt halogen transmitted light; Koehler-type illuminator
 - 100 watt mercury epi-fluorescence
 - Blue (FITC) and green (TRITC) fluorescence filter sets
 - 5-year mechanical & optical/2-year electrical warranty

Meiji Access www.meijitechno.com
and call us to place an order

Polarized Light



Meiji Access www.meijitechno.com
and call us to place an order

Metallurgical



Jenco Access www.jencointernational.com
and call us to place an order

Meiji Access www.meijitechno.com
and call us to place an order