



Rapidly decontaminate metal dissecting instruments with this dry heat sterilizer. Glass beads are heated to a constant 500°F (260°C) in a 2" diameter x 4" deep stainless steel well. Sterilize instrument tips in 15 seconds. 1-year warranty.

**Sterilizer, Hot Bead, Germinator™ 500, USA. 5-1450**..... \$419.00  
 Glass Beads (replacement), 1.5 mm, Bg/10 oz. 5-1455..... \$21.00



Safely incinerate microorganisms in 5-7 seconds at 815°C (1500°F) within this ceramic, stainless steel-encased chamber. Compact and flameless operation allows this unit to be used in anaerobic and aerobic chambers. Ideal for sterilizing loops and needles, glass tube/pipette mouths, and various glass and metal instruments. 1-year warranty.



**Sterilizer, Infrared Heat, BactiZapper™. B1000** ..... \$312.00



- Blue film/white paper
- EO gas/steam
- Latex-free
- Self-sealing


**Sterilization Pouches**

3½ x 9" (8.9 x 22.9 cm), Bx/200. **440166423** ..... \$19.95  
 5¼ x 10" (13.3 x 25.4 cm), Bx/200. **440166424** ..... \$29.95  
 7½ x 13" (19.1 x 33 cm), Bx/200. **440166425** ..... \$52.15



**Sterilizing Trays, plastic, autoclavable, PST, USA**

B/L/M (base, lid & mat) included.  
 6 x 2½ x ¾" (15.2 x 6.4 x 1.9 cm). **2600A** ..... \$20.00  
 (500252)  
 ★ 6½ x 4 x ¾" (16.2 x 10.2 x 1.9 cm). **4600A** ..... \$29.00  
 10 x 6 x ¾" (25.4 x 15.2 x 1.9 cm). **6100A** ..... \$60.00  
 (500254)



B//L/2M (base, insert tray, lid & 2 mats) included.  
 15 x 10 x 1½" (38.1 x 25.4 x 3.8 cm). **1520A** ..... \$398.00




**Instrument Trays**

Polypropylene, Deep w/Flat Cover, autoclavable  
 8½ x 6 x 2¾" (22 x 15 x 7 cm). **60081736**..... \$9.80



**Stainless Steel**

Shallow (Mayo), USA  
 13¾ x 9¾ x ¾" (34.8 x 24.9 x 1.9 cm). **13F** ..... \$42.65  
 15½ x 10½ x ¾" (38.4 x 27 x 1.9 cm). **15F** ..... \$50.65  
 Deep w/Dome Cover  
 8¾ x 3¼ x 1½" (22.2 x 8.4 x 3.8 cm). **35-309**..... \$31.25  
 Deep - Order trays and covers separately, USA  
 8⅞ x 5 x 2" (22.7 x 12.9 x 5.1 cm). **952**..... \$27.85  
 Cover, Fits 952 Tray, dome w/handle. **952D-2**..... \$44.00  
 17⅞ x 4½ x 2⅞" (43.6 x 11.2 x 5.4 cm). **170-0**..... \$49.75  
 Cover, Fits 170-0 Tray, dome w/handle. **170-2** ..... \$39.30



Arkansas stones outlast and produce a better finish than man-made stones.

**Sharpening Stone, natural Arkansas**

Fine Grit (hard), mounted on wood. **ABR-800.10** ..... \$23.95  
 Dimensions: 4 x 2 x ½"



**Ultrasonic Cleaner**

Mechanical Timer, 600 mL capacity. **CLN-520.20** ..... \$45.95  
 Tank dimensions: 6 x 3½ x 2"(D). 42 kHz of power.  
 Includes a plastic cleaning basket.

**Accessories**

Cleaning Basket, fine mesh, 5 x 4 x 1¾". **CLN-650.10** ..... \$14.95  
 Cleaning Solution, 1 Quart. **CLN-850.10**..... \$9.95  
 Biodegradable concentrate mixes 1:30.