

20-Slide
Cardboard
Folder



10-Slide/
10-Block
Consulting
Folder



Slide Folders

- 10-Slide/10-Block (consulting), each. **CSF-101** \$24.75
 12-Slide, holds 75 x 50 mm (3 x 2") slides
 Cardboard, black tape, USA. **SF12-BLACK** \$7.20
 20-Slide, holds 75 x 25 mm (3 x 1") slides
 Cardboard, USA
 Folders made in the USA resist heat-induced warping.
 Black Tape. **SF20-BLACK** \$6.25
 Red Tape. **SF20-RED** \$6.25
 Green Tape. **SF20-GREEN** \$6.25
 Navy Blue Tape. **SF20-BLUE** \$6.25
 Cardboard, Black Tape, China. **HS9902** \$4.70
 (12-587-10)
 Plastic, gray tray w/2 clear doors. **M750-20GY** \$9.20
 Transparent covers hold visible slides securely in place,
 but can fold flat against back of tray when opened.
 Slides are slightly elevated for easy retrieval.




Slide Holders, perfect for slide drying or staining

- 12-Slide (horizontally positioned), polystyrene. **445570000** \$7.80
 52-Slide (approx.), 12 x 4", HDPE, USA. **L500-C** \$27.20
 91-Slide (approx.), 12 x 7", HDPE, USA. **L500-D** \$36.35



2-Slide



5-Slide

Slide Mailers, Plastic

- 2-Slide, Cs/1000. **880041** \$116.00
 5-Slide, end opening, Cs/200. USA. **62706** \$87.30
 62706 mailers can also be used as staining jars.
 (1258717B)

Aluminum



Plastic




Slide Trays

- Aluminum, Eberbach, USA. **2870** \$19.60
 Plastic, white, 20-slide, Bel-Art. **445810000** \$12.30
 Wood, 50-slide \$42.00

Slide & Cassette Pens, Markers & Scribes

Write a 0.4 mm stroke on glass slides with this permanent marker.
 Pilot pens are ideal for cytologists and pathologists.



Pens

- Pilot®, permanent, extra fine point ea bx/12
 Black. **44102** \$1.70 17.45
 Blue. **44103** \$1.70 17.45
 Green. **44105** \$1.70 17.45
 Red. **44104** \$1.70 17.45



Confine immunohistochemical reagents within a light blue, hydrophobic barrier, or separate multiple specimens on a slide/dish. Resists alcohol and acetone, but removes with xylene.

- Super^{HT} PAP (Peroxidase Anti-Peroxidase), hydrophobic
 Regular Tip (4 mm tapered nib). **K-600** \$59.00
 (71310, 195505)
 Fine Tip (2.5 mm tapered nib). **K-650** \$44.00
 (71312, 195506)

Label any glossy surface with this water-resistant wax pencil. Wipes off easily with a cloth. Simply pull the string to remove wrapping and expose more tip.



Markers, USA

- China Marker (Wax Pencil). **02089** ea bx/12
 \$1.10 11.50

Ink resists xylene, alcohol, acetone and formalin.
 Remove with acetic acid, chloroform, or ethanol/xylene mixture.



- Securline® Marker II®, permanent, PDC ea bx/12
 Black, 0.6 mm tip. **1450-20-PDC** \$2.50 20.30
 (14-905-30, 72107)

Fine Tip



Ultra-Fine Tip



Sharpie Retractable, permanent

- Fine Tip ea bx/12
 Black. **36701** — 29.65
 Blue. **36703** \$2.85 29.65
 Red. **36702** \$2.85 29.65
 Ultra-Fine Tip ea bx/12
 Black. **1735790** \$2.85 29.65
 Blue. **1735792** \$2.85 29.65
 Red. **1735791** \$2.85 29.65

Continued



Scribes, retractable

Carbide-Tipped (resharpenable). RSC	\$15.30
Carbide Tip (replacement). RCT-1	\$6.20
Diamond-Tipped. RSD	\$21.85
Diamond Tip (replacement). RDT-1	\$12.90

Evenly warm your slides from room temperature to 70°C. Black anodized surface provides excellent contrast.



Slide Warmers


Small, 26 x 18 cm (10 x 7"). XH-2002	\$237.00
Large, 63 x 20 cm (25 x 8"). XH-2001	\$318.00

Staining

Glass staining jar with ground glass lid holds 8 coverslips, 18-22 mm, back-to-back. Jar is 45 mm wide and tall.



Coverslip Staining Jar, each. **1257200011**..... \$7.40



Staining/Moisture Chambers

10-Slide, polystyrene, USA	<u>ea</u>	<u>pk/3</u>
Clear. 240-9000-010	\$50.10	141.00
Black. 240-9020-Z10	\$50.10	141.00




20-Slide, Simport

Clear Lid (PETG), black base (ABS), each. M920-1	\$158.00
Black Lid (ABS), black base (ABS), each. M920-2	\$158.00




672



Staining Dishes

Glass (flint)

20-Slide (back-to-back), holds slides horizontally in grooves	
Set (grooved dish & lid), USA. 672	\$12.30
Bottom only. 672-BTM	\$8.65
Lid only. 672-LID	\$6.00



20-Slide, holds slides horizontally in metal rack	
Approx. inside dimensions: 9.6 x 7.6 x 5.7 cm(H)	
Set (dish, rack & lid). SDS-20	\$25.65
Dish & Lid. SD-20	\$15.65
Rack only, stainless steel. SR-20	\$12.90
30-Slide, holds slides vertically in metal rack	
Approx. inside dimensions: 9.6 x 8.5 x 8.3 cm(H)	
Set (dish, rack & lid). SDS-30	\$35.00
Dish & Lid. SD-30	\$17.50
Rack only, stainless steel. SR-30	\$21.50
60-Slide, holds slides vertically in metal rack	
Approx. inside dimensions: 15.5 x 8.3 x 8.5 cm(H)	
Set (dish, rack & lid). SDS-60	\$49.90
Rack only, stainless steel. SR-60	\$28.15



Plastic rack and dish resist alcohol and xylene; can be microwaved or frozen. Lightproof staining dish has a tight lid with a perforated slide partition for supporting a drying agent. Small size of dish saves on reagents, since only 150 mL is needed (holds 250 mL).

Plastic (polyoxymethylene [POM]), holds 25 slides	
Set (dish, rack & lid). 0990-99	\$50.90