

20-Slide
Cardboard
Folder



10-Slide/
10-Block
Consulting
Folder



Slide Folders

- 10-Slide/10-Block (consulting), each. **CSF-101** \$24.75
- 12-Slide, holds 75 x 50 mm (3 x 2") slides
Cardboard, black tape, USA. **SF12-BLACK** \$7.20
- 20-Slide, holds 75 x 25 mm (3 x 1") slides
Cardboard, USA
Folders made in the USA resist heat-induced warping.
- Black Tape. **SF20-BLACK** \$6.25
- Red Tape. **SF20-RED** \$6.25
- Green Tape. **SF20-GREEN** \$6.25
- Navy Blue Tape. **SF20-BLUE** \$6.25
- Cardboard, Black Tape, China. **HS9902** \$4.70
(12-587-10)
- Plastic, gray tray w/2 clear doors. **M750-20GY** \$9.20
Transparent covers hold visible slides securely in place,
but can fold flat against back of tray when opened.
Slides are slightly elevated for easy retrieval.




Slide Holders, perfect for slide drying or staining

- 12-Slide (horizontally positioned), polystyrene. **445570000** \$7.80
- 52-Slide (approx.), 12 x 4", HDPE, USA. **L500-C** \$27.20
- 91-Slide (approx.), 12 x 7", HDPE, USA. **L500-D** \$36.35



2-Slide



5-Slide

Slide Mailers, Plastic

- 2-Slide, Cs/1000. **880041** \$116.00
- 5-Slide, end opening, Cs/200. USA. **62706** \$87.30
62706 mailers can also be used as staining jars.
(1258717B)

Aluminum



Plastic



Slide Trays

- Aluminum, Eberbach, USA. **2870** \$19.60
- Plastic, white, 20-slide, Bel-Art. **445810000** \$12.30
- Wood, 50-slide \$42.00

Slide & Cassette Pens, Markers & Scribes

Write a 0.4 mm stroke on glass slides with this permanent marker.
Pilot pens are ideal for cytologists and pathologists.



Pens

- | | | |
|-------------------------------------|-----------|--------------|
| Pilot®, permanent, extra fine point | <u>ea</u> | <u>bx/12</u> |
| ★ Black. 44102 | \$1.70 | 17.45 |
| Blue. 44103 | \$1.70 | 17.45 |
| Green. 44105 | \$1.70 | 17.45 |
| Red. 44104 | \$1.70 | 17.45 |



Confine immunohistochemical reagents within a light blue, hydrophobic barrier, or separate multiple specimens on a slide/dish. Resists alcohol and acetone, but removes with xylene.

- Super^{HT} PAP (Peroxidase Anti-Peroxidase), hydrophobic
Regular Tip (4 mm tapered nib). **K-600**..... \$59.00
(71310, 195505)
- ★ Fine Tip (2.5 mm tapered nib). **K-650** \$44.00
(71312, 195506)

Label any glossy surface with this water-resistant wax pencil. Wipes off easily with a cloth. Simply pull the string to remove wrapping and expose more tip.



Markers, USA

- | | | |
|---|-----------|--------------|
| | <u>ea</u> | <u>bx/12</u> |
| China Marker (Wax Pencil). 02089 | \$1.10 | 11.50 |

Ink resists xylene, alcohol, acetone and formalin.
Remove with acetic acid, chloroform, or ethanol/xylene mixture.



- | | | |
|---|-----------|--------------|
| Securline® Marker II®, permanent, PDC | <u>ea</u> | <u>bx/10</u> |
| Black, 0.6 mm tip. 1450-20-PDC | \$2.50 | 20.30 |
| (14-905-30, 72107) | | |

Fine Tip



Ultra-Fine Tip



Sharpie Retractable, permanent

- | | | |
|-----------------------------|-----------|--------------|
| Fine Tip | <u>ea</u> | <u>bx/12</u> |
| Black. 36701 | — | 29.65 |
| Blue. 36703 | \$2.85 | 29.65 |
| Red. 36702 | \$2.85 | 29.65 |
| Ultra-Fine Tip | <u>ea</u> | <u>bx/12</u> |
| Black. 1735790 | \$2.85 | 29.65 |
| Blue. 1735792 | \$2.85 | 29.65 |
| Red. 1735791 | \$2.85 | 29.65 |

Continued

Rectangular. Recommended for ice pans

36-Hole, 1.5 mL tubes, 23.5 x 11 x 1 cm.....	\$34.40
Black/188480112 Green/188480114 Red/188480113	
24-Hole, 14/15 mL tubes, 23.5 x 11 x 1.7 cm.....	\$36.05
Black/188480122 Green/188480124 Red/188480123	




Incubators, myTemp™ Mini Digital, Benchmark

- 20 L capacity (0.8 cu ft)
- Digital temperature control with LED display
- Two models: Heat Only, or Heat/Cool
- Accepts bottles and flasks up to 2 L
- Two adjustable/removable shelves
- 9.3"(D) x 10.3"(W) x 12.8"(H) interior;
- 14.5"(D) x 13.2"(W) x 18.7"(H) exterior
- 2-year warranty

Heat Only, 1° above ambient to 60°C. H2200-H	\$407.00
Heat/Cool, 15° below ambient to 60°C. H2200-HC	\$637.00
Heat/Cool/CO ₂ , 10° below ambient to 55°C. H2300-HC2	\$2,998.00

Label plastics, glass and metals with these chemically inert and temperature resistant (-196 to 150°C) polyolefin tags.



Labels

Cryo-Babies®, white, Diversified Biotech	
24 x 13 mm (0.94 x 0.5"), Roll/1000. CRY-0500	\$29.75
Fits 0.5 mL tubes.	
33 x 13 mm (1.28 x 0.5"), Roll/1000. CRY-1000	\$29.75
Fits 1.5-2.0 mL tubes.	

Markers, USA



Sharpie Retractable, permanent

Fine Tip	ea	bx/12
Black. 36701	\$2.85	29.65
Blue. 36703	\$2.85	29.65
Red. 36702	\$2.85	29.65
Ultra-Fine Tip	ea	bx/12
Black. 1735790	\$2.85	29.65
Blue. 1735792	\$2.85	29.65
Red. 1735791	\$2.85	29.65

Label any glossy surface with this water-resistant wax pencil. Wipes off easily with a cloth. Simply pull the string to remove wrapping and expose more tip.




China Marker (Wax Pencil). 02089	ea	bx/12
	\$1.10	11.50

Mortar and Pestle Sets

Glass, USA

59 mL (2 oz) Capacity. 826-2	\$15.25
118 mL (4 oz) Capacity. 826-3	\$16.00
237 mL (8 oz) Capacity. 826-4	\$19.00



Porcelain, outside wall and rim glazed

30 mL Capacity, 63 mm top OD. LPC-154	\$4.60
80 mL Capacity, 80 mm top OD. LPC-154-1	\$6.15
100 mL Capacity, 105 mm top OD. LPC-155	\$7.60
320 mL Capacity. LPC-156	\$10.30

Petri Dishes



Petri Dishes

PS (polystyrene)

90 x 20 mm, economy	
Pk/10, nonsterile. 51651	\$5.00
35 x 10 mm, vented lid, stacking ring, Simport	
Cs/500, sterile sleeve/20. D210-15	\$130.00
60 X 15 mm, stacking ring, Simport	
Cs/500, sterile sleeve/20. D210-13	\$102.00
100 X 15 mm, Kord-Valmark	
Slippable (rimless lid for automation)	
Standard, Cs/500, sterile sleeve/25. 2917.. ★	\$103.00
ISO Mark, Cs/500, sterile sleeve/25. 2910P	\$103.00
Stackable	
Standard, Cs/500, sterile sleeve/25. 2900	\$104.00
100 X 20 mm, stacking ring	
Cs/500, sterile sleeve/20, Simport. D210-7	\$163.00

Glass (borosilicate), reinforced rims, economy

90 x 15 mm, each. 31015003	\$2.75
100 x 15 mm, each. 31015004	\$3.20



Store (54) 100 mm petri dishes in this clear, 6-column rack. 34(L) x 20(W) x 23 cm(H).

Petri Dish Rack, 100 mm, acrylic, Bel-Art. 189850100	\$116.00
--	----------

Continued