



**Scribes, retractable**

Carbide-Tipped (resharpenable). <b>RSC</b> .....	\$15.30
Carbide Tip (replacement). <b>RCT-1</b> .....	\$6.20
Diamond-Tipped. <b>RSD</b> .....	\$21.85
Diamond Tip (replacement). <b>RDT-1</b> .....	\$12.90

Evenly warm your slides from room temperature to 70°C. Black anodized surface provides excellent contrast.



**Slide Warmers**


Small, 26 x 18 cm (10 x 7"). <b>XH-2002</b> .....	\$237.00
Large, 63 x 20 cm (25 x 8"). <b>XH-2001</b> .....	\$318.00

**Staining**

Glass staining jar with ground glass lid holds 8 coverslips, 18-22 mm, back-to-back. Jar is 45 mm wide and tall.



Coverslip Staining Jar, each. **1257200011**..... \$7.40



**Staining/Moisture Chambers**

10-Slide, polystyrene, USA	<u>ea</u>	<u>pk/3</u>
Clear. <b>240-9000-010</b> .....	\$50.10	141.00
Black. <b>240-9020-Z10</b> .....	\$50.10	141.00




**20-Slide, Simport**

Clear Lid (PETG), black base (ABS), each. <b>M920-1</b> .....	\$158.00
Black Lid (ABS), black base (ABS), each. <b>M920-2</b> .....	\$158.00



672



**Staining Dishes**

Glass (flint)

20-Slide (back-to-back), holds slides horizontally in grooves	
Set (grooved dish & lid), USA. <b>672</b> .....	\$12.30
Bottom only. <b>672-BTM</b> .....	\$8.65
Lid only. <b>672-LID</b> .....	\$6.00



20-Slide, holds slides horizontally in metal rack	
Approx. inside dimensions: 9.6 x 7.6 x 5.7 cm(H)	
Set (dish, rack & lid). <b>SDS-20</b> .....	\$25.65
Dish & Lid. <b>SD-20</b> .....	\$15.65
Rack only, stainless steel. <b>SR-20</b> .....	\$12.90
30-Slide, holds slides vertically in metal rack	
Approx. inside dimensions: 9.6 x 8.5 x 8.3 cm(H)	
Set (dish, rack & lid). <b>SDS-30</b> .....	\$35.00
Dish & Lid. <b>SD-30</b> .....	\$17.50
Rack only, stainless steel. <b>SR-30</b> .....	\$21.50
60-Slide, holds slides vertically in metal rack	
Approx. inside dimensions: 15.5 x 8.3 x 8.5 cm(H)	
Set (dish, rack & lid). <b>SDS-60</b> .....	\$49.90
Rack only, stainless steel. <b>SR-60</b> .....	\$28.15



Plastic rack and dish resist alcohol and xylene; can be microwaved or frozen. Lightproof staining dish has a tight lid with a perforated slide partition for supporting a drying agent. Small size of dish saves on reagents, since only 150 mL is needed (holds 250 mL).

Plastic (polyoxymethylene [POM]), holds 25 slides	
Set (dish, rack & lid). <b>0990-99</b> .....	\$50.90