



Scribes, retractable

Carbide-Tipped (resharpenable). RSC	\$15.30
Carbide Tip (replacement). RCT-1	\$6.20
Diamond-Tipped. RSD	\$21.85
Diamond Tip (replacement). RDT-1	\$12.90

Evenly warm your slides from room temperature to 70°C. Black anodized surface provides excellent contrast.



Slide Warmers


Small, 26 x 18 cm (10 x 7"). XH-2002	\$237.00
Large, 63 x 20 cm (25 x 8"). XH-2001	\$318.00

Staining

Glass staining jar with ground glass lid holds 8 coverslips, 18-22 mm, back-to-back. Jar is 45 mm wide and tall.



Coverslip Staining Jar, each. **1257200011**..... \$7.40



Staining/Moisture Chambers

10-Slide, polystyrene, USA	<u>ea</u>	<u>pk/3</u>
Clear. 240-9000-010	\$50.10	141.00
Black. 240-9020-Z10	\$50.10	141.00




20-Slide, Simport

Clear Lid (PETG), black base (ABS), each. M920-1	\$158.00
Black Lid (ABS), black base (ABS), each. M920-2	\$158.00



672



Staining Dishes

Glass (flint)

20-Slide (back-to-back), holds slides horizontally in grooves	
Set (grooved dish & lid), USA. 672	\$12.30
Bottom only. 672-BTM	\$8.65
Lid only. 672-LID	\$6.00




20-Slide, holds slides horizontally in metal rack	
Approx. inside dimensions: 9.6 x 7.6 x 5.7 cm(H)	
Set (dish, rack & lid). SDS-20	\$25.65
Dish & Lid. SD-20	\$15.65
Rack only, stainless steel. SR-20	\$12.90
30-Slide, holds slides vertically in metal rack	
Approx. inside dimensions: 9.6 x 8.5 x 8.3 cm(H)	
Set (dish, rack & lid). SDS-30	\$35.00
Dish & Lid. SD-30	\$17.50
Rack only, stainless steel. SR-30	\$21.50
60-Slide, holds slides vertically in metal rack	
Approx. inside dimensions: 15.5 x 8.3 x 8.5 cm(H)	
Set (dish, rack & lid). SDS-60	\$49.90
Rack only, stainless steel. SR-60	\$28.15



Plastic rack and dish resist alcohol and xylene; can be microwaved or frozen. Lightproof staining dish has a tight lid with a perforated slide partition for supporting a drying agent. Small size of dish saves on reagents, since only 150 mL is needed (holds 250 mL).

Plastic (polyoxymethylene [POM]), holds 25 slides	
Set (dish, rack & lid). 0990-99	\$50.90



Staining Jars, hold slides (up to 76 x 26 mm) vertically
Coplin, holds 10 slides back-to-back

	<u>ea</u>	<u>cs/12</u>
Glass w/Ground Glass Cover, USA. 673	\$10.75	119.00
Glass Cover, Fits 673 Jar. 673-LID	\$4.20	—
	<u>ea</u>	<u>cs/12</u>
Glass w/Screw-On Phenolic Cap, USA. 674	\$10.75	119.00
Phenolic Cap, Fits 674 Jar, black. 674-LID	\$4.25	—
	<u>ea</u>	<u>cs/12</u>
Polypropylene w/Screw-On Cap. 600P90101	\$4.60	45.20
Plastic jar is unbreakable, microwave proof, stain resistant, and chemically inert. Cap is flat.		
	<u>ea</u>	<u>cs/12</u>
Wide Top, Glass w/Square Lid, USA. 675	\$14.55	160.00
holds 16 slides back-to-back		
Glass Lid, Fits 675 Jar. 675-LID	\$4.45	—



Move slides, without forceps, between jars linked together in a workstation. Swing out the rack's integral lid for easy accessibility to slides. Stable, wide-based jar requires only 80 mL of reagent. Polyoxymethylene (POM) resists alcohol and xylene; can be microwaved or frozen.

Wide Top, Plastic (POM)
12-Slide, EasyDip™, Simport

	<u>ea</u>	<u>cs/6</u>
Jar only, Blue. M900-12B	\$10.30	56.70
Jar only, Green. M900-12G	\$10.30	56.70
Jar only, Pink. M900-12P	\$10.30	56.70
Jar only, White. M900-12W	\$10.30	56.70
Jar only, Yellow. M900-12Y	\$10.30	56.70
Rack/Lid only, gray. M905-12DGY	\$11.45	63.00



Anodized aluminum staining jar holder


(Jars sold separately)
Staining Jar Holder, 6-Place. **M906** \$65.60



25-Slide, light green, microwaveable

Jar & Lid. 1030	\$11.35
Rack only. 1042	\$18.90
Staining Jar Station, 12-Place. 4501011	\$259.00
Includes metal frame, (12) 25-slide dishes w/ lids (1030) & (1) rack (1042).	
Dimensions: 29(L) x 5(W) x 5"(H)	



Tissue (Embedding) & Biopsy Cassettes



Base Molds

Plastic, disposable, Simport


15 x 15 x 5 mm, Cs/1000. M475-2	\$126.00
24 x 24 x 5 mm, Cs/1000. M475-3	\$126.00
★ 30 x 24 x 5 mm, Cs/1000. M475-4	\$126.00
37 x 24 x 5 mm, Cs/1000. M475-5	\$126.00
37 x 24 x 10 mm, Cs/500. M475-10	\$65.45
Designed for Simport Macrosette™.	
Metal, reusable	
15 x 15 x 5 mm, Pk/12. HB-15	\$18.00
24 x 24 x 5 mm, Pk/12. HB-24	\$18.00
30 x 24 x 5 mm, Pk/12. HB-30	\$18.00
37 x 24 x 5 mm, Pk/12. HB-37	\$18.00
66 x 45 x 15 mm, Simport, Pk/10. M474-7	\$156.00
Designed for Simport Jumbosette™.	



Biopsy Bags

Nylon, 30 x 50 mm, Pk/1000. 540478-1	\$235.00
Paper (Tea Bags), 2¼ x 2¼", Pk/500. 56100	\$38.00

M476-1
Made of 60 PPI cellular polyester urethane foam. Use rectangles with cassettes, circles with capsules.



Biopsy Pads

Rectangular	
Black Foam, 30.2 x 25.4 mm, Bg/1000. 339014	\$40.00
Blue Foam, 30.2 x 25.4 x 2 mm, Cs/10,000. M476-1	\$185.00
Circular, Simport	
Blue Foam, 25.4 x 2.7 mm, Cs/10,000. M476-2	\$225.00