

GENERAL LAB SUPPLIES & EQUIPMENT

Alcohol Prep Pads, Sterile, 70% Isopropanol
Large, 4.5 x 9 cm (1.8 x 3.5"), Cs/1000. **44058-404**..... \$40.00



Applicators, 6" long

Cotton Tipped, wooden shaft

Nonsterile

Bg/100. **95-8702**..... \$1.45

Bx/1000, 10 bg/bx. **95-8702**..... \$11.45

Sterile Pk/1, Bx/100. **4404304**..... \$7.80

Sterile Pk/2, Bx/200. **44024-106-2S**..... \$7.90

Plain, wood, Bg/1000. **95-8701**..... \$6.90



Rubberized Cloth

Aprons

Disposable, polyethylene, 46"(L) x 28"(W)

0.8 mil Thick (0.02 mm), Bx/100. **105DA2846**..... \$10.95

1.75 mil Thick (0.044 mm), Bx/100. **59010485**..... \$38.00

Rubberized Cloth, black

36"(L) x 27"(W), each. **31012760**..... \$7.00

42"(L) x 27"(W), each. **31012761**..... \$8.00


★ 42"(L) x 36"(W), each. **31012763**..... \$10.00

46"(L) x 36"(W), each. **31012762**..... \$10.70

Bags



Collect biohazardous waste in this space-saving, autoclavable bag. Printed with a biohazard symbol and sterilization indicator patch. Measures 20.3 x 30.5 cm (8 x 12") and is 1.57 mil (0.04 mm) thick polypropylene.



HS10044H

Bag, Biohazard, Benchtop, Bx/500. **HS10320**..... \$68.10

Biohazard Bag Holder, Benchtop. **HS10044H**..... \$20.56

Epoxy-coated wire holder measures 9.1 x 12.8 x 21.4 cm.



Bags display a biohazard symbol and precautionary procedures. Sterilization indicator patch darkens when autoclaved. Durable, extra heavy-duty bags provide maximum protection against dangerous tears and leaks. Bel-Art

Bags, Biohazard, orange-red, autoclavable, polypropylene

Standard, 1.5 mil (0.038 mm) thick

14 x 19" (36 x 48 cm), Bx/200. **131641419**..... \$117.00

19 x 23" (48 x 58 cm), Bx/200. **131641923**..... \$170.00

25 x 35" (64 x 89 cm), Bx/200. **131642535**..... \$304.00

Bags, Biohazard (cont)

Heavy-Duty, 2 mil (0.05 mm) thick

14 x 19" (36 x 48 cm), Bx/200. **131651419**..... \$178.00

19 x 23" (48 x 58 cm), Bx/200. **131651923**..... \$260.00

25 x 35" (64 x 89 cm), Bx/200. **131652535**..... \$447.00



Bags, Reclosable, polyethylene

2 mil Thick, Clear

3 x 7", Pk/100. **17730-100**..... \$7.20

2 mil Thick, White Block for labeling

3 x 4", Pk/100. **2012-100**..... \$7.40

6 x 9", Pk/100. **1309-100**..... \$9.20

6 x 12", Pk/100. **10850-100**..... \$15.40

4 mil Thick, clear

3 x 4", Pk/100. **494-100**..... \$12.00

Safely transport patient specimens with these OSHA and NCCLS compliant bags. Each resealable, 2 mil thick polyethylene bag has an attached document pouch and imprinted biohazard symbol.



Bag, Specimen Transport, 6 x 9", Bx/1000. **4919**..... \$56.70



Bags, Waste Disposal, HDPE, Cleanware™, Bel-Art

8 x 10", Pk/50. **131741008**..... \$18.90

12 x 16", Pk/50. **131741012**..... \$31.80

Sliding blade cuts excess material 19 mm from the seal.



Bag Sealer, electric, 12" (30.5 cm) Wide. **HS28272**..... \$165.00

Bag Sealer Replacement Part Kit. **HS28280**..... \$30.65

Poly Tubing

4 mil

4" x 1500'. **1143**..... \$110.00

6" x 1500'. **1144**..... \$162.00

6 mil

4" x 1000'. **1658**..... \$132.00

6" x 1000'. **1659**..... \$162.00

Poly Tubing Dispenser. **H-28**..... \$34.00

Continued