



Polypropylene Beakers

Economy



Premium

Beakers

Glass (borosilicate), Griffin low form, single scale, economy

50 mL, each. 7-440050.....	\$2.15
100 mL, each. 7-440100.....	\$2.35
250 mL, each. 7-440250.....	\$3.30
600 mL, each. 61020229600.....	\$4.95
1000 mL, each. 7-441000.....	\$7.30
2000 mL, each. 7-442000.....	\$15.45
4000 mL, each. 7-444000.....	\$33.15

Polypropylene

Economy, Triple-Spout, Tri-Stir®	ea	cs/100
50 mL, 5 mL graduations. 80092.....	—	\$27.95
100 mL, 10 mL graduations. 80093.....	—	\$33.95
250 mL, 20 mL graduations. 80094.....	—	43.95
1000 mL, 50 mL graduations. 80098.....	\$1.90	84.40

Premium, Griffin low form, single scale, Bel-Art, USA

50 mL, 10 mL graduations, each. 262100000.....	\$2.85
100 mL, 25 mL graduations, each. 262110000.....	\$2.70
250 mL, 50 mL graduations, each. 262130000.....	\$4.75
600 mL, 100 mL graduations, each. 262150000.....	\$7.55
1000 mL, 200 mL graduations, each. 262160000.....	\$7.70
2000 mL, 500 mL graduations, each. 262170000.....	\$17.60
4000 mL, 1000 mL graduations, each. 262190000.....	\$36.40

Bench Protectors

Protect your bench surface using this 40 mil absorbent liner backed with 0.5 mil polyethylene film from Bel-Art.



Bench Liners, white

Roll, 50' x 20", Labmat™. 246750000.....	\$59.95
Sheets,	
Bel-Art, 20 x 24", Pk/100. 246750020.....	\$205.00
Ahlstrom, 46 x 57 cm (18 x 22½"), Pk/50. 4002-4657.....	\$78.60

Bottles



Store immersion oil or mounting media in this borosilicate glass bottle with an applicator rod and ground glass cap.

Bottles

Balsam, glass, 45 mL (empty), USA. 616.....	\$17.50
---	---------

Biochemical Oxygen Demand (BOD) Bottle



BOD, 300 mL, clear, unnumbered, each. 310BTL-BOD \$11.55



BOD Bottle Carrier, holds 8 bottles, polypropylene, Bel-Art
Bx/2. 169680000 \$87.45




Deliver controlled drops from standard 0.030" tip opening with these LDPE bottles. Sealer cap snaps on easily after use.

Dispensing, squeeze bottles

LDPE (low-density polyethylene), Yorker cap, Bel-Art

30 mL (1 oz), Bg/12. 116370001.....	\$16.90
60 mL (2 oz), Bg/12. 116370002.....	\$19.65
125 mL (4 oz), Bg/12. 116370004.....	\$32.80
185 mL (6 oz), Bg/12. 116370006.....	\$32.80

A long nozzle, textured sides, and slim profile make this chemically inert, HDPE container ideal for dispensing immersion oil.



HDPE (high-density polyethylene)

20 mL, each. 18145.....	\$6.90
-------------------------	--------

Amber glass protects light-sensitive liquids. Screw-cap dropper features a black rubber bulb and glass pipet.



Dropping

Amber Glass (flint), Boston Round

30 mL (1 oz), Bx/12. 12574053.....	\$14.00
Bottle Dropper Assembly (replacement)	
20/400, 7 x 76 mm Pipet, Pk/12. 12574053D.....	\$6.75
60 mL (2 oz), Bx/12. 12574054.....	\$22.00
Bottle Dropper Assembly (replacement)	
20/400, 7 x 89 mm Pipet, Pk/12. 12574054D.....	\$8.65

Continued