



Benchtop, Digital, 0-999 range
8-Key (9-unit). **08KP** \$244.00



Colony, handheld electronic. **HS120000** \$238.00



Handheld Tally, 0-9999 range. **600A**..... \$8.50



Crucibles, Tall Form with Lid, porcelain
5 mL Capacity, Pk/6. **600JCT005** \$20.65
10 mL Capacity, Pk/6. **600JCT010** \$23.60
15 mL Capacity, Pk/6. **600JCT015** \$26.60
30 mL Capacity, Pk/6. **600JCT030** \$33.25
50 mL Capacity, Pk/6. **600JCT050** \$38.90
100 mL Capacity, Pk/6. **600JCT100** \$53.20

Cylinders, graduated

Glass (borosilicate), single scale, Class B, economy
10 mL, 0.2 mL graduations, each. **7-711LG** \$2.65
25 mL, 0.5 mL graduations, each. **7-712LG** \$3.40
50 mL, 1.0 mL graduations, each. **7-713LG** \$3.90
100 mL, 1.0 mL graduations, each. **7-714LG** \$4.85
250 mL, 2.0 mL graduations, each. **310CSS250** \$7.80
500 mL, 5.0 mL graduations, each. **310CSS500** \$11.90
1000 mL, 10.0 mL graduations, each. **310CSS1000** \$18.60

Polypropylene, TC/TD double scale, Bel-Art, USA
10 mL, 0.1 mL graduations, each. **284520000** \$4.35
50 mL, 1.0 mL graduations, each. **284540000** \$6.80
100 mL, 1.0 mL graduations, each. **284550000** \$8.40
250 mL, 2.0 mL graduations, each. **284560000** \$9.40
500 mL, 5.0 mL graduations, each. **284570000** \$12.85
1000 mL, 10.0 mL graduations, each. **284580000** \$17.95
2000 mL, 20.0 mL graduations, each. **284590000** \$32.55



Non-Vacuum



Vacuum

Desiccators

Non-Vacuum, glass w/porcelain sample plate, economy
Glazed porcelain plate has 3.5-4.5 cm center hole.
21 cm (8.25") ID. **1257200001** \$41.50
24 cm (10") ID. **1257200002** \$65.75
Vacuum, polycarbonate top/polypropylene bottom
138 mm Plate Diameter. **60055204** \$52.00
190 mm Plate Diameter. **60055205** \$70.65
238 mm Plate Diameter. **60055206** \$89.50

Bel-Art
Cabinets



Clear



Bronze

Cabinet, acrylic

Clear, 30.5 x 30.5 x 30.5 cm. **420660000** \$508.00
Bronze, 30.5 x 30.5 x 30.5 cm. **420660001** \$575.00



367850016

Dippers, 500 mL polypropylene cup with HDPE handle, Bel-Art

3' Handle. **367840016** \$49.55
6' Handle. **367850016** \$64.55
Replacement Cup, screw-on. **367860016** \$11.80

Dish, specimen, glass, stackable, USA. 600 \$9.70
120 mm(OD) x 50 mm



188050000



HS23243A



188190023

Drying Racks

Non-Electric
Bel-Art. **188050000** \$98.35
Dry labware on this epoxy-coated steel wire rack with polystyrene tray. Holds up to 26 pieces on 12 straight pins and 14 hairpin loops. Dimensions: 17 x 43 x 48 cm(H).
Heathrow. **HS23243A** \$109.00
Electric, double-sided, 3 tier, Bel-Art. **188190023** \$917.00
1-year warranty
Dust Filter & Holder (replacement). **999321563** \$12.00
Heating Element (replacement). **999321564** \$193.00