



Eyewash Stations, Bel-Art

- Single
 - 16 oz Bottle (Empty). 248650000 \$41.45
 - 32 oz Bottle (Empty). 248660000 \$44.95
- Double, Two 32 oz Bottles (Empty). 248680000 \$54.75
- Eyewash Solution, sterile saline, 32 oz Bottle. 248790032 \$24.20

Eyewear



Glasses, Safety, integral side shield, absorbs 99.9% UV. 74002... \$5.60



Flexible frame and adjustable elastic strap provides a tight seal.

Goggles, Safety, vented (chemical splash). 97-2001 \$3.55

Film & Foil



- Parafilm**, thermoplastic sealing film, USA
- 2" x 250' (5 cm x 76 m) Roll. HS234526A \$26.10
 - 4" x 125' (10 cm x 38 m) Roll. HS234526B \$26.65
 - 4" x 250' (10 cm x 76 m) Roll. HS234526C \$44.65

Holds one or two 2" rolls, or one 4" roll.



Parafilm Dispenser, ABS plastic, blue. HS234525B \$38.20

Continued



Super-Polyfoil™, polyester-backed aluminum foil, Bel-Art

- 4 x 4" (10.2 x 10.2 cm) Sheets, Pk/100. 132000044 \$32.90
- 5½ x 5½" (13 x 13 cm) Sheets, Pk/100. 132000066 \$42.70

Filter & Blotting Papers



Ahlstrom -Munksjo
Filter Paper

Qualitative, Standard Grades (for gravity filtration)

Slow flow (4 mL/min), 1.5 µm retention, **Grade 610**

- (Whatman 5; S&S 602)
- 9 cm, Pk/100. 6100-0900 \$9.40
- 12.5 cm, Pk/100. 6100-1250 \$11.90
- 15 cm, Pk/100. 6100-1500 \$14.30
- 20.5 cm, Pk/100. 6100-2050 \$22.60
- 24 cm, Pk/100. 6100-2400 \$32.10

Medium flow (30 mL/min), 2 µm retention, **Grade 642**

- (Whatman 2, 2CHR, 20CHR; S&S 2043B, 2045A/B)
- 9 cm, Pk/100. 6420-0900 \$7.35
- 15 cm, Pk/100. 6420-1500 \$13.30

Medium flow (55 mL/min), 2.5 µm retention, **Grade 601**

- (Whatman 1, 1CHR; S&S 593, 597, 2043A)
- 4.25 cm, Pk/100. 6010-0425 \$5.65
- 9 cm, Pk/100. 6010-0900 \$7.55
- 12.5 cm, Pk/100. 6010-1250 \$10.95
- 15 cm, Pk/100. 6010-1500 \$13.35
- 18.5 cm, Pk/100. 6010-1850 \$19.35
- 20.5 cm, Pk/100. 6010-2050 \$20.85
- 33 cm, Pk/100. 6010-3300 \$56.35

Medium flow (60 mL/min), 6 µm retention, **Grade 613**

- (Whatman Student; S&S 595, 0903)
- 9 cm, Pk/100. 6130-0900 \$6.45
- 15 cm, Pk/100. 6130-1500 \$11.30
- 20.5 cm, Pk/100. 6130-2050 \$17.00
- 33 cm, Pk/100. 6130-3300 \$37.75

Medium flow (200 mL/min), 10 µm retention, **Grade 631**

- (Whatman 4, 4CHR; S&S 604, 2040A/B; Fisher P8)
- 11 cm, Pk/100. 6310-1100 \$9.15
- 20.5 cm, Pk/100. 6310-2050 \$21.05
- 33 cm, Pk/100. 6310-3300 \$56.90

Fast flow (235 mL/min), 25 µm retention, creped, **Grade 615**

- 9 cm, Pk/100. 6150-0900 \$6.65
- 12.5 cm, Pk/100. 6150-1250 \$9.70
- 15 cm, Pk/100. 6150-1500 \$11.65
- 20.5 cm, Pk/100. 6150-2050 \$17.50
- 33 cm, Pk/100. 6150-3300 \$38.85



Continued