

Very Fast flow (360 mL/min), 35 µm retention, creped, <b>Grade 617</b>	
9 cm, Pk/50. <b>6170-0900</b> .....	\$6.05
12.5 cm, Pk/50. <b>6170-1250</b> .....	\$8.15
15 cm, Pk/50. <b>6170-1500</b> .....	\$9.50
20.5 cm, Pk/50. <b>6170-2050</b> .....	\$14.45
33 cm, Pk/50. <b>6170-3300</b> .....	\$30.95

**Qualitative, Wet-Strengthened Grades** (for vacuum filtration)

Medium flow (65 mL/min), 3 µm retention, <b>Grade 909</b> ( <i>Whatman 91; S&amp;S 0860, 1574, 3003; Reeve-Angel 201</i> )	
32 cm, Pk/100. <b>9090-3200</b> .....	\$61.40
Fast flow (245 mL/min), 27 µm retention, creped, <b>Grade 962</b> ( <i>Whatman 114; S&amp;S 0856, 0905, 3002; Reeve-Angel 202</i> )	
11 cm, Pk/100. <b>9620-1100</b> .....	\$9.15

**Qualitative, Pleated (Fluted)**

Slow flow (30 mL/min), 2 µm retention, <b>Grade 542</b> ( <i>Whatman 2V</i> )	
12.5 cm, Pk/100. <b>5420-1250</b> .....	\$65.15
18.5 cm, Pk/100. <b>5420-1850</b> .....	\$73.75
Medium flow (60 mL/min), 6 µm retention, <b>Grade 513</b>	
32 cm, Pk/100. <b>5130-3200</b> .....	\$86.65

**Quantitative, Ashless**

(Designed for use in analytical and gravimetric analysis)

Slow flow (4 mL/min), 1.5 µm retention, <b>Grade 94</b> ( <i>Whatman 42; S&amp;S 589/3</i> )	
9 cm, Pk/100. <b>0940-0900</b> .....	\$26.35
15 cm, Pk/100. <b>0940-1500</b> .....	\$45.35
Medium flow (30 mL/min), 2 µm retention, <b>Grade 74</b> ( <i>Whatman 40; S&amp;S 589/2, 589/4</i> )	
9 cm, Pk/100. <b>0740-0900</b> .....	\$26.35
15 cm, Pk/100. <b>0740-1500</b> .....	\$45.35
Fast flow (200 mL/min), 10 µm retention, <b>Grade 54</b> ( <i>Whatman 41; S&amp;S 589/1</i> )	
9 cm, Pk/100. <b>0540-0900</b> .....	\$26.35
12.5 cm, Pk/100. <b>0540-1250</b> .....	\$36.05
15 cm, Pk/100. <b>0540-1500</b> .....	\$45.35
Fast flow (220 mL/min), 15 µm retention, <b>Grade 55</b> ( <i>Whatman 54, 541; S&amp;S 1505</i> )	
9 cm, Pk/100. <b>0550-0900</b> .....	\$29.75
15 cm, Pk/100. <b>0550-1500</b> .....	\$47.80

**Glass Microfiber**

Medium flow (102 mL/min), 1.1 µm retention, <b>Grade 161</b> ( <i>Whatman 934-AH; S&amp;S GF30,GF50; Reeve-Angel 934-AH</i> )	
2.4 cm, Pk/100. <b>1610-0240</b> .....	\$21.95
4.7 cm, Pk/100. <b>1610-0470</b> .....	\$31.10
5.5 cm, Pk/100. <b>1610-0550</b> .....	\$35.90
7 cm, Pk/100. <b>1610-0700</b> .....	\$44.10
9 cm, Pk/100. <b>1610-0900</b> .....	\$57.20
11 cm, Pk/100. <b>1610-1100</b> .....	\$76.15

**Blotting/Chromatography/Electrophoresis Papers**

0.19 mm thick, 125 mm/30 min absorption, <b>Grade 601</b> ( <i>Whatman 1, 1CHR; S&amp;S 593, 597, 2043A</i> )	
46 x 57 cm, Pk/100. <b>6018-4657</b> .....	\$107.00
0.34 mm thick, 140 mm/30 min absorption, <b>Grade 238</b> ( <i>Whatman 3MM, S&amp;S 598,2316</i> )	
7 x 10 cm, Pk/100. <b>2388-0710</b> .....	\$25.30
20 x 20 cm, Pk/100. <b>2388-2020</b> .....	\$55.60
46 x 57 cm, Pk/100. <b>2388-4657</b> .....	\$203.00
0.42 mm thick, 150 mm/30 min absorption, <b>Grade 237</b> ( <i>Whatman 3, 3CHR; S&amp;S 22/2, 591, 2469, 5703</i> )	
46 x 57 cm, Pk/100. <b>2378-4657</b> .....	\$233.00
0.83 mm thick, 220 mm/30 min absorption, <b>Grade 222</b> ( <i>Whatman 31ETCHR; S&amp;S 470, 740E</i> )	
7 x 8 cm, Pk/100. <b>2228-0708</b> .....	\$17.75
7 x 10 cm, Pk/100. <b>2228-0710</b> .....	\$19.75
10 x 15 cm, Pk/100. <b>2228-1015</b> .....	\$36.35
20 x 20 cm, Pk/100. <b>2228-2020</b> .....	\$56.35
0.95 mm thick, 180 mm/30 min absorption, <b>Grade 243</b> ( <i>Whatman 17CHR; S&amp;S 2668</i> )	
11 x 14 cm, Pk/100. <b>2438-1114</b> .....	\$24.05
20 x 20 cm, Pk/100. <b>2438-2020</b> .....	\$44.60
46 x 57 cm, Pk/100. <b>2438-4657</b> .....	\$342.00
2.48 mm thick, 240 mm/30 min absorption, <b>Grade 320</b> ( <i>S&amp;S 2727</i> )	
11 x 14 cm, Pk/50. <b>3208-1114</b> .....	\$38.00
16 x 18 cm, Pk/50. <b>3208-1618</b> .....	\$52.95
20 x 20 cm, Pk/50. <b>3208-2020</b> .....	\$64.60



Forceps, Filter, 11.5 cm. **300-025** ..... \$4.35



<b>Flask Weights, vinyl-coated lead</b>	
C-Shaped lead rings, blue, Heathrow	
Fits 125-500 mL Flasks, 42 mm ID, 227 g. <b>HS8881A</b> .....	\$23.75
Fits 250-1000 mL Flasks, 51 mm ID, 472 g. <b>HS8881B</b> .....	\$27.05
Fits 500-2000 mL Flasks, 54 mm ID, 599 g. <b>HS8881C</b> .....	\$29.80
O-Shaped lead ring, blue, Heathrow	
Fits 1000-4000 mL Flasks, 70 mm ID, 970 g. <b>HS8882D</b> .....	\$29.80
Lead Wire, 300 mm long, 114 g. <b>183240000</b> .....	\$30.10