

## Freezer Racks



175HS-9-3



175HD-442



175V-13-2L

Freezer Racks, stainless steel, USA

### Horizontal (fits upright freezers)

Horizontal, Holds 2" Boxes	1-9	≥10
9-Place (3L x 3H). <b>175HS-9-2</b> .....	\$54.00	48.00
Dimensions: 16.31(L) x 5.5(W) x 6.5"(H)		
★ 12-Place (3L x 4H). <b>175HS-12-2</b> .....	\$72.00	64.00
Dimensions: 16.31(L) x 5.5(W) x 8.75"(H)		
15-Place (3L x 5H). <b>175HS-15-2</b> .....	\$90.00	80.00
Dimensions: 16.31(L) x 5.5(W) x 11"(H)		
16-Place (4L x 4H). <b>175HL-16-2</b> .....	\$96.00	85.00
Dimensions: 22(L) x 5.5(W) x 8.75"(H)		
20-Place (4L x 5H). <b>175HL-20-2</b> .....	\$122.00	108.00
Dimensions: 22(L) x 5.5(W) x 11"(H)		
Horizontal, Holds 3" Boxes, fits upright freezers		
9-Place (3L x 3H). <b>175HS-9-3</b> .....	\$56.00	50.00
Dimensions: 16.31(L) x 5.5(W) x 9.5"(H)		
12-Place (4L x 3H). <b>175HL-12-3</b> .....	\$74.00	66.00
Dimensions: 22(L) x 5.5(W) x 9.5"(H)		

### Horizontal Drawer (fits upright freezers)

Horizontal Drawer, stainless steel drawers		
Holds 2" Boxes		
9-Place (3L x 3H). <b>175HD-332</b> .....	\$79.00	70.00
Dimensions: 16.5(L) x 5.5(W) x 7.125"(H)		
12-Place (3L x 4H). <b>175HD-342</b> .....	\$103.00	92.00
Dimensions: 16.5(L) x 5.5(W) x 9.5"(H)		
15-Place (3L x 5H). <b>175HD-352</b> .....	\$132.00	117.00
Dimensions: 16.5(L) x 5.5(W) x 11.812"(H)		
16-Place (4L x 4H). <b>175HD-442</b> .....	\$141.00	125.00
Dimensions: 22(L) x 5.5(W) x 9.5"(H)		
20-Place (5L x 4H). <b>175HD-542</b> .....	\$175.00	156.00
Dimensions: 27.25(L) x 5.5(W) x 9.5"(H)		
Holds 5" Boxes		
6-Place (3L x 2H). <b>175RD-325</b> .....	\$69.00	61.00
Dimensions: 18.125(L) x 6(W) x 10.75"(H)		
8-Place (4L x 2H). <b>175HD-425</b> .....	\$92.00	82.00
Dimensions: 24(L) x 6(W) x 10.75"(H)		
Holds 15 mL tubes		
160-Place (80L x 2H). <b>175HD-32ML-15</b> .....	\$93.00	82.00
Dimensions: 16.5(L) x 5.5(W) x 9.5"(H)		
220-Place (110L x 2H). <b>175HD-42ML-15</b> .....	\$120.00	107.00
Dimensions: 22(L) x 5.5(W) x 9.5"(H)		
Holds 50 mL tubes		
84-Place (42L x 2H). <b>175HD-42ML-50</b> .....	\$120.00	107.00
Dimensions: 22(L) x 5.5(W) x 9.5"(H)		

## Horizontal Drawer-Tray (fits upright freezers)

Drawer-Tray, aluminum drawer-trays

Holds 2" Boxes		
12-Place (3L x 4H). <b>175HDT-342</b> .....	\$89.00	79.00
Dimensions: 16.63(L) x 5.5(W) x 8.5"(H)		
18-Place (3L x 6H). <b>175HDT-362</b> .....	\$133.00	118.00
Dimensions: 16.63(L) x 5.5(W) x 12.75"(H)		

### Vertical (fits chest freezers)

Vertical, Holds 2" Boxes		
13-Place. <b>175V-13-2L</b> .....	\$81.00	72.00
Dimensions: 5.5(L) x 5.625(W) x 28.562"(H)		

## Funnels



### Funnels

Disposable, polypropylene, short stems, Simport		
55 mm Top ID, takes 11 cm filter, Cs/100. <b>F490-3</b> .....	\$31.15	
65 mm Top ID, takes 12.5 cm filter, Cs/100. <b>F490-4</b> .....	\$31.15	
Glass (borosilicate), long stem (150 mm)		
50 mm Top ID, each. <b>600GF614050</b> .....	\$5.95	
65 mm Top ID, each. <b>600GF614065</b> .....	\$7.15	
75 mm Top ID, each. <b>600GF614075</b> .....	\$7.35	
Heavy-Duty, polypropylene, Bel-Art		
152 mm Top ID, 480 mL capacity, Pk/3. <b>146970000</b> .....	\$37.85	
184 mm Top ID, 720 mL capacity, Pk/3. <b>146980000</b> .....	\$42.15	
267 mm Top ID, 3845 mL capacity, each. <b>146990000</b> .....	\$33.60	
Heavy-Duty, Molded-In Handle, polypropylene		
250 mm Top ID, 2485 mL capacity, each. <b>60057232</b> .....	\$19.70	
High-Speed, TPX (polymethylpentene), Urbanti®, Bel-Art		
Internal helicoid ribs increase filter speed, and vertical outer ribs prevent air lock.		
100 mm Top ID, 224 mL capacity, Pk/4. <b>146460000</b> .....	\$24.80	
140 mm Top ID, 586 mL capacity, Pk/3. <b>146480000</b> .....	\$69.20	
196 mm Top ID, 1532 mL capacity, Pk/2. <b>146500000</b> .....	\$60.70	
Micro, polypropylene, Bel-Art		
24 mm Top ID, 2.6 mL capacity, Pk/12. <b>146850024</b> .....	\$29.70	
28 mm Top ID, 5.2 mL capacity, Pk/12. <b>146850028</b> .....	\$24.30	
35 mm Top ID, 8.1 mL capacity, Pk/12. <b>146850035</b> .....	\$22.40	
43 mm Top ID, 15.7 mL capacity, Pk/12. <b>146850043</b> .....	\$51.80	

Continued