

— Hot Plate, Magnetic Stirrer & Hot Plate Stirrers —



Hot Plates

- 6" Diameter Cast Iron Burner, 1000W. 1257225 \$28.05
- 7.5 x 7.5" Ceramic Platform, 80-380°C. H4000-H \$208.00




- Ceramic platform is 7.5 x 7.5" (19 x 19 cm)
- Stir speed is adjustable from 60-1500 rpm
- 2-year warranty

Magnetic Stirrer, Benchmark. H4000-S \$208.00

Hot Plate Stirrer, 50-380°C, Benchmark. H4000-HS \$313.00

— Homogenizers & Tissue Grinder —



Rapidly homogenize, emulsify, suspend and disrupt 0.1 mL to 250 mL samples. Includes a 5 x 50 mm flat-bottom probe & 7 x 50 mm saw tooth probe for use with micro tubes.

Homogenizers

Handheld. D1000 \$919.00



Homogenize a variety of sample types, from soft tissue to bone, in 35 seconds or less.


Benchtop, Microtube, BeadBlaster™ 24. D2400 \$6,796.00

- Holds (24) 2.0 mL tubes, or (6) 5 mL tubes w/optional tube holder
- Speed adjusts from 4.00 to 7.00 m/s in 0.05 m/s increments
- Stores 50 programs in memory, each with a maximum of 10 cycles
- 2-year warranty

Homogenizer Prefilled Tubes

2.0 mL screw-capped tubes are prefilled with your choice of the following beads:

- Glass, acid washed
 - 0.1 mm, Pk/50. D1031-01 \$97.00
 - 0.5 mm, Pk/50. D1031-05 \$97.00
- Zirconium Beads, acid washed/heat treated/nuclease-free
 - 1.5 mm, Pk/50. D1032-15 \$136.00
 - 3.0 mm, Pk/50. D1032-30 \$136.00
- Stainless Steel Beads, acid washed
 - 2.8 mm, Pk/50. D1033-28 \$136.00



Tissue Grinder, Dual® Glass Pestle

- 1 mL Capacity, LG10640100 \$78.50
- 3 mL Capacity, LG10640102 \$88.35
- 5 mL Capacity, LG10640104 \$100.00



Hydrometers Access www.belart.com and call us to place an order




Temperature Range: -196 to 100°C

Ice Buckets with lids, expanded PVC foam, Bel-Art, USA

- 2.5 L Capacity (Handy)
 - Black. 188482002 \$77.35
 - Blue. 188482001 \$77.35
 - Purple. 188482006 \$77.35
- 4 L Capacity (N'Icer)
 - Black. 188484002 \$94.40
 - Blue. 188484001 \$94.40

Ice Pans w/o lids, expanded PVC foam, Bel-Art, USA

- 1 L Capacity (Mini), 15.5 x 15.5 x 8.5 cm(H)
 - Black. 188481102 \$55.70
 - Blue. 188481101 \$55.70
 - Purple. 188481106 \$55.70
- 4 L Capacity (Midi), 33 x 24.5 x 12.5 cm(H)
 - Black. 188484102 \$91.85
 - Blue. 188484101 \$91.85
 - Purple. 188484106 \$91.85
- 9 L Capacity (Maxi), 44 x 35 x 13 cm(H)
 - Black. 188489102 \$110.00
 - Blue. 188489101 \$110.00
 - Purple. 188489106 \$110.00



Tube Racks, Floating, PVC Foam, Bel-Art, USA

- Circular. Recommended for ice buckets
 - 44-Hole, 1.5 mL tubes, 22 cm dia. x 1 cm \$34.40
 - Black/188480012 Green/188480014 Red/188480013
 - 30-Hole, 14/15 mL tubes, 22 cm dia. x 1.7 cm \$36.05
 - Black/188480022 Green/188480024 Red/188480023

Continued