

— Hot Plate, Magnetic Stirrer & Hot Plate Stirrers —




Hot Plates
 6" Diameter Cast Iron Burner, 1000W. 1257225 \$28.05
 7.5 x 7.5" Ceramic Platform, 80-380°C. H4000-H..... \$208.00



- Ceramic platform is 7.5 x 7.5" (19 x 19 cm)
- Stir speed is adjustable from 60-1500 rpm
- 2-year warranty

Magnetic Stirrer, Benchmark. H4000-S..... \$208.00
 Hot Plate Stirrer, 50-380°C, Benchmark. H4000-HS..... \$313.00

— Homogenizers & Tissue Grinder —



Rapidly homogenize, emulsify, suspend and disrupt 0.1 mL to 250 mL samples. Includes a 5 x 50 mm flat-bottom probe & 7 x 50 mm saw tooth probe for use with micro tubes.

Homogenizers
 Handheld. D1000 \$919.00



Homogenize a variety of sample types, from soft tissue to bone, in 35 seconds or less.

Benchtop, Microtube, BeadBlaster™24. D2400..... \$6,796.00


- Holds (24) 2.0 mL tubes, or (6) 5 mL tubes w/optional tube holder
- Speed adjusts from 4.00 to 7.00 m/s in 0.05 m/s increments
- Stores 50 programs in memory, each with a maximum of 10 cycles
- 2-year warranty

Homogenizer Prefilled Tubes
 2.0 mL screw-capped tubes are prefilled with your choice of the following beads:

Glass, acid washed
 0.1 mm, Pk/50. D1031-01.....\$97.00
 0.5 mm, Pk/50. D1031-05.....\$97.00

Zirconium Beads, acid washed/heat treated/nuclease-free
 1.5 mm, Pk/50. D1032-15.....\$136.00
 3.0 mm, Pk/50. D1032-30.....\$136.00

Stainless Steel Beads, acid washed
 2.8 mm, Pk/50. D1033-28.....\$136.00



Tissue Grinder, Dual® Glass Pestle
 1 mL Capacity, LG10640100 \$78.50
 3 mL Capacity, LG10640102 \$88.35
 5 mL Capacity, LG10640104 \$100.00



Hydrometers Access www.belart.com
 and call us to place an order



Temperature Range: -196 to 100°C

Ice Buckets with lids, expanded PVC foam, Bel-Art, USA
 2.5 L Capacity (Handy)
 Black. 188482002..... \$77.35
 Blue. 188482001 \$77.35
 Purple. 188482006 \$77.35
 4 L Capacity (N'icer)
 Black. 188484002..... \$94.40
 Blue. 188484001 \$94.40

Ice Pans w/o lids, expanded PVC foam, Bel-Art, USA
 1 L Capacity (Mini), 15.5 x 15.5 x 8.5 cm(H)
 Black. 188481102..... \$55.70
 Blue. 188481101 \$55.70
 Purple. 188481106 \$55.70
 4 L Capacity (Midi), 33 x 24.5 x 12.5 cm(H)
 Black. 188484102..... \$91.85
 Blue. 188484101 \$91.85
 Purple. 188484106 \$91.85
 9 L Capacity (Maxi), 44 x 35 x 13 cm(H)
 Black. 188489102..... \$110.00
 Blue. 188489101 \$110.00
 Purple. 188489106 \$110.00



Tube Racks, Floating, PVC Foam, Bel-Art, USA
 Circular. Recommended for ice buckets
 44-Hole, 1.5 mL tubes, 22 cm dia. x 1 cm \$34.40
 Black/188480012 Green/188480014 Red/188480013
 30-Hole, 14/15 mL tubes, 22 cm dia. x 1.7 cm \$36.05
 Black/188480022 Green/188480024 Red/188480023

Continued

Rectangular. *Recommended for ice pans*

36-Hole, 1.5 mL tubes, 23.5 x 11 x 1 cm.....	\$34.40
Black/188480112 Green/188480114 Red/188480113	
24-Hole, 14/15 mL tubes, 23.5 x 11 x 1.7 cm.....	\$36.05
Black/188480122 Green/188480124 Red/188480123	




Incubators, myTemp™ Mini Digital, Benchmark

- 20 L capacity (0.8 cu ft)
- Digital temperature control with LED display
- Two models: Heat Only, or Heat/Cool
- Accepts bottles and flasks up to 2 L
- Two adjustable/removable shelves
- 9.3"(D) x 10.3"(W) x 12.8"(H) interior;
- 14.5"(D) x 13.2"(W) x 18.7"(H) exterior
- 2-year warranty

Heat Only, 1° above ambient to 60°C. H2200-H	\$407.00
Heat/Cool, 15° below ambient to 60°C. H2200-HC	\$637.00
Heat/Cool/CO ₂ , 10° below ambient to 55°C. H2300-HC2	\$2,998.00

Label plastics, glass and metals with these chemically inert and temperature resistant (-196 to 150°C) polyolefin tags.



Labels

Cryo-Babies®, white, Diversified Biotech	
24 x 13 mm (0.94 x 0.5"), Roll/1000. CRY-0500	\$29.75
Fits 0.5 mL tubes.	
33 x 13 mm (1.28 x 0.5"), Roll/1000. CRY-1000	\$29.75
Fits 1.5-2.0 mL tubes.	

Markers, USA



Sharpie Retractable, permanent

Fine Tip	<u>ea</u>	<u>bx/12</u>
Black. 36701	\$2.85	29.65
Blue. 36703	\$2.85	29.65
Red. 36702	\$2.85	29.65
Ultra-Fine Tip	<u>ea</u>	<u>bx/12</u>
Black. 1735790	\$2.85	29.65
Blue. 1735792	\$2.85	29.65
Red. 1735791	\$2.85	29.65

Label any glossy surface with this water-resistant wax pencil. Wipes off easily with a cloth. Simply pull the string to remove wrapping and expose more tip.



China Marker (Wax Pencil). 02089	<u>ea</u>	<u>bx/12</u>
	\$1.10	11.50

Mortar and Pestle Sets

Glass, USA

59 mL (2 oz) Capacity. 826-2	\$15.25
118 mL (4 oz) Capacity. 826-3	\$16.00
237 mL (8 oz) Capacity. 826-4	\$19.00



Porcelain, outside wall and rim glazed

30 mL Capacity, 63 mm top OD. LPC-154	\$4.60
80 mL Capacity, 80 mm top OD. LPC-154-1	\$6.15
100 mL Capacity, 105 mm top OD. LPC-155	\$7.60
320 mL Capacity. LPC-156	\$10.30

Petri Dishes



Petri Dishes

PS (polystyrene)

90 x 20 mm, economy	
Pk/10, nonsterile. 51651	\$5.00
35 x 10 mm, vented lid, stacking ring, Simport	
Cs/500, sterile sleeve/20. D210-15	\$130.00
60 X 15 mm, stacking ring, Simport	
Cs/500, sterile sleeve/20. D210-13	\$102.00
100 X 15 mm, Kord-Valmark	
Slippable (rimless lid for automation)	
Standard, Cs/500, sterile sleeve/25. 2917.. ★	\$103.00
ISO Mark, Cs/500, sterile sleeve/25. 2910P	\$103.00
Stackable	
Standard, Cs/500, sterile sleeve/25. 2900	\$104.00
100 X 20 mm, stacking ring	
Cs/500, sterile sleeve/20, Simport. D210-7	\$163.00

Glass (borosilicate), reinforced rims, economy

90 x 15 mm, each. 31015003	\$2.75
100 x 15 mm, each. 31015004	\$3.20



Store (54) 100 mm petri dishes in this clear, 6-column rack. 34(L) x 20(W) x 23 cm(H).

Petri Dish Rack, 100 mm, acrylic, Bel-Art. 189850100	\$116.00
--	----------

Continued




Tube Racks, Floating, HDPE Foam w/Handle, Heathrow

- 18-Hole, 0.2/0.5/1.5-2.0 mL tubes, blue circle, Pk/5. **HS2166** \$14.80
- 24-Hole, 1.5-2.0 mL tubes, yellow rectangle, Pk/5. **HS2165C** \$14.55
- 8-Hole, 15 mL tubes, blue parallelogram, Pk/5. **HS2165B** \$12.80
- 4-Hole, 50 mL tubes, green diamond, Pk/5. **HS2165A** \$14.30




Tube Racks, Grip, polypropylene, No-Wire™, Bel-Art

- 10-13 mm Tubes, 90-place, each. **187490000** \$24.40
- 15-16 mm Tubes, 60-place, each. **187490001** \$24.40
- 18-20 mm Tubes, 40-place, each. **187490002** \$24.40



Tube Racks, Interlocking, polypropylene, Heathrow

- Cube, mixed colors, Pk/5. **HS29050A** \$39.80
Holds (4) 50 mL tubes, (10) 15 mL tubes, (12) 12 x 75 or 12 x 100 tubes, or (16) 1.5-2.0 mL tubes.
- 4-Way, mixed colors, Pk/5. **HS29022G** \$36.85
Holds (4) 50 mL tubes, (12) 15 mL tubes, (32) 1.5 mL tubes, or (32) 0.5 mL tubes.
- 4-Way Rotatable, Rota-Rack™, four connected modules
Large Tubes, 22.4 x 11 x 11 cm, each. **HS29040A** \$38.85
Each module holds (2) 50 mL tubes, (8) 15 mL tubes, (6) 20 mm tubes, or (10) 10-13 mm tubes. Mixed colors.
- Microtubes, 21 x 10 x 10 cm, each. **HS29040B** \$38.85
Each module holds (6) 15 mL tubes, (9) 1.5/2.0 mL tubes, (12) 0.5/0.6 mL tubes, or (32) 0.2 mL PCR tubes or (4) 8-tube strips. Mixed colors.




Tube Racks, Mega Rack, polypropylene, Heathrow
10-13 mm Tubes

- 216-Place, 33.7 x 26.7 x 9.1 cm(H)
Blue, each. **HS24322A** \$35.55
White, each. **HS24322E** \$35.55
- 432-Place, 67.4 x 26.7 x 9.1 cm(H)
Blue, each. **HS120113** \$64.90
White, each. **HS120112** \$64.90



Tube Racks, Microtube, polypropylene, Heathrow

- 80-Well (16 x 5), molded-in coordinates, 22.5 x 6.7 x 2.8 cm
Holds 1.5-2.0 mL tubes in conical wells.
Blue, Pk/5. **HS29025B** \$23.75
Pink, Pk/5. **HS29025D** \$23.75
Mixed, Pk/5. **HS29025G** \$23.75
- 96-Well reversible (12 x 8), snap-on cover, 24.6 x 12.1 x 5 cm
Holds 1.5-2.0 mL tubes on one side or 0.5 mL on reverse.
Blue, each. **HS2345C** \$11.05
Green, each. **HS2345D** \$11.05
Pink, each. **HS2345E** \$11.05



Tube Racks, PCR, polypropylene, mixed colors, Heathrow
(Colors: blue, green, orange, purple & yellow)

- Reversible, 96-Well, Pk/5. **HS2344A** \$49.75
Holds 0.2, 0.5 and 1.5 mL tubes.
- 96-Well, embossed numbers and letters, Pk/5. **HS23461A** \$35.70
Holds 0.2 mL individual tubes and/or strips of 8-12 tubes.




Tube Racks, Peg, polypropylene, Heathrow

- 10-13 mm Tubes, 96-place
Blue, Pk/2. **HS24312B** \$23.45
White, Pk/2. **HS24312A** \$23.45
- 14-17 mm Tubes, 50-place
Blue, Pk/2. **HS24311B** \$21.30
White, Pk/2. **HS24311A** \$21.30



Tube Racks, Slant, polypropylene, Bel-Art
 13-16 mm Tubes, 40-place, each. **188900016**..... \$55.05
 16-20 mm Tubes, 40-place, each. **188900020**..... \$55.05




Tilt your test tube rack 25°, 35°, or 45°. Create agar slant tubes, or adjust the angle of tube racks when pipetting to reduce hand strain.

Tube Rack Slanter, Bel-Art. **188900100**..... \$19.50
 Platform size: 19.6 cm(W) x 12.7 cm(D)



Tube Rack, Stackable, polypropylene, Bel-Art
 1.5 mL Tubes, 72-place, each. **188410000**..... \$38.25



Tube Rack, SmoothRack™, polypropylene, blue, USA
 10-17 mm Tubes, 72-place, each. **SR-B** \$14.75




Remains stable in water bath.

Tube Racks, Fold & Snap, polypropylene, blue, Heathrow
 13 mm Tubes, 90-place, each. **HS243073B** \$11.85
 16 mm Tubes, 60-place, each. **HS243077B** \$11.85



Half-Size




Full-Size

Tube Racks, Teardrop, OneRack™, POM, Heathrow
 Half-Size
 16 mm Tubes, 36-place, 12.7 x 12.7 x 7 cm(H)
 Blue, Pk/8. **HS27502B** \$94.50
 Red, Pk/8. **HS27502D** \$94.50
 Full-Size
 13 mm Tubes, 72-place, 20.2 x 10.4 x 5.9 cm(H)
 Blue, Pk/8. **HS27511B** \$122.00
 Red, Pk/8. **HS27511D** \$122.00
 16 mm Tubes, 72-place, 25 x 12.7 x 7 cm(H)
 Blue, Pk/8. **HS27512B** \$122.00
 Red, Pk/8. **HS27512D** \$122.00
 20 mm Tubes, 40-place, 25.2 x 10 x 8.3 cm(H)
 Blue, Pk/8. **HS27513B** \$122.00
 Red, Pk/8. **HS27513D** \$122.00



Polypropylene cube and lid has two snap-in, gridded inserts to hold either (16) 15 mL tubes or (9) 50 mL tubes.

Tube Racks, Transformer Cube, Pk/5. **HS23271A**..... \$50.20




120-Place Wire Rack

Tube Racks, Wire, epoxy-coated, blue, Poxgrid®, Bel-Art
 10-13 mm Tubes, 120-place (6 x 20), each. **187860750**..... \$40.60
 Dimensions: 33.1 x 11 x 6 cm(H)
 15-16 mm Tubes
 100-Place (5 x 20), each. **187860780**..... \$41.25
 Dimensions: 38.6 x 10.8 x 6.4 cm(H)
 180-Place (9 x 20), each. **187860770**..... \$41.45
 Dimensions: 38.6 x 18.3 x 6.4 cm(H)



HS23215



17515ML-25


Tube Racks, 15 mL
 15-Place, epoxy-coated wire, ea. **HS23215**..... \$25.20
 25-Place (5 x 5), stainless steel, USA, ea. **17515ML-25**..... \$14.80
 80-Place (5 x 16), stainless steel, USA, ea. **17515ML-80**..... \$33.40
 Dimensions: 40.7 x 12.7 x 10.2 cm(H)
 110-Place (5 x 22), stainless steel, USA, ea. **17515ML-110**..... \$43.00
 Dimensions: 54.6 x 12.7 x 10.2 cm(H)



Wire



Friction-Fit



Interlocking

Tube Racks, 50 mL

1-Place, Heathrow

- Wire, epoxy-coated, each. **HS23051** \$11.45
- Friction-Fit, polystyrene, Pk/5. **HS23053A** \$16.45
- Interlocking, polypropylene, Pk/5. **HS23052** \$32.00



HS232461



17550ML-9

- 8-Place (2 x 4), epoxy-coated wire, ea. **HS232461** \$24.50
- 9-Place (3 x 3), stainless steel, USA, ea. **17550ML-9** \$14.80
- 16-Place (4 x 4), epoxy-coated wire, ea. **187940001** \$42.85
- 24-Place (3 x 8), polypropylene, ea. **189150000** \$34.20
- Dimensions: 29 x 12.7 x 7.5 cm(H)
- 30-Place (3 x 10), stainless steel, USA, ea. **17550ML-30** \$33.40
- Dimensions: 40.7 x 12.7 x 10.2 cm(H)
- 42-Place (3 x 14), stainless steel, USA, ea. **17550ML-42** \$43.00
- Dimensions: 54.6 x 12.7 x 10.2 cm(H)

Weighing Dishes



Hexagonal



Standard



Pour

Weigh Boats, Heathrow

Polystyrene
Standard

- Small, 4.6 x 4.6 x 0.8 cm(D), Bx/500. **HS1420A** \$19.85
- ★ Medium, 8.5 x 8.5 x 2.4 cm(D), Bx/500. **HS1420B** \$32.80
- Large, 14 x 14 x 2.2 cm(D), Bx/500. **HS1420C** \$53.40

Hexagonal, anti-static

- Small, 3.8 cm top dia., Bx/500. **HS14251A** \$14.35
- Medium, 6.4 cm top dia., Bx/500. **HS14251B** \$25.45
- Large, 12 x 14 cm top OD, Bx/500. **HS14251C** \$66.50
- (25433-104)
- X-Large, 14 x 16 cm top OD, Bx/500. **HS14251D** \$84.90
- (25433-106)

Pour (Canoes)

- Small, 5.8 x 4.3 x 1.3 cm(D), Bx/500. **HS1419A** \$33.30
- Medium, 13.2 x 8.3 x 2.6 cm(D), Bx/500. **HS1419B** ... \$65.75
- Large, 18.3 x 10.8 x 2.6 cm(D), Bx/500. **HS1419C** \$93.70



Aluminum, crimped sides w/tab

- Small, 43 mm bottom ID, Bx/100. **HS14521A** \$10.95
- Medium, 57 mm bottom ID, Bx/100. **HS14521B** \$12.30
- Large, 70 mm bottom ID, Bx/100. **HS14521C** \$16.85



199170001

Polypropylene vacuum bottles are autoclavable and withstand 700 mm Hg systems.



EV432

Stand sold separately.

Vacuum Waste Collection Systems

- 1.9 L (0.5 gal), Bel-Art. **199170001** \$191.00
- Attached trap protects vacuum source from overflow.
- 4 L (1 gal), Argos. **EV432** \$122.00
- Cap has traditional barbed fitting.



Square



Syracuse

Watch Glasses, flint glass, USA

- Square, ground front, 41 mm. **910-B** \$7.90
- Syracuse, grooved for stacking, 67 mm OD. **913-A** \$5.85

Wipes



Extra low-lint, nonabrasive, highly absorbent, 1-ply wipes are perfect for microscope maintenance.

- | | | | |
|------------------------------------|---------|---------------|-----------------|
| Kimwipes® | | <u>bx/280</u> | <u>cs/60 bx</u> |
| 4½ x 8½" Wipes. 34155 | \$3.55 | 190.00 | |
| | | <u>bx/140</u> | <u>cs/15 bx</u> |
| 15 x 17" Wipes. 34256 | \$10.75 | 149.00 | |

Continued