


glass w/rubber bulb; 75 mm (3")

Medicine (Eye) Dropper, Pk/12. 97-4301 \$3.50



Pipet Fillers

- 3-Valve, rubber, black. 207B \$14.40
- Wheel-Type, Pipette Pump™, Bel-Art
 - Fits Pipets Up To 2 mL, blue. 378970000 \$23.15
 - Fits Pipets Up To 10 mL, green. 378980000 \$24.20
 - Fits Pipets Up To 25 mL, red. 378990000 \$33.30
- Wheel-Type, Pipette Pump™ II, quick-release ea pk/12
 - Fits Pipets Up To 2 mL, blue. 379111002 \$25.05 265.00
 - Fits Pipets Up To 10 mL, green. 379111010 \$25.05 265.00
 - Fits Pipets Up To 25 mL, red. 379111025 \$33.30 353.00




Motorized. HS3000 \$235.00



Dimensions: 29 x 8.8 x 40 cm

Pipet Rack, magnetized, blue, Heathrow. HS20612M \$71.40



Store up to 4 sterilization cans with this epoxy-coated wire rack. (Pipet cans not included)

Pipet Sterilizing Box Rack, Bel-Art. 189420000 \$44.60



Pipet Tip Disposal Box, Heathrow. HS234634 \$83.65




Dispense fluids accurately and safely with these color-coded and plugged pipets. Negative graduations provide extra volume (not available on 25 mL). RNase-/DNase-/Pyrogen-free, and noncytotoxic. Calibrated "To Deliver" (TD).

Pipets, Serological, polystyrene, UniPlast™

- Nonsterile
 - 5 mL, blue band, pk/25, Cs/375. 1755 \$73.15
 - 10 mL, orange band, pk/25, Cs/250. 1775 \$50.00
 - 25 mL, red band, pk/10, Cs/100. 1784 \$46.55
- Sterile Pk/1 (individually paper/plastic wrapped)
 - 1 mL, yellow band, Cs/500. 1700 \$73.45
 - 2 mL, green band, Cs/500. 1720 \$78.80
 - 5 mL, blue band, reverse grad., Cs/250. 1740 \$58.10
 - 10 mL, orange band, reverse grad., Cs/200. 1760 \$50.65
 - 25 mL, red band, reverse grad., Cs/100. 1780 \$54.75
 - 50 mL, black band, reverse grad., Cs/100. 1820 \$139.00



5 mL Capacity, 1 mL graduated



7 mL Capacity, 3 mL graduated

Pipets, Transfer, LDPE (low-density polyethylene)

- Graduated, Nonsterile, 15 cm long
 - 5 mL capacity, 1 mL graduated, 25 drops/mL
 - Bg/25. PTP-02 \$1.00
 - Bx/500. PTP-02 \$15.65
 - Cs/10 Bx. PTP-02 \$145.00
 - 7 mL capacity, 3 mL graduated, 14 drops/mL
 - Bg/25. PTP-01 \$1.00
 - Bx/500. PTP-01 \$15.65
 - Cs/10 Bx. PTP-01 \$145.00
- Graduated, Sterile Pk/1 (individually wrapped), 15 cm long
 - 5 mL capacity, 1 mL graduated, Simport
 - Bx/400. P200-521S \$39.65
 - Cs/10 Bx. P200-521S \$364.00
- Graduated, Sterile Pk/20, 15 cm long
 - 5 mL capacity, 1 mL graduated
 - Bx/500. HS206371B \$32.85
 - 5 mL capacity, 2 mL graduated
 - Bx/400. P200-5220S (HS206372B) \$23.80
- Long Stem, Sterile Pk/1 (individually wrapped)
 - 22.5 cm (9" long), 6 mL capacity
 - Cs/400. 139038 \$47.70
- Long Stem, Nonsterile
 - 30 cm (12" long), 23 mL capacity
 - Bx/100. 139050 \$17.20



Made of borosilicate glass in compliance with ASTM specification E-969, Class A requirements. Individually serialized and certified. Calibrated "To Deliver" (TD).



Pipets, Volumetric, serialized, Class A

- 2 mL (±0.006 mL), orange band, Pk/5. 600PT7100B \$64.75
- 5 mL (±0.010 mL), white band, Pk/5. 600PT7100C \$68.05
- 10 mL (±0.020 mL), red band, Pk/5. 600PT7100D \$77.75
- 20 mL (±0.030 mL), yellow band, Pk/5. 600PT7100E \$110.00

Continued