


Rockers, Shakers & Mixers



Rockers, Fixed Speed, Benchmark (2-year warranty)

2-D (seesaw), 12 x 12" platform w/flat mat

BenchBlotter™. **BR1000** \$413.00

- Fixed speed of 12 rocks/minute
- Platform angle is $\pm 8^\circ$
- Cold room & incubator safe (4-45°C)
- Load capacity is 2.3 kg (5 lbs)

3-D (nutating)

Small, 10.5 x 7.5" platform w/dimpled mat

MiniMixer™. **B3D1020** \$350.00

Large, 12 x 12" platform w/dimpled mat

BioMixer™. **B3D1320** \$422.00

- Fixed speed of 24 rocks/minute
- Platform angle is $\pm 20^\circ$
- Cold room & incubator safe (4-45°C)
- Load capacity is 1.8 kg (4 lbs)



Rockers, Variable Speed

Benchmark (2-year warranty)

2-D (seesaw), 14 x 12" platform w/flat mat

BenchRocker™ 2D. **BR2000** \$679.00

- Variable speed of 2-30 rocks/minute
- Platform tilt angle adjusts from 0-30°
- Cold room & incubator safe (4-45°C)
- Load capacity is 3.6 kg (8 lbs)

3-D (nutating), 12 x 12" platform w/dimpled and flat mats

BenchRocker™ 3D. **B3D2300** \$722.00

Stacking Platform. **BR1000-STACK** \$94.00

12 x 12" platform w/flat mat has 3" of clearance.

Stacking Platform Extenders, Pk/4. **BR2000-SP** \$36.10

Adds 1.25" to platform separation.

- Variable speed of 2-30 rocks/minute
- Platform tilt angle adjusts from 0-30°
- Cold room & incubator safe (4-45°C)
- Load capacity is 2.3 kg (5 lbs)



Rockers, Variable Speed

Reliable Scientific (5-year warranty), USA

2-D (seesaw)

Standard, 11 x 14" platform w/rubber mat

Single Platform. **55S1114** \$576.00

Double Platform, 3/4" apart. **55D1114** \$634.00

Large, 12 x 16" platform w/rubber mat

Single Platform. **55S1216** \$615.00

Double Platform, 3/4" apart. **55D1216** \$672.00

- Variable speed of 40-160 rocks/minute
- Platform tilt angle adjusts from 0-20°
- Cold room & incubator safe (0-65°C)
- Load capacity is 4.5 kg (10 lbs)



Rocker, Tube, 12 x 4.5" grooved mat, Benchmark. **M2100** \$309.00

- Fixed speed of 24 rocks/minute
- Platform tilt angle is $\pm 30^\circ$
- Coldroom & incubator safe (4-45°C)
- Accommodates 1.5 to 50 mL tubes up to 2.3 kg (5 lbs)
- 2-year warranty

Orbital Shakers



Aerate and mix with these low to high capacity orbital shakers.



Orbital Shakers, BenchMark (2-year warranty)

Low Speed, 14 x 12" platform w/flat mat

Orbi-Blotter™. **BT30** \$710.00

- Speed adjusts from 3-70 rpm
- Orbit diameter measures 19 mm
- Cold room & incubator safe (4-65°C)
- Load capacity is 3 x 2 L, or 4 kg (9 lbs)

Small, 9.5 x 8" platform w/flat mat

Orbi-Shaker™ JR. **BT302** \$809.00

- Speed adjusts from 50-300 rpm in 10 rpm increments
- Orbit diameter measures 19 mm
- Cold room & incubator safe (4-65°C)
- Load capacity is 4 x 1 L, or 4 kg (9 lbs)
- Timer is programmable from 1 min - 99 hrs, 59 min

Continued