


Rockers, Shakers & Mixers



Rockers, Fixed Speed, Benchmark (2-year warranty)

2-D (seesaw), 12 x 12" platform w/flat mat

BenchBlotter™. **BR1000** \$413.00

- Fixed speed of 12 rocks/minute
- Platform angle is $\pm 8^\circ$
- Cold room & incubator safe (4-45°C)
- Load capacity is 2.3 kg (5 lbs)

3-D (nutating)

Small, 10.5 x 7.5" platform w/dimpled mat

MiniMixer™. **B3D1020** \$350.00

Large, 12 x 12" platform w/dimpled mat

BioMixer™. **B3D1320** \$422.00

- Fixed speed of 24 rocks/minute
- Platform angle is $\pm 20^\circ$
- Cold room & incubator safe (4-45°C)
- Load capacity is 1.8 kg (4 lbs)



Rockers, Variable Speed

Benchmark (2-year warranty)

2-D (seesaw), 14 x 12" platform w/flat mat

BenchRocker™ 2D. **BR2000** \$679.00

- Variable speed of 2-30 rocks/minute
- Platform tilt angle adjusts from 0-30°
- Cold room & incubator safe (4-45°C)
- Load capacity is 3.6 kg (8 lbs)

3-D (nutating), 12 x 12" platform w/dimpled and flat mats

BenchRocker™ 3D. **B3D2300** \$722.00

Stacking Platform. **BR1000-STACK** \$94.00

12 x 12" platform w/flat mat has 3" of clearance.

Stacking Platform Extenders, Pk/4. **BR2000-SP** \$36.10

Adds 1.25" to platform separation.

- Variable speed of 2-30 rocks/minute
- Platform tilt angle adjusts from 0-30°
- Cold room & incubator safe (4-45°C)
- Load capacity is 2.3 kg (5 lbs)



Rockers, Variable Speed

Reliable Scientific (5-year warranty), USA

2-D (seesaw)

Standard, 11 x 14" platform w/rubber mat

Single Platform. **55S1114** \$576.00

Double Platform, 3/4" apart. **55D1114** \$634.00

Large, 12 x 16" platform w/rubber mat

Single Platform. **55S1216** \$615.00

Double Platform, 3/4" apart. **55D1216** \$672.00

- Variable speed of 40-160 rocks/minute
- Platform tilt angle adjusts from 0-20°
- Cold room & incubator safe (0-65°C)
- Load capacity is 4.5 kg (10 lbs)



Rocker, Tube, 12 x 4.5" grooved mat, Benchmark. **M2100** \$309.00

- Fixed speed of 24 rocks/minute
- Platform tilt angle is $\pm 30^\circ$
- Coldroom & incubator safe (4-45°C)
- Accommodates 1.5 to 50 mL tubes up to 2.3 kg (5 lbs)
- 2-year warranty

Orbital Shakers



Aerate and mix with these low to high capacity orbital shakers.



Orbital Shakers, Benchmark (2-year warranty)

Low Speed, 14 x 12" platform w/flat mat

Orbi-Blotter™. **BT30** \$710.00

- Speed adjusts from 3-70 rpm
- Orbit diameter measures 19 mm
- Cold room & incubator safe (4-65°C)
- Load capacity is 3 x 2 L, or 4 kg (9 lbs)

Small, 9.5 x 8" platform w/flat mat

Orbi-Shaker™ JR. **BT302** \$809.00

- Speed adjusts from 50-300 rpm in 10 rpm increments
- Orbit diameter measures 19 mm
- Cold room & incubator safe (4-65°C)
- Load capacity is 4 x 1 L, or 4 kg (9 lbs)
- Timer is programmable from 1 min - 99 hrs, 59 min

Continued



BT3000

Large, 13 x 12" platform w/flat mat

Orbi-Shaker™. **BT3000** \$1,277.00

- Speed adjusts from 30-300 rpm in 1 rpm increments
- Orbit diameter measures 19 mm
- Cold room & incubator safe (4-65°C)
- Load capacity is 9 x 500 mL, or 5 x 1 L, or 2 x 2 L,
- or 6.4 kg (14 lbs)
- Timer is programmable from 1 min - 48 hrs/continuous

Optional Magic Clamp™ platform & accessories sold separately - specify configuration

Steel Platform, Magic Clamp™, 14 x 12". **BT3000-MR**... \$208.00

250 mL Flask Magnetic Clamp. **H1000-MR-250** \$28.05

500 mL Flask Magnetic Clamp. **H1000-MR-500** \$29.45

1000 mL Flask Magnetic Clamp. **H1000-MR-1000** \$45.60

2000 mL Flask Magnetic Clamp. **H1000-MR-2000** \$60.35



Tube Rotator



Tube Rotator, Reliable Scientific, USA. **Model 24** \$349.00

- Speed is adjustable from 4-60 rpm
- Cold room & incubator safe (0-65°C)
- 5-year warranty

Rotator disks (carousels) sold separately - specify configuration

24-Place Disk, holds 9-15 mm tubes. **24/9-15** \$89.00

12-Place Disk, holds 15-23 mm tubes. **12/15-23** \$89.00

6-Place Disk, holds 23-29 mm tubes. **6/23-29** \$89.00

Vortexers



BV101



BV1000

Vortex Mixers, Benchmark (2-year warranty)

Vornado™, mini vortexer, fixed 2800 rpm. **BV101** ___* \$179.00

*Specify Gray, Blue or Purple cup head

BenchMixer™, single cup head, 200-3200 rpm. **BV1000** \$216.00

Flat Head Attachment, 3". **BV1000-FLAT** \$28.00

Rulers, Meter Sticks & Tape Measure

Stainless Steel



Rulers, 15 cm/6" long, scaled in both metric and English units

Vinyl, flexible

Transparent. **4461** \$0.95

White. **4460** \$0.55

Stainless Steel, 12 mm wide, USA. **605-065** \$5.50

Meter Sticks, lacquered wood

One side is graduated on both edges in centimeters with 1 mm divisions, while the other side is graduated on one edge in inches with 1/8" divisions.

Plain Ends. **605-050** \$4.90

Metal Ends. **605-051** \$6.15

Tape Measure, 152.4 cm/60", reinforced fiberglass. **11036** \$3.70

One side is graduated on the top edge in inches with 1/16" divisions and on the bottom edge in cm with 1 mm divisions. The opposite side is graduated in inches with 1/16" divisions.



C300-1R



4402269

Sharps Container, red polypropylene

1 L, Biodisposer™, Simport, each. **C300-1R** \$3.15

Dimensions: 8.9 x 8.9 x 15.25 cm (3 x 3 x 6" H)

5.1 L, horizontal entry lid, Prevent™, each. **4402269** \$7.95

Dimensions: 12" x 4.75" x 11"(H)



Screw-Cap



Snap-Cap

Specimen Containers, polypropylene (PP)

Screw-Cap

60 mL, nonsterile, Pk/500. **840095**..... \$102.00

90 mL, sterile, Pk/300. **840175** \$80.70

Snap-Cap, white PP bottom/white LDPE lid, Simport

300 mL (10 oz), Cs/100. **C580-3** \$42.90

11.1 cm opening x 4.5 cm high

500 mL (17 oz), Cs/100. **C580-4** \$55.00

11.1 cm opening x 7.5 cm high

1000 mL (34 oz), Cs/100. **271032W**..... \$82.25

11.1 cm opening x 14 cm high.