



BT3000

Large, 13 x 12" platform w/flat mat

- Orbi-Shaker™. **BT3000** \$1,277.00
- Speed adjusts from 30-300 rpm in 1 rpm increments
- Orbit diameter measures 19 mm
- Cold room & incubator safe (4-65°C)
- Load capacity is 9 x 500 mL, or 5 x 1 L, or 2 x 2 L,
- or 6.4 kg (14 lbs)
- Timer is programmable from 1 min - 48 hrs/continuous

Optional Magic Clamp™ platform & accessories sold separately - specify configuration

- Steel Platform, Magic Clamp™, 14 x 12". **BT3000-MR**... \$208.00
- 250 mL Flask Magnetic Clamp. **H1000-MR-250** \$28.05
- 500 mL Flask Magnetic Clamp. **H1000-MR-500** \$29.45
- 1000 mL Flask Magnetic Clamp. **H1000-MR-1000** \$45.60
- 2000 mL Flask Magnetic Clamp. **H1000-MR-2000** \$60.35



Tube Rotator



- Tube Rotator, Reliable Scientific, USA. **Model 24** \$349.00
- Speed is adjustable from 4-60 rpm
- Cold room & incubator safe (0-65°C)
- 5-year warranty

Rotator disks (carousels) sold separately - specify configuration

- 24-Place Disk, holds 9-15 mm tubes. **24/9-15** \$89.00
- 12-Place Disk, holds 15-23 mm tubes. **12/15-23** \$89.00
- 6-Place Disk, holds 23-29 mm tubes. **6/23-29** \$89.00

Vortexers



BV101



BV1000

Vortex Mixers, Benchmark (2-year warranty)

- Vornado™, mini vortexer, fixed 2800 rpm. **BV101** ___* \$179.00
- *Specify Gray, Blue or Purple cup head
- BenchMixer™, single cup head, 200-3200 rpm. **BV1000** \$216.00
- Flat Head Attachment, 3". **BV1000-FLAT** \$28.00

Rulers, Meter Sticks & Tape Measure

Stainless Steel



Rulers, 15 cm/6" long, scaled in both metric and English units

- Vinyl, flexible
- Transparent. **4461** \$0.95
- White. **4460** \$0.55
- Stainless Steel, 12 mm wide, USA. **605-065** \$5.50

Meter Sticks, lacquered wood

One side is graduated on both edges in centimeters with 1 mm divisions, while the other side is graduated on one edge in inches with 1/8" divisions.

- Plain Ends. **605-050** \$4.90
- Metal Ends. **605-051** \$6.15

Tape Measure, 152.4 cm/60", reinforced fiberglass. **11036** \$3.70

One side is graduated on the top edge in inches with 1/16" divisions and on the bottom edge in cm with 1 mm divisions. The opposite side is graduated in inches with 1/16" divisions.



C300-1R



4402269

Sharps Container, red polypropylene

- 1 L, Biodisposer™, Simport, each. **C300-1R** \$3.15
- Dimensions: 8.9 x 8.9 x 15.25 cm (3 x 3 x 6" H)
- 5.1 L, horizontal entry lid, Prevent™, each. **4402269** \$7.95
- Dimensions: 12" x 4.75" x 11"(H)



Screw-Cap



Snap-Cap

Specimen Containers, polypropylene (PP)

- Screw-Cap
- 60 mL, nonsterile, Pk/500. **840095**..... \$102.00
- 90 mL, sterile, Pk/300. **840175** \$80.70
- Snap-Cap, white PP bottom/white LDPE lid, Simport
- 300 mL (10 oz), Cs/100. **C580-3** \$42.90
- 11.1 cm opening x 4.5 cm high
- 500 mL (17 oz), Cs/100. **C580-4** \$55.00
- 11.1 cm opening x 7.5 cm high
- 1000 mL (34 oz), Cs/100. **271032W**..... \$82.25
- 11.1 cm opening x 14 cm high.



Jewelry, sterling silver

Charms


- Caduceus, 7/8" long. **88501** \$11.00
- Chimpanzee, 1 1/16" long. **88503** \$12.00
- Microscope, 3/4" long. **88505** \$11.00
- Mouse, 1" wide. **88507** \$12.00

Earrings, pair

- Stud, microscope, 3/4" long. **88521** \$34.00
- Wire Hoops**
- Microscope, 3/4" long. **88531** \$25.00
- Mouse, 1" wide. **88533** \$27.00
- Pin, Microscope. **88541** \$17.00

Micro Bible. 89585 \$15.00

View with a 10X objective all 1,245 pages of the King James version of the Holy Bible on this approximately 1" square microfilm.



Ruler, 1 Decimeter, stainless steel, Germany. 605-060 \$10.20

Ruler is 1 dm long, 1 cm wide, and 1 mm thick. Its volume is 1 cm³ and its weight in grams is the specific gravity. One edge of one side is graduated in centimeters with 1 mm divisions.



This beautiful, transparent ornament is a collage of cytological color photomicrographs creatively captured on a 3" glass disk. Suspended from a gold braid.

Suncatcher, "Cytological Crash Course." CCC01 \$22.50



Tube-Wringer®, Heavy-Duty, USA. **295401** \$24.80
3" aluminum rollers with steel handle