

Illuminators

Please call regarding availability of items not listed



Fiber Optics, Dolan-Jenner

- Light Source, 150W EKE lamp, rheostat. **FL150/115** \$392.00
(MI-150)
- Gooseneck, Dual, chrome plated, 23" long. **FL150/05**..... \$230.00
(EEG2823M)
- Focusing Lens. **FL150/50**..... \$39.00
(LH-759)



- LED Illuminator, 5W LED. MA264/LED**..... \$431.00

Accessories & Parts

Cushion and support your elbows/forearms with a pair of these ergonomic armrests. Pliable, latex bag stays firmly in position on the workbench, where it conforms to surface and your arm.



- Armrests, latex rubber bag, black, USA** ea pair
- ★ Low-Profile, 1.25 lbs each. **14146** \$12.00 22.00
Recommended for newer 'scopes with low controls
- High-Profile, 2 lbs each. **14141** \$13.00 23.00



Boxes, Plastic, hinged polystyrene, USA

Catalog #	Inside Dimensions	Foam Insert	Suggested Application	Price
H-1	1 x 1 x 1/4"	-	Reticles, Pointers	\$0.30
H-1/11B	1 x 1 x 1/4"	1	Reticles, Pointers	\$0.35
H-54	2 9/16 x 2 9/16 x 1/2"	-	Filters	\$1.20
H-54/2916G	2 9/16 x 2 9/16 x 1/2"	1	Filters	\$1.35
H-28/334B2	3 3/4 x 1 3/16 x 5/8"	2	Slides, Stage Micrometers	\$1.40




Transport and store your microscope in these protective cases. Velcro® straps and accessory pockets are provided.

Carrying Cases, hard-sided, anodized aluminum w/fabric lining (Inside dimensions provided)

- 10.5"(D) x 7.75"(W) x 14.9"(H). **975-001** \$129.00
- 11.75"(D) x 8"(W) x 17.25"(H). **975-002** \$139.00

Counting Chambers

Bright-Line®
Hemocytometer



Fuchs-Rosenthal Chamber, USA

- Chamber with 2 Coverslips, Hausser. **3720** \$375.00
- Coverslips only, 20 x 26 x 0.4 mm, Pk/2, Erie. **2026H** \$10.80

Hemocytometer, bright-lined, V-slash, Germany

- Chamber with 2 Coverslips, Marienfeld. **0650030** \$134.00

Hemocytometers, Bright-Line®, molded borosilicate glass, USA

- Standard Thickness (5.5 mm thick)
- Chamber with 2 Coverslips, Hausser. **1492** \$283.00
- Chamber only, Hausser. **1490** \$255.00
- Coverslips, 20 x 26 x 0.4 mm, Pk/2, Erie. **2026H** \$10.80
- Reduced Thickness/phase counting chamber (3.5 mm thick)
- Chamber with 2 Coverslips, Hausser. **1475** \$302.00
- Coverslips, 20 x 26 x 0.4 mm, Pk/2, Erie. **2026H** \$10.80

Howard Mold Chamber, Hausser, USA


Requires Howard Disc - see RETICLES & STAGE MICROMETERS

- Chamber with 2 Coverslips (5080 & 5090). **3820** \$375.00
- Coverslip only, 28 x 33 x 0.5 mm (thin). **5080** \$46.80
- Coverslip only, 28 x 33 x 1.0 mm (thick). **5090** \$46.80

Petroff-Hausser Bacteria Counter, Hausser, USA

- Chamber with 2 Coverslips (5051 & 1 grid). **3900** \$578.00
- Coverslip only, 20 x 26 x 0.2 mm. **5051** \$46.80

- England Field Finder, 1 mm thick, USA. KR-834** \$159.00

- Eyepiece Pointer, wire. 14301**  \$2.00



- Eyetube Sleeve, 23.2 mm increased to 30 mm. 14320** \$35.00