

Stoppers

Stoppers, rubber, black



Stopper Sizing Guide			
Size #	Top (mm)	Bottom (mm)	Fits ID Openings (mm)
1	19	14	15-17
2	20	16	16-18.5
3	24	18	18-21
4	26	20	20-23
5	27	23	23-25
6	32	26	26-27
6½	34	27	27-31.5
7	37	30	30-34
8	41	33	33-37
9	45	37	37-41
10	50	42	42-46
12	64	54	54-59

All stoppers are 25 mm long. Stopper size #1 has a 4 mm hole; sizes #2 and larger have 5 mm holes.

Solid	pk/5	lb
Size 2, ~39 pcs/lb. 50020	\$2.95	12.85
Size 3, ~29 pcs/lb. 50030	\$4.00	12.85
Size 4, ~23 pcs/lb. 50040	\$5.00	12.85
Size 5, ~20 pcs/lb. 50050	—	12.85
Size 6, ~15 pcs/lb. 50060	—	12.85
Size 6.5, ~14 pcs/lb. 500650	—	12.85
Size 7, ~11 pcs/lb. 50070	—	12.85
Size 8, ~9 pcs/lb. 50080	—	12.85
1-Hole	pk/5	lb
Size 1, ~51 pcs/lb. 50011	\$2.25	12.85
Size 2, ~43 pcs/lb. 50021	\$2.70	12.85
★ Size 3, ~30 pcs/lb. 50031	\$3.85	12.85
Size 4, ~24 pcs/lb. 50041	\$4.80	12.85
Size 5, ~21 pcs/lb. 50051	—	12.85
Size 6, ~18 pcs/lb. 50061	—	12.85
Size 9, ~8 pcs/lb. 50091	—	12.85
Size 12, ~3 pcs/lb. 50121	—	12.85
2-Hole	lb	
Size 6, ~17 pcs/lb. 50062	\$12.85	
Size 7, ~12 pcs/lb. 50072	\$12.85	
Size 8, ~9 pcs/lb. 50082	\$12.85	
Size 9, ~7 pcs/lb. 50092	\$12.85	
Size 10, ~5 pcs/lb. 50102	\$12.85	
Assortment, 1 lb. 50012	\$17.10	

Solid, 1-hole, and 2-hole stoppers in sizes 1-6.



Accurately determine the correct size of rubber stopper, size 00-12, or cork, size 0-26, for any vessel using this durable plastic tool. Reverse side has a circular scale for sizing 4-52 mm OD tubing.

Stopper/Tubing Sizer. 133920000 \$12.55



Surgical Drape, nonferrosted, 18 x 26", Bx/50. 44025-515..... \$25.50
sterile, poly-lined

Tape

Tape



Autoclave, steam
¾" x 60 Yd Roll. 44073ST048 \$8.50
1" x 60 Yd Roll. 44073ST036 \$10.70

Kraft paper with rubber resin adhesive.



Label, Write-On™, Bel-Art
1" Core, 15 yd (540") roll
½" Wide, each \$5.35
Green. 134623005 Yellow. 134622005
White. 134620005
¾" Wide, Pk/4 \$23.20
Green. 134820075 Yellow. 134810075
White. 134800075
3" Core, 40 yd roll
½" Wide, each \$9.80
Green. 134633005 White. 134630005
½" Wide, Pk/6 \$55.75
Red. 134840050 Yellow. 134850050
¾" Wide, each \$13.00
Green. 134633075 White. 134630075
Yellow. 134632075
¾" Wide, Pk/4 \$49.35
White. 134830075 Red. 134840075
¾" Wide, Pk/6. 134630600 \$75.25
Mixed: blue, green, orange, red, white & yellow



Teflon®, ½ x 300", USA, each. 240200001 \$4.20