

**Stoppers**

Stoppers, rubber, black



Stopper Sizing Guide			
Size #	Top (mm)	Bottom (mm)	Fits ID Openings (mm)
1	19	14	15-17
2	20	16	16-18.5
3	24	18	18-21
4	26	20	20-23
5	27	23	23-25
6	32	26	26-27
6½	34	27	27-31.5
7	37	30	30-34
8	41	33	33-37
9	45	37	37-41
10	50	42	42-46
12	64	54	54-59

All stoppers are 25 mm long. Stopper size #1 has a 4 mm hole; sizes #2 and larger have 5 mm holes.

Solid	pk/5	lb
Size 2, ~39 pcs/lb. 50020	\$2.95	12.85
Size 3, ~29 pcs/lb. 50030	\$4.00	12.85
Size 4, ~23 pcs/lb. 50040	\$5.00	12.85
Size 5, ~20 pcs/lb. 50050	—	12.85
Size 6, ~15 pcs/lb. 50060	—	12.85
Size 6.5, ~14 pcs/lb. 500650	—	12.85
Size 7, ~11 pcs/lb. 50070	—	12.85
Size 8, ~9 pcs/lb. 50080	—	12.85
1-Hole	pk/5	lb
Size 1, ~51 pcs/lb. 50011	\$2.25	12.85
Size 2, ~43 pcs/lb. 50021	\$2.70	12.85
★ Size 3, ~30 pcs/lb. 50031	\$3.85	12.85
Size 4, ~24 pcs/lb. 50041	\$4.80	12.85
Size 5, ~21 pcs/lb. 50051	—	12.85
Size 6, ~18 pcs/lb. 50061	—	12.85
Size 9, ~8 pcs/lb. 50091	—	12.85
Size 12, ~3 pcs/lb. 50121	—	12.85
2-Hole	lb	
Size 6, ~17 pcs/lb. 50062	\$12.85	
Size 7, ~12 pcs/lb. 50072	\$12.85	
Size 8, ~9 pcs/lb. 50082	\$12.85	
Size 9, ~7 pcs/lb. 50092	\$12.85	
Size 10, ~5 pcs/lb. 50102	\$12.85	
Assortment, 1 lb. 50012	\$17.10	

Solid, 1-hole, and 2-hole stoppers in sizes 1-6.



Accurately determine the correct size of rubber stopper, size 00-12, or cork, size 0-26, for any vessel using this durable plastic tool. Reverse side has a circular scale for sizing 4-52 mm OD tubing.

Stopper/Tubing Sizer. 133920000 ..... \$12.55



Surgical Drape, nonfenestrated, 18 x 26", Bx/50. 44025-515..... \$25.50  
sterile, poly-lined

**Tape**

Tape



Autoclave, steam

¾" x 60 Yd Roll. 44073ST048 ..... \$8.50  
1" x 60 Yd Roll. 44073ST036 ..... \$10.70

Kraft paper with rubber resin adhesive.



Label, Write-On™, Bel-Art

1" Core, 15 yd (540") roll  
½" Wide, each ..... \$5.35  
Green. 134623005 Yellow. 134622005  
White. 134620005  
¾" Wide, Pk/4 ..... \$23.20  
Green. 134820075 Yellow. 134810075  
White. 134800075  
3" Core, 40 yd roll  
½" Wide, each ..... \$9.80  
Green. 134633005 White. 134630005  
½" Wide, Pk/6 ..... \$55.75  
Red. 134840050 Yellow. 134850050  
¾" Wide, each ..... \$13.00  
Green. 134633075 White. 134630075  
Yellow. 134632075  
¾" Wide, Pk/4 ..... \$49.35  
White. 134830075 Red. 134840075  
¾" Wide, Pk/6. 134630600 ..... \$75.25  
Mixed: blue, green, orange, red, white & yellow



Teflon®, ½ x 300", USA, each. 240200001 ..... \$4.20

**Fuses**

¼ x 1¼"



Fast-Acting, AGC	
0.1A/250V. 312-100 .....	\$2.00
0.2A/250V. 312-200 .....	\$2.00
1A/250V. 312-001 .....	\$2.00
2A/250V. 312-002 .....	\$2.00
3A/250V. 312-003 .....	\$2.00
4A/250V. 312-004 .....	\$2.00

Time-Delay (Slow-Blow), MDL	
0.7A/250V. 313-700 .....	\$2.00
1A/250V. 313-001 .....	\$2.00
2.5A/250V .....	\$2.00

5 x 20 mm

Fast-Acting, GMA	
0.5A/250V. 235-500 .....	\$2.00
1A/250V. 235-001 .....	\$2.00
3A/250V. 235-003 .....	\$2.00
5A/125V. 235-005 .....	\$2.00

Time-Delay (Slow-Blow), Glass tube	
0.5A/250V. 239-500 .....	\$2.00
1A/250V. 239-001 .....	\$2.00
2A/250V. 239-002 .....	\$2.00
3A/250V. 239-003 .....	\$2.00
5A/250V. 239-005 .....	\$2.00

Time-Delay (Slow-Blow), Ceramic	
1A/250V. 215-001 .....	\$2.00



Cone Point



Cup Point

**Set Screws, steel**

Cone Point, DIN 914

M4 x 0.7 x 5 mm length, 2 mm hex key size. 77702 .....	\$0.50
M4 x 0.7 x 10 mm length, 2 mm hex key size. 77704 .....	\$0.50
M4 x 0.7 x 12 mm length, 2.5 mm hex key size. 77706 .....	\$0.50
M5 x 0.7 x 10 mm length, 2.5 mm hex key size. 77708 .....	\$0.50

Cup Point, DIN 916

M3 x 0.5 x 3 mm length, 1.5 mm hex key size. 77722 .....	\$0.50
M3 x 0.5 x 5 mm length, 1.5 mm hex key size. 77724 .....	\$0.50

**Tools**



Cutter, diagonal, 5". PLR-766.00 ..... \$10.95




Demagnetize or magnetize parts and tools.

Demagnetizer/Magnetizer, 1½ x 3" opening. DEM-345.00 ..... \$80.95

**Gauges**



Vernier Caliper



Compasses

Vernier Caliper, digital, 3-button. 600VCD001 ..... \$28.80  
Measures in millimeters (0-150 mm with 0.01 mm precision) and inches (0-6"). Provided with a plastic case.

Compasses, replaceable tips, Germany

3" Long. DIV-320.00 .....	\$11.95
6" Long. DIV-326.00 .....	\$16.95
Tips (replacement), Pair. DIV-320.05 .....	\$1.95



Reproducibility enhanced by rapid stop. Provided with padded, plastic case.

Micrometer, 0-25 mm in 0.01 mm div, Germany. 601A ..... \$39.00



Measure inside and outside diameters of glass and plastic tubing, connecting pieces, rods, test tubes, etc. between 1.5 and 16 mm. Bel-Art

Tubing Sizer, metric diameter gauge, PP. 133900000 ..... \$33.35

*Continued*