

## Illuminators

*Please call regarding availability of items not listed*



### Fiber Optics, Dolan-Jenner

Light Source, 150W EKE lamp, rheostat. **FL150/115** ..... \$392.00  
(MI-150)

Gooseneck, Dual, chrome plated, 23" long. **FL150/05**..... \$230.00  
(EEG2823M)

Focusing Lens. **FL150/50**..... \$39.00  
(LH-759)




**LED Illuminator, 5W LED. MA264/LED**..... \$431.00

## Accessories & Parts

Cushion and support your elbows/forearms with a pair of these ergonomic armrests. Pliable, latex bag stays firmly in position on the workbench, where it conforms to surface and your arm.



**Armrests, latex rubber bag, black, USA** \_ea\_    \_pair\_  
★ Low-Profile, 1.25 lbs each. **14146** ..... \$12.00    22.00  
*Recommended for newer 'scopes with low controls*  
High-Profile, 2 lbs each. **14141** ..... \$13.00    23.00



### Boxes, Plastic, hinged polystyrene, USA

Catalog #	Inside Dimensions	Foam Insert	Suggested Application	Price
H-1	1 x 1 x 1/4"	-	Reticles, Pointers	\$0.30
H-1/11B	1 x 1 x 1/4"	1	Reticles, Pointers	\$0.35
H-54	2 9/16 x 2 9/16 x 1/2"	-	Filters	\$1.20
H-54/2916G	2 9/16 x 2 9/16 x 1/2"	1	Filters	\$1.35
H-28/334B2	3 3/4 x 1 3/16 x 5/8"	2	Slides, Stage Micrometers	\$1.40



Transport and store your microscope in these protective cases. Velcro® straps and accessory pockets are provided.

### Carrying Cases, hard-sided, anodized aluminum w/fabric lining


(Inside dimensions provided)

10.5"(D) x 7.75"(W) x 14.9"(H). **975-001** ..... \$129.00

11.75"(D) x 8"(W) x 17.25"(H). **975-002** ..... \$139.00

## Counting Chambers

Bright-Line®  
Hemocytometer



### Fuchs-Rosenthal Chamber, USA

Chamber with 2 Coverslips, Hausser. **3720** ..... \$375.00

Coverslips only, 20 x 26 x 0.4 mm, Pk/2, Erie. **2026H** ..... \$10.80

### Hemocytometer, bright-lined, V-slash, Germany

Chamber with 2 Coverslips, Marienfeld. **0650030** ..... \$134.00

### Hemocytometers, Bright-Line®, molded borosilicate glass, USA

Standard Thickness (5.5 mm thick)

Chamber with 2 Coverslips, Hausser. **1492** ..... \$283.00

Chamber only, Hausser. **1490** ..... \$255.00

Coverslips, 20 x 26 x 0.4 mm, Pk/2, Erie. **2026H** ..... \$10.80

Reduced Thickness/phase counting chamber (3.5 mm thick)

Chamber with 2 Coverslips, Hausser. **1475** ..... \$302.00

Coverslips, 20 x 26 x 0.4 mm, Pk/2, Erie. **2026H** ..... \$10.80

### Howard Mold Chamber, Hausser, USA

*Requires Howard Disc - see RETICLES & STAGE MICROMETERS*

Chamber with 2 Coverslips (5080 & 5090). **3820** ..... \$375.00

Coverslip only, 28 x 33 x 0.5 mm (thin). **5080** ..... \$46.80


Coverslip only, 28 x 33 x 1.0 mm (thick). **5090** ..... \$46.80

### Petroff-Hausser Bacteria Counter, Hausser, USA

Chamber with 2 Coverslips (5051 & 1 grid). **3900** ..... \$578.00

Coverslip only, 20 x 26 x 0.2 mm. **5051** ..... \$46.80

**England Field Finder, 1 mm thick, USA. KR-834** ..... \$159.00

**Eyepiece Pointer, wire. 14301** .....  ..... \$2.00




**Eyetube Sleeve, 23.2 mm increased to 30 mm. 14320** ..... \$35.00

# SLIDE & CASSETTE PROCESSING

## Microscope Slides

### Slides, general, 75 x 25 x 1 mm, Bx/72

Plain, ground edges		<u>1 bx</u>	<u>cs/20 bx</u>
Premier®, safety corners. <b>9101</b> .....		\$3.75	72.00
Frosted, beveled edges			
Premier®, sandblasted, safety corners. <b>8205</b> ....		\$5.70	108.00
White End, clipped corners, "Premiere"			
Premier®. <b>9108W-E</b> .....		—	108.00
Superfrost®, "SPECIMEN". <b>4951-001T</b> .....		—	257.00
(12-550-143, 48311-600)			




### Slides, specialty


Charged (silanated), 75 x 25 x 1 mm			
Recommended for adhering tissue sections and cell preparations.			
Diamond™, white end, Globe			
Cs/20 Bx. <b>1358W</b> .....			\$258.00
Superfrost® Plus, white end, Erie, USA			
Cs/20 Bx. <b>4951PLUS-001</b> .....			\$597.00
Preferred slide for Ventana instrumentation. (12-550-15, 48311-703)			



Concavity			
75 x 25 x 1.2 mm, economy		<u>ea</u>	<u>1 bx/12</u>
1-Well, 15 mm dia., 0.5 mm deep. <b>7103</b> .....	\$0.70		6.30
2-Well, 15 mm dia., 0.5 mm deep. <b>7104</b> .....	\$0.95		8.45
76 x 26 x 1.5 mm, Erie		<u>1 bx/12</u>	
1-Well, 18 mm dia., 0.8 mm deep. <b>1518</b> .....			\$43.10
2-Well, 18 mm dia., 0.8 mm deep. <b>1519</b> .....			\$53.75
76 x 26 x 3.2 mm, Erie		<u>1 bx/12</u>	
1-Well, 18 mm dia., 1.75 mm deep. <b>1526</b> .....			\$64.95
2-Well, 18 mm dia., 1.75 mm deep. <b>1527</b> .....			\$72.55
Fused Quartz, NSG Type OZ®, 3" x 1" x 1 mm, USA			
<b>QMS3X1X1</b> , each.....			\$20.00
Fused Silica, Corning 7980, 3" x 1" x 1 mm, USA			
<b>FSMS3X1X1</b> , each.....			\$29.00
3" x 2", economy, 1 gross. <b>6101</b> .....			\$16.75


## Coverslips (Cover Glasses)


### Coverslips, premium, borosilicate glass, 1 oz/pk, Cs/10 oz

#0 Thickness (0.08-0.12 mm)			
Circles, Germany*			
25 mm, ~222 pcs/oz, Cs/10 oz. <b>0-25</b> .....			\$214.00

#1 Thickness (0.13-0.17 mm)			
Squares & Rectangles, Erie, USA		<u>1 oz</u>	<u>cs/10 oz</u>
22 x 22 mm, ~162 pcs/oz. <b>1-2222</b> .....		\$7.45	62.70
22 x 30 mm, ~119 pcs/oz. <b>1-2230</b> .....		—	57.30
22 x 50 mm, ~71 pcs/oz. <b>1-2250</b> .....		—	57.30
24 x 50 mm, ~65 pcs/oz. <b>1-2450</b> .....		\$6.85	57.30
24 x 60 mm, ~54 pcs/oz. <b>1-2460</b> .....		—	57.30

#1 Thickness (0.13-0.17 mm)			
Circles, Germany*		<u>cs/5 oz</u>	<u>cs/10 oz</u>
12 mm, ~663 pcs/oz. <b>1-12</b> .....		\$268.00	—
18 mm, ~296 pcs/oz. <b>1-18</b> .....		\$150.00	—
25 mm. <b>1-25</b> .....		—	\$150.00

#1.5 Thickness (0.16-0.19 mm). Recommended for most users			
Squares & Rectangles, Erie, USA		<u>1 oz</u>	<u>cs/10 oz</u>
18 x 18 mm, ~209 pcs/oz. <b>1.5-1818</b> .....		—	61.40
22 x 22 mm, ~140 pcs/oz. <b>1.5-2222</b> .....		\$7.45	62.70
22 x 40 mm, ~77 pcs/oz. <b>1.5-2240</b> .....		—	57.30
22 x 50 mm, ~62 pcs/oz. <b>1.5-2250</b> .....		—	57.30
24 x 30 mm, ~94 pcs/oz. <b>1.5-2430</b> .....		—	57.30
24 x 40 mm, ~71 pcs/oz. <b>1.5-2440</b> .....		—	57.30
24 x 50 mm, ~56 pcs/oz. <b>1.5-2450</b> .....		\$6.85	57.30
24 x 60 mm, ~47 pcs/oz. <b>1.5-2460</b> .....		—	57.30

#1.5 Thickness (0.16-0.19 mm). Recommended for most users			
Circles, Germany*		<u>cs/5 oz</u>	<u>cs/10 oz</u>
12 mm, ~552 pcs/oz. <b>1.5-12</b> .....		\$223.00	—
18 mm, ~245 pcs/oz. <b>1.5-18</b> .....		\$250.00	—
22 mm, ~164 pcs/oz. <b>1.5-22</b> .....		—	\$223.00
25 mm, ~127 pcs/oz. <b>1.5-25</b> .....		—	\$252.00
Circles, Erie, USA		<u>cs/1 oz</u>	<u>cs/10 oz</u>
40 mm, ~54 pcs/oz. <b>1.5-40</b> .....		\$29.00	244.00
40 mm required for Biotech's FCS2 chamber.			

\*German glass is preferred for cell growth/adhesion

### Coverslips, specialty

Counting Chamber - see Counting Chamber Coverslips in MICROSCOPE Accessories & Parts section

#### Fused Quartz, USA

0.15-0.25 mm Thick, NSG Type OZ®			
1" Square, each. <b>QCS1X1</b> .....			\$19.00
1" Circle, each. <b>QCS1X0</b> .....			\$19.00
0.16 ± 0.01 mm Thick, A1			
25 mm Square, each. <b>R525025</b> .....			\$32.40
25 mm Circle, each. <b>R525000</b> .....			\$23.40

#### Fused Silica, USA

0.15-0.25 mm Thick, Corning 7980			
1" Square, each. <b>FSCS1X1</b> .....			\$21.50
1" Circle, each. <b>FSCS1X0</b> .....			\$21.50
0.16 ± 0.01 mm Thick, S1-UV			
25 mm Square, each. <b>R425025</b> .....			\$32.00
25 mm Circle, each. <b>R425000</b> .....			\$28.80