



134730000




134610000

Test Strips..... Call regarding availability

**Label Tape Dispensers, Bel-Art**

Multi-Rol, 134730000 ..... \$39.30  
 Single-Roll, 134610000 ..... \$87.05

**Thermometers, Digital, Traceable®**



Long Stem Ultra™/  
X-Long Stem Ultra™

4352/4354



4371/  
4378

Lollipop/  
Ultra Lollipop




Refrigerator



Refrigerator/  
Freezer


4127



4374 Dial

Model	Range	Resolution	Accuracy	Probe Length	Unit Size	Replacement Battery	Price
4352	-58 to 302°F/ -50 to 150°C	0.1°	±0.2°C	8"	11 x 1 <sup>3</sup> / <sub>16</sub> " dia.	*357	\$58.10
4354	-58 to 572°F/ -50 to 300°C	0.1°	±0.5°C	11 3/8"	14 3/8 x 1 <sup>3</sup> / <sub>16</sub> " dia.	*357	\$58.75
4371	-58 to 572°F/ -50 to 300°C	0.1°	±1.0°C	8"	10 1/4 x 1 3/4" dia.	*LR1130	\$27.90
4378	-58 to 572°F/ -50 to 300°C	0.1°	±0.4°C	8"	10 1/4 x 1 3/4" dia.	* LR1130	\$63.40
4158	14.2 to 131°F/ -9.9 to 55°C	0.1°	±0.5°C	N/A	2 1/2 x 1 x 3/8"	* LR1130	\$29.55
4127	-58 to 158°F/ -50 to 70°C	1°	±1.0°C	2 1/2"	4 1/4 x 2 3/4 x 3/4"	*AA	\$52.85
4374	-58 to 302°F/ -50 to 150°C	0.1°	±0.3°C	5 1/4"	6 1/2 x 2 1/8" dia.	*357	\$52.20

\*Battery (replacement): 357 / \$4.00; LR1130/ \$2.25. Warranty: 1 year



**Thermometers, Glass, student grade**

Model	Range	Divisions	Length	Features	Price
97-3104	-20 to 230°F/ -30 to 110°C	2°F/1°C	5 3/4"	Red Alcohol/Aluminum V-Back.	\$2.15
97-3110	-10 to 102°C	1°C	6"	Red Alcohol/White Background. Partial immersion.	\$3.15
97-3108	0-230°F/ -20 to 110°C	2°F/1°C	12"	Red Alcohol/White Background. Partial immersion.	\$3.15

