



Tube Coolers, benchtop

0°C, CoolCube™, Benchmark

Microtube and PCR Plate. **R1000** \$104.00
 Each cube holds (36) 1.5/2.0 mL tubes, or (96) 0.2 mL tubes, or (12) 0.2 mL strips, or (1) 96-well PCR plate and (10) 1.5/2.0 mL tubes. Maintains 0°C for 4 hours.

0°C, Cryo-Safe™, Bel-Art

Mini, holds (12) 0.5/1.5 mL tubes. **188460000** \$124.00
 Lasts 2.5 hours.

Maxi, holds (32) 0.5/1.5 mL tubes. **188462000** \$219.00
 Lasts 3 hours.

-20°C, Cryo-Safe™, Bel-Art

Mini, holds (12) 0.5/1.5/2.0 mL tubes. **188460020** \$95.15
 Lasts 1 hour.

Maxi, holds (32) 0.5/1.5/2.0 mL tubes. **188462020** \$189.00
 Lasts 1.5 hours.

Dividers only, fit standard 2 & 3" boxes

- 100-Cell (10 x 10), Pk/12. **HS2860E** \$25.45
 Holds 1.5-2.0 mL tubes and 1.2-2.0 mL cryovials.
- 81-Cell (9 x 9), Pk/12. **HS2860H** \$25.45
 Holds 13 mm vials and 4.0-5.0 mL cryovials.
- 64-Cell (8 x 8), Pk/12. **HS2860F** \$24.15
 Holds 14 mm vials.
- 25-Cell (5 x 5), 25 mm cell, Pk/12. **D-25** \$21.60
- 16-Cell (4 x 4), 31 mm cell, Pk/12. **D-16** \$21.60
- 5" Box Only, 5¼ x 5¼ x 5", Pk/10. **HS2860M** \$43.90
- 5" Box w/dividers
- 36-Cell (6 x 6), 15 mL tubes, Pk/10. **HS2860M/T** \$65.10
- 16-Cell (4 x 4), 50 mL tubes, Pk/10. **HS2860M/P** \$62.00



Polycarbonate boxes store vials from -196°C to 121°C.

Arctic Squares®, polycarbonate, Heathrow

- 25-Well (5 x 5), friction-fit lid, Pk/8. **HS23325A** \$52.05
 Holds vials up to 2 mL; red. 7.6 x 7.6 x 5.3 cm
- 81-Well (9 x 9), hinged lid, Pk/4. **HS23381B** \$42.55
 Holds vials up to 2 mL; mixed. 13.3 x 13.3 x 5.3 cm
- Cryo-Safe™, polycarbonate, 13.3 x 13.3 x 9.5 cm, Bel-Art
- 81-Well (9 x 9), Pk/4. **188490002** \$56.75
 Holds 3.0/5.0 mL vials.



HS2860A/H

Tube Storage Boxes (cryogenic)

Cardboard, vented

2" Box Only, 5¼ x 5¼ x 2", Pk/12. **HS2860A** \$26.60
 2" box w/dividers

White

- 100-Cell (10 x 10), Pk/12. **HS2860A/E** \$46.90
- ★ 81-Cell (9 x 9), Pk/12. **HS2860A/H** \$46.90
- 64-Cell (8 x 8), Pk/12. **HS2860A/F** \$46.25

Mixed colors (four each of blue, red & yellow)

- 100-Cell (10 x 10), Pk/12. **HS2860C/E** \$47.95
- 81-Cell (9 x 9), Pk/12. **HS2860C/H** \$47.95
- 64-Cell (8 x 8), Pk/12. **HS2860C/F** \$47.25


2" Mini box w/dividers, 3 x 3 x 2"

25-Cell (5 x 5), Pk/12. **MB-2/MD-25** \$50.40
 Holds 13 mm vials and 4.0-5.0 mL cryovials.

3" Box Only, 5¼ x 5¼ x 3", Pk/12. **HS2860B** \$28.80
 3" box w/dividers

- 100-Cell (10 x 10), Pk/12. **HS2860B/E** \$48.30
- 81-Cell (9 x 9), Pk/12. **HS2860B/H** \$48.30
- 64-Cell (8 x 8), Pk/12. **HS2860B/F** \$47.60

HS120035/7



HS120044



Polypropylene, Heathrow

- 81-Well (9 x 9), 13 x 13 x 4.7 cm
 Holds 1.5-2.0 mL tubes and cryovials.
 Features imprinted coordinates on a colored base and a translucent, friction-fit lid with frosted panels.
- Blue base, each. **HS120035** (HS23255B) \$6.95
- Pink base, each. **HS120037** (HS23255D) \$6.95
- Black base & lid, Pk/5. **HS120041** \$32.30
- 100-Well (10 x 10), mixed, Pk/5. **HS120044** (HS1065) \$33.65
 Holds 1.5-2.0 mL tubes and cryovials.
 Hinged lid has molded grid lines; base has imprinted and raised coordinates for each well. 14.1 x 15.1 x 5.7 cm.

Holds 1.5-2.0 mL tubes in a two-piece rack.



Polystyrene Foam (styrofoam), Heathrow, USA

- 50-Well, Pk/10. **HS4283P** \$40.10
- 100-Well, Pk/10. **HS4282** \$50.45

Continued