


Microscope Stage Laminate, USA. 16251 \$8.00
6.25 x 7.75" sheet provided with a #10 scalpel.
White Teflon® (PTFE) on self-adhesive vinyl is 0.010" thick.



Nosepiece Dust Plug, RMS thread (20.32 mm Dia.). 15420..... \$0.50



Mark slides, dishes, flasks, and chambers with a 1.8 mm ring for later referral. Works on upright and inverted microscopes.

Object Marker, RMS thread (20.32 mm Dia.)
Black. **14361** \$172.00
Refill Ink, black, 10 mL, Shachihata. **22112** \$7.00
Adapters
RMS/45 to CFI60, 15 mm long. **14111** \$43.00
RMS ID to 25 mm OD Thread, brass. **14115** \$52.00

Microscopes with Nikon CFI60 infinity-corrected (∞) objectives require a 14111 adapter.
Microscopes with Leica 25 mm thread/45 mm parfocal length objectives require a 14115 adapter.
Microscopes with 160 mm tube length objectives do not require an adapter.



Objectives, 160 mm tube length, DIN, Jenco, Meiji, National
Achromats, Motic/National
4X, N.A. 0.10. **704-160** \$32.00
10X, N.A. 0.25. **710-160** \$33.00
40X, N.A. 0.65. **740-160** \$42.00
100X Oil, N.A. 1.25. **799-160** \$69.00
Super High Contrast Achromats, National
4X, N.A. 0.10. **704-160ASC** \$28.00
10X, N.A. 0.25. **710-160ASC** \$36.00
40X, N.A. 0.65. **740-160ASC** \$66.00
100X Oil, N.A. 1.25. **799-160ASC** \$99.00
Semi-Plan Achromats, National (unless noted)
2X, N.A. 0.05, Accu-Scope. **3193SP-00** \$70.00
4X, N.A. 0.10. **704-160SP** \$32.00
10X, N.A. 0.25. **710-160SP** \$33.00
40X, N.A. 0.65. **740-160SP** \$42.00
50X Oil, N.A. 0.95, Meiji. **MA694** \$295.00
60X, N.A. 0.85. **MA645** \$315.00
100X Oil, N.A. 1.25. **799-160SP** \$69.00

Continued

Plan Achromats, National (unless noted)
2.5X, N.A. 0.08, Meiji. **MA432** \$285.00
4X, N.A. 0.10. **704-160P** \$33.00
10X, N.A. 0.25. **710-160P** \$55.00
20X, N.A. 0.40. **720-160P** \$125.00
40X, N.A. 0.65. **740-160P** \$104.00
50X Oil, N.A. 0.95, Accu-Scope. **3190PL-00** \$315.00
60X, N.A. 0.85. **760-160P** \$199.00
100X Oil, N.A. 1.25. **799-160P** \$168.00

Protect spare microscope objectives.



Objective Case, DIN, polystyrene. 15601 \$1.00




Extender tubes do not contain lenses.

Objective Extender Tubes, RMS thread, brass, USA
10 mm Long. **14105** \$44.20
18 mm Long. **14106** \$42.75

Objective Shims, fit RMS thread (20.32 mm Dia.), USA
Parfocalize objectives with these stainless steel rings.
25 µm (0.001") Thick, Pk/10. **10032061** \$16.15
50 µm (0.002") Thick, Pk/10. **10032062** \$16.15
75 µm (0.003") Thick, Pk/10. **10032063** \$16.15
100 µm (0.004") Thick, Pk/10. **10032064** \$16.15
125 µm (0.005") Thick, Pk/10. **10032065** \$16.15
Assortment Pack, pk/10. **100320630** \$19.80
Includes two each of the following thicknesses: 25, 50, 75, 100, and 125 µm.

Highly efficient, linear polarizing film has excellent color reproducibility. Transmission is 38% for a single polarizer, while extinction is better than 99.99% for crossed polars. Film is 0.76 mm thick.



Polarizing Film, 50 x 50 mm Pair w/Plastic Case, USA. 15751 ... \$12.00