


Reticles & Stage Micrometers



Reticles, eyepiece, 1.5 mm thick, Klarmann Rulings, USA

*Specify reticle diameter using 2-digit suffix,
for example, use KR-207-21 for 21 mm diameter reticle

★ 5 mm/100 Parts. KR-209-_*	\$77.00
10 mm/100 Parts. KR-207-_*	\$77.00
Crossline 90°. KR-301-_*	\$77.00
Crossline 90° with 10 mm/100 Parts Scale. KR-257-_*	\$77.00
Howard Disc. KR-832-_*	\$77.00
Indexed Square, 100 sqs, each 1 x 1 mm. KR-406B-_*	\$77.00
1-10 on x-axis; A-J on y-axis	
Miller Disc	
Corner Square. KR-802-_*	\$77.00
Central Square. KR-802A-_*	\$77.00
Walton Beckett, G-22. KR-827-1-_*	\$189.00
Whipple Disc. KR-413-_*	\$77.00




Stage Micrometers, Transmitted Light

1 mm/100 Parts Scale (0.01 mm), Meiji. MA285	\$105.00
1 mm-0.01 mm/div on 1 mm thick slide, protected with microcoverslip. Preferred for compound microscopes.	
10 mm/100 Parts Scale (0.1 mm), Meiji. MA187	\$105.00
10 mm-0.1 mm/div on 1 mm thick slide, protected with 1 mm thick coverslip. Designed for stereo microscopes.	
25 mm Scale, USA. KR-812	\$173.00
25mm-1 mm/div, 5 mm-0.1 mm/div, 2 mm-0.02 mm/div on 3 mm thick slide. Preferred at magnifications below 400X. Ideal for large field stereo microscopes.	
0.2 mm Square/400 Squares, Meiji. MA287	\$105.00
20 vertical lines, 1 mm long and spaced 0.01 mm apart, intersect 20 similar horizontal lines to create a central 0.2 mm square consisting of 400 squares, each 0.01 mm. 1 mm thick slide is protected with microcoverslip. Useful for x-y calibration.	
0.04"/40 Parts Scale (0.001"), Meiji. MA286	\$105.00
0.04"-0.001"/div on 1 mm thick slide, protected with microcoverslip.	

Stage Micrometer, Reflected Light

1 mm/100 Parts Scale (0.01 mm), Meiji. MA292	\$210.00
1 mm-0.01 mm/div on 3 mm thick metal slide	



Deter microscope theft. Each kit includes 3 steel mini plates, super bonding adhesive, steel padlock, and 5' steel cable. Additional cables and plates are available for expanding the basic kit.

Security Kit, Super Bond, USA. 95267	\$39.00
Extension Cable, 5'. 5859	\$11.00
Security Plate, 2 x 4". 5862S	\$25.00

Cameras

Analog (Video)




Color, 480 TVL(H), Japan

1/3" IT CCD, 2:1 interlace, CS/C-mount. STC-630PWJ	\$490.00
1/2" IT CCD, 2:1 interlace, CS/C-mount. STC-620PWJ	\$576.00
<ul style="list-style-type: none"> • NTSC (Composite) output. • Y/C (S-VHS) output. • 0.11 lux minimum illumination (AGC on) • Push to set white balance • 3-way image flip • Includes power supply and 12' BNC cable • 51(L) x 51(W) x 60.5(H) mm, 190 g 	

Cameras require a video coupler for microscope attachment - see VIDEO COUPLERS

Digital



Digital Camera, Sentech. STC-MC83USB-AS	\$415.00
<ul style="list-style-type: none"> • 0.8 megapixel USB 2.0 output [30 fps (live) at 1024 x 768 pixels] • 1/3" IT color CCD, progressive scan, XGA • 21 lux minimum illumination • Sentech Viewing software • CS/C-mount • 28(L) x 28(W) x 37(H) mm 	

Continued



- Digital Camera, National. Moticam 2** \$332.00
- 2.0 megapixel (1600 x 1200 pixels) USB 2.0 output with card
 - 1/3" CMOS (Complementary Metal Oxide Semiconductor) chip
 - Adapters for attaching camera to microscope include a C-mount ring and couplers to fit 30-38 mm OD eyepieces
 - 12 mm focusable lens and viewing tube for macro viewing
 - 4-dot calibration slide
 - Motic Images Plus 2.0 imaging software

Digital SLR



- Digital SLR Camera, 18 megapixel, Canon. EOS Rebel T5i** \$779.00
 Includes Canon 18-55mm EF-S IS II lens, SD-16 GB card and shoulder bag.
 AC Adapter Kit. **ACK-E8** \$104.00


Accessories



- Adapter, CS-Mount to C-Mount, 5 mm spacer** \$15.00
Do not use this adapter ring on an existing C-mount camera, as it will negatively impact the quality and parfocality of the imaging system.

Cables

- BNC, male-to-male, 75Ω, RG59/U
- 6' \$8.80
 - 12' \$10.60
- S-VHS, 4-pin male to 4-pin male
- 6' \$10.00
 - 12' \$14.00



Cable Connectors

- BNC Female to RCA Male. **8110** \$5.65
- BNC Female to BNC Female. **25-7430** \$4.20
- BNC Tee, female-male-female. **27-8140** \$9.00
- RCA Female to BNC Male. **7510** \$4.50

- Cable Release. AR-3** \$24.95



Couplers, Diagnostic Instruments, USA

- 1X C-Mount for 38 mm ISO Photoport. **D10NLC** \$135.00
Fits newer Nikon and Leica microscopes
- Nikon SLR Lens to C-Mount TV. **NIKL-CMTC** \$116.00

Couplers, National

- 1X C-Mount for 23.2 mm ID Eyetube. **930-023** \$29.00
Lens free - this coupler is not parfocal

Magnifiers




ELP-174.10



ELP-750.01

Loupes

- Eye, plastic, 10X at 1", black. **ELP-174.10** \$2.60
- Triplet, 10X at 1", 18 mm, black. **ELP-750.01** \$15.95
Three optical glass lenses bonded together create a single achromatic (color-corrected) lens. Leather case provided.



Large diameter (4 1/8") glass lens is supported on a sturdy plastic, foldout stand for hands-free viewing. Base is graduated in both cm and inches for measuring objects.

- Magnifier, Folding (Linen Tester), 3X. ELP-800.03** \$9.95



A 4X magnifying lens is attached to fine-tipped, stainless steel forceps. Forceps double as a lens support, freeing viewer's hands. Folds for storage.

- Magnifier, Forceps. M15** \$4.00

Continued