



- Digital Camera, National. Moticam 2** ..... \$332.00
- 2.0 megapixel (1600 x 1200 pixels) USB 2.0 output with card
  - 1/3" CMOS (Complementary Metal Oxide Semiconductor) chip
  - Adapters for attaching camera to microscope include a C-mount ring and couplers to fit 30-38 mm OD eyepieces
  - 12 mm focusable lens and viewing tube for macro viewing
  - 4-dot calibration slide
  - Motic Images Plus 2.0 imaging software

**Digital SLR**



- Digital SLR Camera, 18 megapixel, Canon. EOS Rebel T5i** .... \$779.00  
 Includes Canon 18-55mm EF-S IS II lens, SD-16 GB card and shoulder bag.  
 AC Adapter Kit. **ACK-E8** ..... \$104.00


**Accessories**



- Adapter, CS-Mount to C-Mount, 5 mm spacer** ..... \$15.00  
*Do not use this adapter ring on an existing C-mount camera, as it will negatively impact the quality and parfocality of the imaging system.*

**Cables**

- BNC, male-to-male, 75Ω, RG59/U
- 6' ..... \$8.80
  - 12' ..... \$10.60
- S-VHS, 4-pin male to 4-pin male
- 6' ..... \$10.00
  - 12' ..... \$14.00



**Cable Connectors**

- BNC Female to RCA Male. **8110** ..... \$5.65
- BNC Female to BNC Female. **25-7430** ..... \$4.20
- BNC Tee, female-male-female. **27-8140** ..... \$9.00
- RCA Female to BNC Male. **7510** ..... \$4.50

- Cable Release. AR-3** ..... \$24.95



**Couplers, Diagnostic Instruments, USA**

- 1X C-Mount for 38 mm ISO Photoport. **D10NLC** ..... \$135.00  
*Fits newer Nikon and Leica microscopes*
- Nikon SLR Lens to C-Mount TV. **NIKL-CMTC** ..... \$116.00

**Couplers, National**

- 1X C-Mount for 23.2 mm ID Eyetube. **930-023** ..... \$29.00  
*Lens free - this coupler is not parfocal*

**Magnifiers**



ELP-174.10



ELP-750.01

**Loupes**

- Eye, plastic, 10X at 1", black. **ELP-174.10** ..... \$2.60
- Triplet, 10X at 1", 18 mm, black. **ELP-750.01** ..... \$15.95  
*Three optical glass lenses bonded together create a single achromatic (color-corrected) lens. Leather case provided.*



Large diameter (4 1/8") glass lens is supported on a sturdy plastic, foldout stand for hands-free viewing. Base is graduated in both cm and inches for measuring objects.

- Magnifier, Folding (Linen Tester), 3X. ELP-800.03** ..... \$9.95



A 4X magnifying lens is attached to fine-tipped, stainless steel forceps. Forceps double as a lens support, freeing viewer's hands. Folds for storage.

- Magnifier, Forceps. M15** ..... \$4.00

*Continued*



Lenses are lightweight, water-clear, and shatterproof. Molded handles protect the lens and offer comfort and balance.

**Magnifier, Handle, Ultra-Optix**


- 2" Lens, 3X. **SV-2P**..... \$5.45
- 3" Lens, 2X with 5X bifocal inset. **SV-3P**..... \$7.45
- 4" Lens, 2X with 4X bifocal inset. **SV-4P**..... \$12.55



Provides precise binocular viewing. Visor tilts out of the way when optical glass magnifiers are not needed. Adjustable headband is leather lined. Other magnifications and replacement lenses are available.

**Magnifier, Headband, Optivisor® , USA**

- 1.75X at 14", DA3. **ELP-503.00**..... \$48.95
- 2.5X at 8", DA5. **ELP-505.00**..... \$48.95
- 2.75X at 6", DA7. **ELP-507.00**..... \$50.95
- Light Attachment, Quasar™ LED. **ELP-558.00**..... \$29.95



View small details in low light with this compact, illuminated magnifier. Protective case keeps the lens scratch-free and serves as a handle when lens is extended.

**Magnifier, Lighted, 2 x 2" lens, 2X. SV-BPL**..... \$7.45



Confirm patient identification and preview the slide before examining microscopically. This macro inspection magnifier mounts onto the microscope's arm, where the lens angle can be adjusted using the double-ball joint arm. Large 3 x 1¾" lens, 2.25X, permits viewing of entire 3 x 1" glass slide.

**Magnifier, Magnistitch™, USA. 15401**..... \$12.50