


SLIDE & CASSETTE PROCESSING

Microscope Slides

Slides, general, 75 x 25 x 1 mm, Bx/72

Plain, ground edges		<u>1 bx</u>	<u>cs/20 bx</u>
Premier®, safety corners. 9101		\$3.75	72.00
Frosted, beveled edges			
Premier®, sandblasted, safety corners. 8205		\$5.70	108.00
White End, clipped corners, "Premiere"			
Premier®. 9108W-E		—	108.00
Superfrost®, "SPECIMEN". 4951-001T		—	257.00
(12-550-143, 48311-600)			





Slides, specialty



Charged (silanated), 75 x 25 x 1 mm			
Recommended for adhering tissue sections and cell preparations.			
Diamond™, white end, Globe			
Cs/20 Bx. 1358W		\$258.00	
Superfrost® Plus, white end, Erie, USA			
Cs/20 Bx. 4951PLUS-001		\$597.00	
Preferred slide for Ventana instrumentation. (12-550-15, 48311-703)			
Concavity			
75 x 25 x 1.2 mm, economy		<u>ea</u>	<u>1 bx/12</u>
1-Well, 15 mm dia., 0.5 mm deep. 7103		\$0.70	6.30
2-Well, 15 mm dia., 0.5 mm deep. 7104		\$0.95	8.45
76 x 26 x 1.5 mm, Erie			<u>1 bx/12</u>
1-Well, 18 mm dia., 0.8 mm deep. 1518		\$43.10	
2-Well, 18 mm dia., 0.8 mm deep. 1519		\$53.75	
76 x 26 x 3.2 mm, Erie			<u>1 bx/12</u>
1-Well, 18 mm dia., 1.75 mm deep. 1526		\$64.95	
2-Well, 18 mm dia., 1.75 mm deep. 1527		\$72.55	
Fused Quartz, NSG Type OZ®, 3" x 1" x 1 mm, USA			
QMS3X1X1 , each.....		\$20.00	
Fused Silica, Corning 7980, 3" x 1" x 1 mm, USA			
FSMS3X1X1 , each.....		\$29.00	
3" x 2", economy, 1 gross. 6101		\$16.75	


Coverslips (Cover Glasses)


Coverslips, premium, borosilicate glass, 1 oz/pk, Cs/10 oz

#0 Thickness (0.08-0.12 mm)			
Circles, Germany*			
25 mm, ~222 pcs/oz, Cs/10 oz. 0-25		\$214.00	

#1 Thickness (0.13-0.17 mm)			
Squares & Rectangles, Erie, USA		<u>1 oz</u>	<u>cs/10 oz</u>
22 x 22 mm, ~162 pcs/oz. 1-2222		\$7.45	62.70
22 x 30 mm, ~119 pcs/oz. 1-2230		—	57.30
22 x 50 mm, ~71 pcs/oz. 1-2250		—	57.30
24 x 50 mm, ~65 pcs/oz. 1-2450		\$6.85	57.30
24 x 60 mm, ~54 pcs/oz. 1-2460		—	57.30

#1 Thickness (0.13-0.17 mm)			
Circles, Germany*		<u>cs/5 oz</u>	<u>cs/10 oz</u>
12 mm, ~663 pcs/oz. 1-12		\$268.00	—
18 mm, ~296 pcs/oz. 1-18		\$150.00	—
25 mm. 1-25		—	\$150.00

#1.5 Thickness (0.16-0.19 mm). Recommended for most users			
Squares & Rectangles, Erie, USA		<u>1 oz</u>	<u>cs/10 oz</u>
18 x 18 mm, ~209 pcs/oz. 1.5-1818		—	61.40
22 x 22 mm, ~140 pcs/oz. 1.5-2222		\$7.45	62.70
22 x 40 mm, ~77 pcs/oz. 1.5-2240		—	57.30
22 x 50 mm, ~62 pcs/oz. 1.5-2250		—	57.30
24 x 30 mm, ~94 pcs/oz. 1.5-2430		—	57.30
24 x 40 mm, ~71 pcs/oz. 1.5-2440		—	57.30
24 x 50 mm, ~56 pcs/oz. 1.5-2450		\$6.85	57.30
24 x 60 mm, ~47 pcs/oz. 1.5-2460		—	57.30

#1.5 Thickness (0.16-0.19 mm). Recommended for most users			
Circles, Germany*		<u>cs/5 oz</u>	<u>cs/10 oz</u>
12 mm, ~552 pcs/oz. 1.5-12		\$223.00	—
18 mm, ~245 pcs/oz. 1.5-18		\$250.00	—
22 mm, ~164 pcs/oz. 1.5-22		—	\$223.00
25 mm, ~127 pcs/oz. 1.5-25		—	\$252.00
Circles, Erie, USA		<u>cs/1 oz</u>	<u>cs/10 oz</u>
40 mm, ~54 pcs/oz. 1.5-40		\$29.00	244.00
40 mm required for Biotech's FCS2 chamber.			

*German glass is preferred for cell growth/adhesion

Coverslips, specialty

Counting Chamber - see Counting Chamber Coverslips in MICROSCOPE Accessories & Parts section

Fused Quartz, USA

0.15-0.25 mm Thick, NSG Type OZ®

1" Square, each. QCS1X1	\$19.00
1" Circle, each. QCS1X0	\$19.00

0.16 ± 0.01 mm Thick, A1

25 mm Square, each. R525025	\$32.40
25 mm Circle, each. R525000	\$23.40

Fused Silica, USA


0.15-0.25 mm Thick, Corning 7980

1" Square, each. FSCS1X1	\$21.50
1" Circle, each. FSCS1X0	\$21.50

0.16 ± 0.01 mm Thick, S1-UV

25 mm Square, each. R425025	\$32.00
25 mm Circle, each. R425000	\$28.80

Accessories



Bibulous Paper, 10 x 15 cm, Bk/50 sheets. **BIB-46-50** bk cs/48 bk
 \$3.65 169.00

Forceps

Spade Tips



- Angled (Kuehne)
 - Stainless Steel, 11 cm (4¼"), economy. **3270**..... \$4.15
 - Teflon® Coated, 11 cm (4¼"), Bel-Art. **379420000**..... \$20.90
 - Straight, stainless steel, economy. **57.926**..... \$4.90

Square Point Tips, straight



- Stainless Steel, 12 cm (4¾"), economy. **1244**..... \$2.25
- Teflon® Coated, 11 cm (4¼"), Bel-Art. **379430000**..... \$22.90



Minimize contamination during slide & coverslip handling. Will not scratch; corrosion and heat resistant.

Hemostat, Polypropylene, 12.5 cm. **HEM-300.00** \$1.50

Cytology Funnels

Single,
Metal Clips
Required



Double,
All Plastic



White filter cards will accept up to 0.5 mL of sample, while tan filters will accept volumes up to 0.4 mL.

Cytology Funnels, Simport CytoSep™

Shandon (Thermo) Cytospin® Cytocentrifuges

Single. *Metal slide clips required*

White Filter & Cap, Cs/500. **M964-10FW** #1 \$382.00
 (59 91 040, 110 25 48)

Tan Filter & Cap, Cs/500. **M964-10FT** \$382.00

Single. *All plastic - metal slide clips not required*

White Filter & Cap, Cs/480. **M965-10FW** \$503.00
 (A787 100003)

Tan Filter & Cap, Cs/480. **M965-10FT** \$503.00

Double. *All plastic - metal slide clips not required*

White Filter & Cap, Cs/480. **M965-20FW** \$626.00
 (A787 100005)

Immersion Oils & Mounting Media



Cargille immersion oils A, B, and NVH can be mixed with each other for custom viscosities. USA

Immersion Oil

Type A, low viscosity (150 cSt). *Short W.D. optics*

- 7.5 mL in LDPE dispensing bottle. **18101** \$5.35
- 20 mL in HDPE dispensing bottle. **18102** \$10.00
- 120 mL in amber bottle. **18103** \$22.45
- 480 mL in amber bottle. **18104** \$49.60

Type B, high viscosity (1250 cSt). *Long W.D. optics*

- 7.5 mL in LDPE dispensing bottle. **18105** \$5.35
- 20 mL in HDPE dispensing bottle. **18106** \$10.00
- 120 mL in amber bottle. **18107** \$22.45
- 480 mL in amber bottle. **18108** \$49.60


Type NVH, very high viscosity (21,000 cSt)

- 480 mL in amber bottle. **18113** \$57.75

Type FF, fluorescence free (170 cSt)

- 7.5 mL in LDPE dispensing bottle. **18121** \$7.35
- 480 mL in amber bottle. **18124** \$81.75

A long nozzle, textured sides, and slim profile make this chemically inert, high-density polyethylene (HDPE) container ideal for dispensing immersion oil.



Bottle, Dispensing, HDPE, 20 mL, each. **18145** \$6.90



Store immersion oil or mounting media in this borosilicate glass bottle with an applicator rod and ground glass cap.

Bottle, Balsam, glass, 45 mL (empty), USA. **616**..... \$17.30

K-750

Suitable for fluorescence microscopy.



K-700



Takes only 10 minutes to dry.


Mounting Media

Mount-Quick, Kiyota

- 30 mL in Dropper Bottle, aqueous-based. **K-750**..... \$22.40
- 30 mL in Tube, solvent-based, 30 cps. **K-700**..... \$15.40

Slide Storage

PLASTIC SLIDE BOXES



Slide Boxes

- 25-Slide, Holds 75 x 50 mm (3 x 2") Slides
White, ABS plastic, China. **97-0625** \$2.85
squeeze-release cover, cork bottom
- 25-Slide, holds 75 x 25 mm (3 x 1") slides
Squeeze-Release cover, cork bottom, 9.5 x 8.5 x 3 cm
- Blue, ABS plastic, China. **97-0026** \$1.95
- Green, ABS plastic, China. **97-0027** \$1.95
- #1** Red, ABS plastic, China. **97-0028** \$1.95
- White, ABS plastic, China. **97-0025** \$1.95
- One-Piece construction, beige, USA. **PP-25** \$3.40
Resistant to freezer damage and light penetration.
UV protected polypropylene. Dimensions: 12 x 9.5 x 4 cm
Friction-Fit cover, cork bottom, 14.1 x 8.8 x 3.5 cm
White, ABS plastic, China. **HS15989E** \$5.90
- 50-Slide, hinged cover, cork bottom, holds 75 x 25 mm slides
Dimensions: 20 x 9.5 x 3.5 cm
- Blue, ABS plastic, China. **97-0056** \$5.55
- Red, ABS plastic, China. **97-0058** \$5.55
- 100-Slide, holds 75 x 25 mm (3 x 1") slides
One-Piece construction, Unlined
- Blue, polypropylene, China. **HS15991A** \$5.10
- Steel-Hinged cover, Cork bottom ea cs/48
- Gray, polystyrene, USA. **PS-100** \$10.55 452.00
PS-100 high-impact slide box offers superior resistance
to freezer distortion and light penetration. UV protected.
- Blue, ABS plastic, China. **97-0126** \$7.25
- Green, ABS plastic, China. **97-0127** \$7.25
- #1** Red, ABS plastic, China. **97-0128** \$7.25
- White, ABS plastic, China. **97-0125** \$7.25
- Yellow, ABS plastic, China. **HS15994D** \$7.25
- Steel-Hinged cover, Foam bottom
- White, ABS plastic, China. **HS15994R** \$6.95
- Wood Frame Covered w/Black Paper, USA. **44321** \$12.30
A plastic, numbered insert holds slides.

Steel Slide Cabinet



Slide Cabinets

Steel, Paragon, USA

*Specify Gray, Green or Tan for cabinet and/or base

6-Drawer, stackable

Holds 4500 75 x 25 slides in 6 divided drawers.

Each cabinet comes with a sample of 3 compressor follow blocks and 100 memorandum markers.

Gray. **C101Gray** \$132.00

Green. **C101Green** \$132.00

★ Tan. C101Tan \$132.00

12 Cabinet System (includes base). **C112B** * \$1,588.00

Base, raises cabinets 4 3/4" off floor. **C103** * \$54.00

Memorandum Markers, paperboard, Pk/100. **C108-M** .. \$7.40

Compressor Follow Block, rubber, each. **C109** \$1.40

Oak with Aluminum Trays, Eberbach, USA

Holds 75 x 25 mm slides

500-Slide Capacity. **2800** \$549.00

1000-Slide Capacity. **2820** \$1,033.00

2000-Slide Capacity. **2830** \$1,915.00

3000-Slide Capacity. **2840** \$2,570.00


Slide Tray only, Aluminum, holds 20 slides. **2870** \$19.60

Holds 480 75 x 25 mm slides, or 240
75 x 51 mm slides in indexed compartments.
Dimensions: 24.1 x 15.2 x 8.3 cm.



Slide Files

Fiberboard, hinged lid, Cargille. **73010** \$34.75



Slide Sleeves, 10 gauge vinyl

4 x 5 1/2" Card Holder Plus 4 Slide Pockets. **SL-3-264** \$0.45

File Box, 4 x 6 x 11 1/2", fiberboard. **40589** \$16.80

5 x 8" Card Holder Plus 5 Slide Pockets. **SL-3-108** \$0.55

File Box, 6 x 9 x 11 1/2", fiberboard. **40591** \$26.25

8 1/2 x 11" Page Holds 18 Slides. **SL-3-182** \$0.75

20-Slide
Cardboard
Folder



10-Slide/
10-Block
Consulting
Folder



Slide Folders

- 10-Slide/10-Block (consulting), each. **CSF-101** \$24.75
 12-Slide, holds 75 x 50 mm (3 x 2") slides
 Cardboard, black tape, USA. **SF12-BLACK** \$7.20
 20-Slide, holds 75 x 25 mm (3 x 1") slides
 Cardboard, USA
 Folders made in the USA resist heat-induced warping.
 Black Tape. **SF20-BLACK** \$6.25
 Red Tape. **SF20-RED** \$6.25
 Green Tape. **SF20-GREEN** \$6.25
 Navy Blue Tape. **SF20-BLUE** \$6.25
 Cardboard, Black Tape, China. **HS9902** \$4.70
 (12-587-10)
 Plastic, gray tray w/2 clear doors. **M750-20GY** \$9.20
 Transparent covers hold visible slides securely in place,
 but can fold flat against back of tray when opened.
 Slides are slightly elevated for easy retrieval.



12-Slide



52-Slide




91-Slide

Slide Holders, perfect for slide drying or staining

- 12-Slide (horizontally positioned), polystyrene. **445570000** \$7.80
 52-Slide (approx.), 12 x 4", HDPE, USA. **L500-C** \$27.20
 91-Slide (approx.), 12 x 7", HDPE, USA. **L500-D** \$36.35



2-Slide



5-Slide

Slide Mailers, Plastic

- 2-Slide, Cs/1000. **880041** \$116.00
 5-Slide, end opening, Cs/200. USA. **62706** \$87.30
 62706 mailers can also be used as staining jars.
 (1258717B)

Aluminum



Plastic



Slide Trays

- Aluminum, Eberbach, USA. **2870** \$19.60
 Plastic, white, 20-slide, Bel-Art. **445810000** \$12.30
 Wood, 50-slide \$42.00

Slide & Cassette Pens, Markers & Scribes

Write a 0.4 mm stroke on glass slides with this permanent marker.
 Pilot pens are ideal for cytologists and pathologists.



Pens

- | | | |
|-------------------------------------|-----------|--------------|
| Pilot®, permanent, extra fine point | <u>ea</u> | <u>bx/12</u> |
| ★ Black. 44102 | \$1.70 | 17.45 |
| Blue. 44103 | \$1.70 | 17.45 |
| Green. 44105 | \$1.70 | 17.45 |
| Red. 44104 | \$1.70 | 17.45 |



Confine immunohistochemical reagents within a light blue, hydrophobic barrier, or separate multiple specimens on a slide/dish. Resists alcohol and acetone, but removes with xylene.

- Super^{HT} PAP (Peroxidase Anti-Peroxidase), hydrophobic
 Regular Tip (4 mm tapered nib). **K-600**..... \$59.00
 (71310, 195505)
 ★ Fine Tip (2.5 mm tapered nib). **K-650** \$44.00
 (71312, 195506)

Label any glossy surface with this water-resistant wax pencil. Wipes off easily with a cloth. Simply pull the string to remove wrapping and expose more tip.



Markers, USA

- | | | |
|---|-----------|--------------|
| | <u>ea</u> | <u>bx/12</u> |
| China Marker (Wax Pencil). 02089 | \$1.10 | 11.50 |

Ink resists xylene, alcohol, acetone and formalin.
 Remove with acetic acid, chloroform, or ethanol/xylene mixture.



- | | | |
|---|-----------|--------------|
| Securline® Marker II®, permanent, PDC | <u>ea</u> | <u>bx/10</u> |
| Black, 0.6 mm tip. 1450-20-PDC | \$2.50 | 20.30 |
| (14-905-30, 72107) | | |

Fine Tip



Ultra-Fine Tip



Sharpie Retractable, permanent

- | | | |
|-----------------------------|-----------|--------------|
| Fine Tip | <u>ea</u> | <u>bx/12</u> |
| Black. 36701 | — | 29.65 |
| Blue. 36703 | \$2.85 | 29.65 |
| Red. 36702 | \$2.85 | 29.65 |
| Ultra-Fine Tip | <u>ea</u> | <u>bx/12</u> |
| Black. 1735790 | \$2.85 | 29.65 |
| Blue. 1735792 | \$2.85 | 29.65 |
| Red. 1735791 | \$2.85 | 29.65 |

Continued



Scribes, retractable

Carbide-Tipped (resharpenable). RSC	\$15.30
Carbide Tip (replacement). RCT-1	\$6.20
Diamond-Tipped. RSD	\$21.85
Diamond Tip (replacement). RDT-1	\$12.90

Evenly warm your slides from room temperature to 70°C. Black anodized surface provides excellent contrast.



Slide Warmers


Small, 26 x 18 cm (10 x 7"). XH-2002	\$237.00
Large, 63 x 20 cm (25 x 8"). XH-2001	\$318.00

Staining

Glass staining jar with ground glass lid holds 8 coverslips, 18-22 mm, back-to-back. Jar is 45 mm wide and tall.



Coverslip Staining Jar, each. **1257200011**..... \$7.40




Staining/Moisture Chambers

10-Slide, polystyrene, USA	<u>ea</u>	<u>pk/3</u>
Clear. 240-9000-010	\$50.10	141.00
Black. 240-9020-Z10	\$50.10	141.00



20-Slide, Simport

Clear Lid (PETG), black base (ABS), each. M920-1	\$158.00
Black Lid (ABS), black base (ABS), each. M920-2	\$158.00



Staining Dishes

Glass (flint)

20-Slide (back-to-back), holds slides horizontally in grooves	
Set (grooved dish & lid), USA. 672	\$12.30
Bottom only. 672-BTM	\$8.65
Lid only. 672-LID	\$6.00





20-Slide, holds slides horizontally in metal rack	
Approx. inside dimensions: 9.6 x 7.6 x 5.7 cm(H)	
Set (dish, rack & lid). SDS-20	\$25.65
Dish & Lid. SD-20	\$15.65
Rack only, stainless steel. SR-20	\$12.90
30-Slide, holds slides vertically in metal rack	
Approx. inside dimensions: 9.6 x 8.5 x 8.3 cm(H)	
Set (dish, rack & lid). SDS-30	\$35.00
Dish & Lid. SD-30	\$17.50
Rack only, stainless steel. SR-30	\$21.50
60-Slide, holds slides vertically in metal rack	
Approx. inside dimensions: 15.5 x 8.3 x 8.5 cm(H)	
Set (dish, rack & lid). SDS-60	\$49.90
Rack only, stainless steel. SR-60	\$28.15



Plastic rack and dish resist alcohol and xylene; can be microwaved or frozen. Lightproof staining dish has a tight lid with a perforated slide partition for supporting a drying agent. Small size of dish saves on reagents, since only 150 mL is needed (holds 250 mL).

Plastic (polyoxymethylene [POM]), holds 25 slides	
Set (dish, rack & lid). 0990-99	\$50.90



Staining Jars, hold slides (up to 76 x 26 mm) vertically
Coplin, holds 10 slides back-to-back

Glass w/Ground Glass Cover, USA. 673	<u>ea</u>	<u>cs/12</u>
Glass Cover, Fits 673 Jar. 673-LID	\$10.75	119.00
	\$4.20	—
	<u>ea</u>	<u>cs/12</u>
Glass w/Screw-On Phenolic Cap, USA. 674	\$10.75	119.00
Phenolic Cap, Fits 674 Jar, black. 674-LID	\$4.25	—
	<u>ea</u>	<u>cs/12</u>
Polypropylene w/Screw-On Cap. 600P90101	\$4.60	45.20
Plastic jar is unbreakable, microwave proof, stain resistant, and chemically inert. Cap is flat.		
	<u>ea</u>	<u>cs/12</u>
Wide Top, Glass w/Square Lid, USA. 675	\$14.55	160.00
holds 16 slides back-to-back		
Glass Lid, Fits 675 Jar. 675-LID	\$4.45	—



Move slides, without forceps, between jars linked together in a workstation. Swing out the rack's integral lid for easy accessibility to slides. Stable, wide-based jar requires only 80 mL of reagent. Polyoxymethylene (POM) resists alcohol and xylene; can be microwaved or frozen.

Wide Top, Plastic (POM)

12-Slide, EasyDip™, Simport	<u>ea</u>	<u>cs/6</u>
Jar only, Blue. M900-12B	\$10.30	56.70
Jar only, Green. M900-12G	\$10.30	56.70
Jar only, Pink. M900-12P	\$10.30	56.70
Jar only, White. M900-12W	\$10.30	56.70
Jar only, Yellow. M900-12Y	\$10.30	56.70
Rack/Lid only, gray. M905-12DGY	\$11.45	63.00




Anodized aluminum staining jar holder

(Jars sold separately)
Staining Jar Holder, 6-Place. **M906** \$65.60



25-Slide, light green, microwaveable	
Jar & Lid. 1030	\$11.35
Rack only. 1042	\$18.90
Staining Jar Station, 12-Place. 4501011	\$259.00
Includes metal frame, (12) 25-slide dishes w/ lids (1030) & (1) rack (1042).	
Dimensions: 29(L) x 5(W) x 5"(H)	



Tissue (Embedding) & Biopsy Cassettes



Base Molds

Plastic, disposable, Simport


15 x 15 x 5 mm, Cs/1000. M475-2	\$126.00
24 x 24 x 5 mm, Cs/1000. M475-3	\$126.00
★ 30 x 24 x 5 mm, Cs/1000. M475-4	\$126.00
37 x 24 x 5 mm, Cs/1000. M475-5	\$126.00
37 x 24 x 10 mm, Cs/500. M475-10	\$65.45
Designed for Simport Macrosette™.	
Metal, reusable	
15 x 15 x 5 mm, Pk/12. HB-15	\$18.00
24 x 24 x 5 mm, Pk/12. HB-24	\$18.00
30 x 24 x 5 mm, Pk/12. HB-30	\$18.00
37 x 24 x 5 mm, Pk/12. HB-37	\$18.00
66 x 45 x 15 mm, Simport, Pk/10. M474-7	\$156.00
Designed for Simport Jumbosette™.	



Biopsy Bags


Nylon, 30 x 50 mm, Pk/1000. 540478-1	\$235.00
Paper (Tea Bags), 2¼ x 2¼", Pk/500. 56100	\$38.00

M476-1
Made of 60 PPI cellular polyester urethane foam. Use rectangles with cassettes, circles with capsules.



Biopsy Pads

Rectangular	
Black Foam, 30.2 x 25.4 mm, Bg/1000. 339014	\$40.00
Blue Foam, 30.2 x 25.4 x 2 mm, Cs/10,000. M476-1	\$185.00
Circular, Simport	
Blue Foam, 25.4 x 2.7 mm, Cs/10,000. M476-2	\$225.00



Molded from a high-density, solvent resistant polymer. Textured writing surface on front and sides.

Biopsy Cassettes, Detached Lids (packed separately)

Economy, White, small square pores, Pk/500. **BC-0109W** \$40.65

Tissue Cassettes, Detached Lids, slotted pores

Economy, Pk/500 \$40.65

Blue. **EC-0109B** Pink. **EC-0109P**
 Green. **EC-0109G** White. **EC-0109W**
 Lilac. **EC-0109L**

Histosette™ II, Simport, Pk/1000 1 pk 3 pk
 \$84.00 241.00
 45° front labeling area.

Blue. **M485-6** Pink. **M485-3**
 Green. **M485-4** White. **M485-2**

Uni-Cassettes™, USA, Pk/500

\$59.30
 Aqua. **6059** Lime. **6062**
 Blue. **6054** Peach. **6064**
 Burnt Orange. **6056** Pink. **6052**
 Chartreuse. **6061** Red. **6060**
 Gold. **6058** Rose. **6063**
 Green. **6053** Tan. **6055**
 Lavender. **6065** White. **6050**
 Lilac. **6057** Yellow. **6051**



Simport
Histosette I



Simport
Micromesh

Biopsy Cassettes, Attached Lids, acetyl polymer 1 pk 3 pk

Histosette™ I, 1 mm square pores, Pk/1500 \$122.00 348.00
 30° front labeling area.

Blue. **M491-6** Pink. **M491-3**
 Green. **M491-4** White. **M491-2**

Unisette™, hinged lids, Pk/1500 1 pk 3 pk
 \$130.00 373.00
 35° front labeling area. (Alternative to Sakura's Uni-Cassette®)

Blue. **M506-6** Pink. **M506-3**
 Green. **M506-4** White. **M506-2**

Micromesh™, 0.38mm pores, Pk/1000 1 pk 3 pk
 \$206.00 552.00
 45° front labeling area.

Blue. **M507-6** Pink. **M507-3**
 Green. **M507-4** White. **M507-2**



Simport



Simport
Macrosette

Tissue Cassettes, Attached Lids, slotted pores 1 pk 3 pk

Histosette™ I, snap-apart lids, Pk/1500 \$122.00 348.00
 30° front labeling area.

Blue. **M490-6** Pink. **M490-3**
 Green. **M490-4** White. **M490-2**

Histosette™ II, integral lids, Pk/1500 1 pk 3 pk
 \$165.00 474.00
 45° front labeling area.



Blue. **M492-6** Pink. **M492-3**
 Green. **M492-4** White. **M492-2**

Unisette™, hinged lids, Pk/1500 1 pk 3 pk
 \$130.00 373.00
 35° front labeling area. (Alternative to Sakura's Uni-Cassette®)

Blue. **M505-6** Pink. **M505-3**
 Green. **M505-4** White. **M505-2**

Macrosette™, 2X normal depth, white, Pk/750. **M512** \$172.00
 Dimensions: 40.1 x 28.5 x 13 mm(H)

Jumbosette™, white, Pk/100. **M513-2** \$140.00
 Dimensions: 75 x 52 x 18 mm(H)



M480-4

Cassettes, without Lids (compatible with all available metal lids)

Tissue, acetyl polymer, premium

Capsettes™, USA, Pk/500 \$40.80

Aqua. **6039** Lime. **6042**
 Blue. **6034** Peach. **6044**
 Burnt Orange. **6036** Pink. **6032**
 Chartreuse. **6041** Red. **6040**
 Gold. **6038** Rose. **6043**
 Green. **6033** Tan. **6035**
 Lavender. **6045** White. **6030**
 Lilac. **6037** Yellow. **6031**

Simport, 45° front labeling area
 Green, Pk/1500. **M480-4** \$92.00

Cassette Storage

Box, Cardboard, white
 14 x 14 x 2", Pk/50. 8342..... \$121.00
 Each box can store approximately 540 paraffin blocks.

Store 3,150 cassettes or 1,950 blocks
 in 4 drawers, each divided into 12 rows.
 Cabinets can be stacked up to 8 high.
 Dimensions: 49.5 x 38 x 21 cm.



Tissue Cassette Cabinet, Steel, 4-drawer. 485400 \$239.00



Store 1500 cassettes or 990 embedding
 rings in 6 stackable drawers.
 High-impact polystyrene resists
 temperatures from -80°C to +80°C.
 Drawer dimensions: 40.5 x 23 x 5.1 cm.

Tissue Cassette Drawers, Polystyrene, Pk/6. M495-6..... \$150.00

Remove paraffin, ice, etc., safely
 from your laboratory equipment
 with these plastic scrapers.



Scrapers, plastic
 5 cm (2") Wide Blade. 311-245 \$3.70
 10 cm (4") Wide Blade. 311-250 \$4.70
 15 cm (6") Wide Blade. 311-255 \$5.70



Tissue Flotation Bath. XH-1001..... \$191.00

Ultra pure and nonflammable;
 chills to -62°F.



Freeze Spray, 10 oz Can. 49185..... \$12.00



Histology Fixative, B-Plus Fix™, 60 mL, Cs/256. MA0102015 \$261.00



Microtome Blades

Accu-Edge®, Sakura Finetek
 Low Profile, 80 x 8.3 mm, Pk/50. 4689..... \$134.00
 High Profile, 75.7 x 14 mm, Pk/50. 4685..... \$134.00
 DURAEDGE®, PTFE coated, USA
 Low Profile, 77 x 8 mm, Pk/50. 7223..... \$64.80
 High Profile, 77 x 14 mm, Pk/50. 7310..... \$64.80
 Encore™, PTFE coated, USA
 Low Profile, 83 x 8 mm, Pk/50. 7380..... \$79.80
 High Profile, 83 x 14 mm, Pk/50. 7390..... \$79.80