

**Slide & Cassette Pens, Markers & Scribes**

Write a 0.4 mm stroke on glass slides with this permanent marker. Pilot pens are ideal for cytologists and pathologists.



**Pens**

Pilot®, permanent, extra fine point	<u>ea</u>	<u>bx/12</u>
Black. <b>44102</b> .....	\$1.70	17.45
★ Blue. <b>44103</b> .....	\$1.70	17.45
Green. <b>44105</b> .....	\$1.70	17.45
Red. <b>44104</b> .....	\$1.70	17.45



Confine immunohistochemical reagents within a light blue, hydrophobic barrier, or separate multiple specimens on a slide/dish. Resists alcohol and acetone, but removes with xylene.

Super <sup>HT</sup> PAP (Peroxidase Anti-Peroxidase), hydrophobic		
Regular Tip (4 mm tapered nib). <b>K-600</b> .....	\$59.00	
(71310, 195505)		
★ Fine Tip (2.5 mm tapered nib). <b>K-650</b> .....	\$44.00	
(71312, 195506)		

Label any glossy surface with this water-resistant wax pencil. Wipes off easily with a cloth. Simply pull the string to remove wrapping and expose more tip.



**Markers, USA**

China Marker (Wax Pencil), Bx/12. <b>02089</b> .....	\$11.50
--	---------

Ink resists xylene, alcohol, acetone and formalin. Remove with acetic acid, chloroform, or ethanol/xylene mixture.



Securline® Marker II®, permanent, PDC	<u>ea</u>	<u>bx/10</u>
Black, 0.6 mm tip. <b>1450-20-PDC</b> .....	\$2.50	20.30
(14-905-30, 72107)		



Sharpie Retractable, permanent		
Fine Tip	<u>ea</u>	<u>bx/12</u>
Black. <b>36701</b> .....	—	29.65
Blue. <b>36703</b> .....	\$2.85	29.65
Red. <b>36702</b> .....	\$2.85	29.65
Ultra-Fine Tip	<u>ea</u>	<u>bx/12</u>
Black. <b>1735790</b> .....	\$2.85	29.65
Blue. <b>1735792</b> .....	—	29.65



**Scribes, retractable**

Carbide-Tipped (resharpenable). <b>RSC</b> .....	\$15.30
Carbide Tip (replacement). <b>RCT-1</b> .....	\$6.20
Diamond-Tipped. <b>RSD</b> .....	\$21.85
Diamond Tip (replacement). <b>RDT-1</b> .....	\$12.90




Evenly warm your slides from room temperature to 70°C. Black anodized surface provides excellent contrast.

**Slide Warmers**

Small, 26 x 18 cm (10 x 7"). <b>XH-2002</b> .....	\$251.00
Large, 63 x 20 cm (25 x 8"). <b>XH-2001</b> .....	\$335.00

**Staining**



**Staining/Moisture Chambers**

10-Slide, polystyrene, <b>USA</b>	<u>ea</u>	<u>pk/3</u>
Clear. <b>240-9000-010</b> .....	\$54.95	154.00
Black. <b>240-9020-Z10</b> .....	\$56.10	157.00



20-Slide, Simport		
★ Clear Lid (PETG), black base (ABS), each. <b>M920-1</b> .....	\$160.00	
Black Lid (ABS), black base (ABS), each. <b>M920-2</b> .....	\$160.00	

Continued



**Incubators**, myTemp™ Mini Digital, Benchmark

- 20 L capacity (0.8 cu ft)
- Digital temperature control with LED display
- Two models: Heat Only, or Heat/Cool
- Accepts bottles and flasks up to 2 L
- Two adjustable/removable shelves
- 9.3"(D) x 10.3"(W) x 12.8"(H) interior;  
14.5"(D) x 13.2"(W) x 18.7"(H) exterior
- 2-year warranty

Heat Only, 1° above ambient to 60°C. **H2200-H** ..... \$413.00  
 Heat/Cool, 15° below ambient to 60°C. **H2200-HC** ..... \$637.00  
 Heat/Cool/CO<sub>2</sub>, 10° below ambient to 55°C. **H2300-HC2** \$3,089.00

Label any glossy surface with this water-resistant wax pencil. Wipes off easily with a cloth. Simply pull the string to remove wrapping and expose more tip.



**Markers, USA**

China Marker (Wax Pencil), Bx/12. **02089** ..... \$11.50



Sharpie Retractable, permanent

Fine Tip	<u>ea</u>	<u>bx/12</u>
Black. <b>36701</b> .....	---	\$29.65
Blue. <b>36703</b> .....	\$2.85	29.65
Red. <b>36702</b> .....	\$2.85	29.65
Ultra-Fine Tip	<u>ea</u>	<u>bx/12</u>
Black. <b>1735790</b> .....	\$2.85	29.65
Blue. <b>1735792</b> .....	---	29.65

**Mortar and Pestle Sets**

Glass, **USA**

59 mL (2 oz) Capacity. **826-2** ..... \$15.50  
 118 mL (4 oz) Capacity. **826-3** ..... \$16.25  
 237 mL (8 oz) Capacity. **826-4** ..... \$19.30



Porcelain, outside wall and rim glazed

30 mL Capacity, 63 mm top OD. **LPC-154** ..... \$4.70  
 80 mL Capacity, 80 mm top OD. **LPC-154-1** ..... \$6.30  
 100 mL Capacity, 105 mm top OD. **LPC-155** ..... \$7.75  
 320 mL Capacity. **LPC-156** ..... \$10.55

**Petri Dishes**



**Petri Dishes**

- PS (polystyrene)  
 90 x 20 mm, economy  
 Pk/10, nonsterile. **51651** ..... \$5.00  
 35 x 10 mm, vented lid, stacking ring, Simport  
 Cs/500, sterile sleeve/20. **D210-15** ..... \$132.00  
 60 X 15 mm, stackable, Kord-Valmark  
 Cs/500, sterile sleeve/20. **2901** ..... \$92.00 ↓  
 100 X 15 mm, Kord-Valmark  
 Slippable (rimless lid for automation)  
 Standard, Cs/500, sterile sleeve/25. **2917** ..... \$103.00  
 ISO Mark, Cs/500, sterile sleeve/25. **2910P** ..... \$103.00  
 ★ Stackable  
 Standard, Cs/500, sterile sleeve/25. **2900** ..... \$106.00  
 100 X 25 mm, stackable, Kord-Valmark  
 Cs/340, sterile sleeve/17. **2906** ..... \$106.00 ↓  
 Glass (borosilicate), reinforced rims, economy  
 90 x 15 mm, each. **31015003** ..... \$2.75  
 100 x 15 mm, each. **31015004** ..... \$3.20



Store (54) 100 mm petri dishes in this clear, 6-column rack. 34(L) x 20(W) x 23 cm(H).

**Petri Dish Rack**, 100 mm, acrylic. **189850100** ..... \$123.00