

Student Microscopes



National

185 Student Microscope, Monocular \$89.00

- Upright and unreversed image
- Widefield 10X eyepiece
- Monocular, 45° inclined viewing head rotates 360°
- 2X objective provides 20X magnification
- Stage has two clips and reversible black/white plate
- Rack and pinion focusing
- 26 cm (10¼") tall
- 5-year warranty
- Grades K-12



C&A Scientific

MFL-05 Student Microscope, Monocular \$130.00

- Widefield 10X eyepiece with pointer
 - Monocular, 45° inclined viewing head rotates 360°
 - Achromat 4X, 10X and 40XR objectives
 - Stage provided with safety stop to protect slides & objectives
 - Coarse focus knob
 - In-stage condenser with 6-place disc diaphragm
 - Adjustable intensity LED illuminator, equivalent to a 20 watt tungsten bulb, operates on rechargeable battery or power cord
 - 27.8 cm (11") tall
 - 1-year warranty
 - Five prepared glass slides
 - Grades K-8
-



National

131 Student Microscope, Monocular, DIN..... \$217.00

- Widefield 10X eyepiece w/ pointer accepts optional reticle
- Monocular, 45° inclined viewing head rotates 360°
- DIN achromat 4X, 10X, and 40XR objectives
- Stage provided with safety stop to protect slides & objectives
- Low positioned coarse and fine focus; tension adjustment
- In-stage 0.65 N.A. condenser with 5-place disc diaphragm
- 20 watt tungsten illuminator with blue filter
- 36 cm (14") tall
- Lifetime mechanical and optical/3-year electrical warranty
- Grades K-12

National

131-CLED-MS Student Microscope, Monocular, DIN \$283.00

- Widefield 10X eyepiece w/ pointer accepts optional reticle
- Monocular, 45° inclined viewing head rotates 360°
- DIN achromat 4X, 10X, and 40XR objectives
- Graduated mechanical stage provided with safety stop to protect slides & objectives
- Low positioned coarse and fine focus; tension adjustment
- In-stage 0.65 N.A. condenser with 5-place disc diaphragm
- LED illuminator, corded (equivalent to a 20 watt tungsten bulb)
- 36 cm (14") tall
- Lifetime mechanical and optical/3-year electrical warranty
- Grades K-12