

## GENERAL LAB SUPPLIES & EQUIPMENT

**Alcohol Prep Pads**, Sterile, 70% Isopropanol  
Large, 4.5 x 9 cm (1.8 x 3.5"), Cs/1000. **44058-404**..... \$40.00



### Applicators, 6" long

Cotton Tipped, wooden shaft

Nonsterile

Bg/100. **95-8702**..... \$1.50

Bx/1000, 10 bg/bx. **95-8702**..... \$12.00

Sterile Pk/1, Bx/100. **4404304**..... \$7.80

Sterile Pk/2, Bx/200. **44024-106-2S**..... \$7.90

Plain, wood, Bg/1000. **95-8701**..... \$7.25



Rubberized Cloth

### Aprons

Disposable, polyethylene, 46"(L) x 28"(W)

0.8 mil Thick (0.02 mm), Bx/100. **105DA2846**..... \$10.95

1.75 mil Thick (0.044 mm), Bx/100. **59010485**..... \$38.00

Rubberized Cloth, black

36"(L) x 27"(W), Pk/10. **3101276010**..... \$66.50

42"(L) x 27"(W), Pk/10. **3101276110**..... \$76.00




42"(L) x 36"(W), Pk/10. **3101276310**..... \$90.00

46"(L) x 36"(W), Pk/10. **3101276210**..... \$96.30

### Bags



Collect biohazardous waste in this space-saving, autoclavable bag. Printed with a biohazard symbol and sterilization indicator patch. Measures 20.3 x 30.5 cm (8 x 12") and is 1.57 mil (0.04 mm) thick polypropylene.



HS10044H

**Bag, Biohazard, Benchtop**, Bx/500. **HS10320**..... \$69.35

Biohazard Bag Holder, Benchtop. **HS10044H**..... \$20.85

Epoxy-coated wire holder measures 9.1 x 12.8 x 21.4 cm.



Bags display a biohazard symbol and precautionary procedures. Sterilization indicator patch darkens when autoclaved. Polypropylene bags are orange-red. Durable, extra heavy-duty bags provide maximum protection against dangerous tears and leaks. Bel-Art

**Bags, Biohazard, Standard**, 1.5 mil (0.038 mm) thick

14 x 19" (36 x 48 cm), Bx/200. **131641419**..... \$125.00

19 x 23" (48 x 58 cm), Bx/200. **131641923**..... \$182.00

25 x 35" (64 x 89 cm), Bx/200. **131642535**..... \$326.00

**Bags, Biohazard, Heavy-Duty**, 2 mil (0.05 mm) thick

14 x 19" (36 x 48 cm), Bx/200. **131651419**..... \$191.00

19 x 23" (48 x 58 cm), Bx/200. **131651923**..... \$279.00

25 x 35" (64 x 89 cm), Bx/200. **131652535**..... \$505.00



**Bags, Reclosable**, polyethylene

2 mil Thick, Clear

3 x 7", Pk/100. **17730-100**..... \$7.20

2 mil Thick, White Block for labeling

3 x 4", Pk/100. **2012-100**..... \$7.40

6 x 9", Pk/100. **1309-100**..... \$9.20

6 x 12", Pk/100. **10850-100**..... \$15.40

4 mil Thick, clear

3 x 4", Pk/100. **494-100**..... \$12.00

Safely transport patient specimens with these OSHA and NCCLS compliant bags. Each resealable, 2 mil thick polyethylene bag has an attached document pouch and imprinted biohazard symbol.



**Bag, Specimen Transport**, 6 x 9", Bx/1000. **4919**..... \$56.70



**Bags, Waste Disposal**, HDPE, Cleanware™

8 x 10", Pk/50. **131741008**..... \$19.90

12 x 16", Pk/50. **131741012**..... \$33.45

Sliding blade cuts excess material 19 mm from the seal.



**Bag Sealer**, electric, 12" (30.5 cm) Wide. **HS28272**..... \$168.00

Bag Sealer Replacement Part Kit. **HS28280**..... \$30.65

Poly Tubing

4 mil

4" x 1500'. **1143**..... \$110.00

6" x 1500'. **1144**..... \$162.00

6 mil

4" x 1000'. **1658**..... \$132.00

6" x 1000'. **1659**..... \$162.00

Poly Tubing Dispenser. **H-28**..... \$34.00

Continued