

Balances

Balance, Portable



Model	DWP-1001
Capacity (g)	1000
Readability (g)	0.1
Modes	g, oz, ct, gn
Pan	14.6 cm (5.75") stainless steel
Power	4 AA/AC adapter
Price	\$72.00

Warranty: 1 year



Benchmark
W3200-1200

Balances, Toploading

Model	W3200-500	W3200-1200
Capacity (g)	500	1200
Readability (g)	0.001	0.01
Modes	g, oz, ct	g, oz, ct
Pan	11.5 cm (4.5")	16 cm (6.3")
Draft Shield	Yes, 3-door glass	No
Price	\$999.00	\$808.00

Warranty: 2 years



Benchmark

Balances, Analytical

Model	★ W3100-120	W3100A-120
Calibration	External (100 g weight provided)	Internal (Quick-Cal™ software)
Capacity (g)	120	120
Readability (g)	0.0001	0.0001
Modes	g, oz, lb, ct	g, oz, lb, ct
Pan	9 cm (3.5")	9 cm (3.5")
Price	\$1,487.00	\$1,580.00

Warranty: 2 years

Baths



BSH1002

Aluminum
Blocks

Dry Baths, digital, Benchmark (2-year warranty)

1-Block Capacity. BSH1001	\$359.00
2-Block Capacity. BSH1002	\$425.00
4-Block Capacity. BSH1004	\$725.00

Blocks sold separately – please specify configuration

Blocks, aluminum

Holds (48) 0.2 mL Tubes. BSW02	\$85.00
Holds (24) 0.5 mL Tubes. BSW05	\$85.00
Holds (24) 1.5 mL Tubes. BSW1500	\$85.00
Holds (24) 1.5/2.0 mL Tubes. BSW1520	\$85.00
Holds (20) 12/13 mL Tubes. BSW13	\$85.00
Holds (12) 15 mL Tubes. BSW15	\$85.00
Holds (5) 50 mL Tubes. BSW50	\$85.00



Racks are only included with model B2000-4.

Water Baths, digital, MyBath™, Benchmark (2-year warranty)

4 L, 3 tube racks. B2000-4	\$739.00
8 L. B2000-8	\$799.00
12 L. B2000-12	\$956.00



Replace aluminum blocks, water, and ice in baths and buckets with this thermal metallic media.



Lab Armor™ Beads, USA

0.75 L (replaces two blocks in a dry bath). 42370-750	\$98.00
2 L. 42370-002	\$191.00
(B-0700-2L)	
4 L. 42370-004	\$314.00
(B-0700-4L)	
8 L (2 x 4 L). 42370-008	\$490.00
(B-0700-8L)	