



Bottle Top Dispensers

- 0.5-5 mL. **600BTDR2**\$249.00
- 5-60 mL. **600BTDR5**\$295.00



Matchbox Style

Box, Cardboard, Matchbox Style (Pill Box), white

- 2 1/8 x 1 1/8 x 5/8", Pk/72. **20047**\$31.10

Box, Cardboard, 14 x 14 x 2", white, Pk/50. 8342\$121.00

Each box can store approximately 540 paraffin blocks.



246530002

246530001

Boxes, Cardboard, Glass Disposal

- Bench Model, 8 x 8 x 10", Pk/6. **246530002**\$65.55
- Floor Model, 12 x 12 x 27", Pk/6. **246530001**\$91.05

Hinged
Polystyrene



Boxes, Plastic, hinged polystyrene (see following chart)

Plastic Boxes				
Catalog #	Inside Dimensions	Foam Insert	Suggested Application	Price
H-1	1 x 1 x 1/4"	-	Reticles, Pointers	\$0.40
H-1/11B	1 x 1 x 1/4"	1	Reticles, Pointers	\$0.50
H-54	2 9/16 x 2 9/16 x 1/2"	-	Filters	\$1.20
H-54/2916G	2 9/16 x 2 9/16 x 1/2"	1	Filters	\$1.35
H-28/334B2	3 3/4 x 1 3/16 x 5/8"	2	Slides, Stage Micrometers	\$1.40
H-10	2 7/8 x 1 3/16 x 3/8"	-	Small Parts	\$0.60
H-21/512B	5 1/2 x 1 3/4 x 1/2"	1	Forceps & Scissors	\$1.80


BRS-969.10



LB-40



LB-38



LB-30



LB-14



Brushes

- Bench, horse hair, 10" long. **BRS-969.10**\$9.40
Bristle head: 4 3/4"(L) x 1 1/2"(W); bristle length: 1".
- Brass Wire, 9 3/4" long. **LB-40**\$8.70
Bristle head: 5"(L) x 3/4"(W); bristle length: 1 1/8".
- Flat Bristles, paint style, camel hair
1/2" Wide, 6 1/4" total length. **LB-36**\$3.90
- 1" Wide, 7" total length. **LB-37**\$6.15
- 2" Wide, 7 3/4" total length. **LB-38**\$5.30
- Round Bristles, camel hair, 8 3/4" long. **LB-30**\$10.25
- Swivel-Head, black nylon bristles. **LB-14**\$4.20
Bristle head: 4 1/4"(L) x 1 1/2"(W). Total length: 12".

Remove debris, rust, or corrosion in hard to reach places.



- Pen Style (scratch), retractable, 11.5 cm (4 1/2") long, Germany
Brass. **BRS-295.00**\$12.05
- Nylon. **BRS-293.00**\$9.55
- Steel. **BRS-292.00**\$12.60



H-5870 vertical mixing tube is equipped with a flame stabilizer and adjustable air supply ports. Serrated gas inlet accepts 5/16" ID tubing.

H-5550 gridded top, 3.8 cm (1 1/2") diameter, produces a uniform blue cone of light.



Burners, Bunsen, natural gas, Humboldt Mfg, USA

- Standard, flame stabilizer. **H-5870**\$22.45
- High Temperature, Meeker-type. **H-5550**\$78.90



- Burner Lighter**, replaceable single flint, round file. **SOL-290.00** ...\$3.40
- Flints (replacement), Pk/5. **SOL-290.01**\$2.00