

## Filter & Blotting Papers



Ahlstrom - Munksjo  
Filter Paper

### Qualitative, Standard Grades (for gravity filtration)

Slow flow (4 mL/min), 1.5 µm retention, <b>Grade 610</b> (Whatman 5; S&S 602)	
9 cm, Pk/100. <b>6100-0900</b> .....	\$9.70
12.5 cm, Pk/100. <b>6100-1250</b> .....	\$12.25
15 cm, Pk/100. <b>6100-1500</b> .....	\$14.75
20.5 cm, Pk/100. <b>6100-2050</b> .....	\$23.30
24 cm, Pk/100. <b>6100-2400</b> .....	\$33.05
Medium flow (30 mL/min), 2 µm retention, <b>Grade 642</b> (Whatman 2, 2CHR, 20CHR; S&S 2043B, 2045A/B)	
9 cm, Pk/100. <b>6420-0900</b> .....	\$7.55
15 cm, Pk/100. <b>6420-1500</b> .....	\$13.70
Medium flow (55 mL/min), 2.5 µm retention, <b>Grade 601</b> (Whatman 1, 1CHR; S&S 593, 597, 2043A)	
4.25 cm, Pk/100. <b>6010-0425</b> .....	\$5.80
9 cm, Pk/100. <b>6010-0900</b> .....	\$7.80
12.5 cm, Pk/100. <b>6010-1250</b> .....	\$11.30
15 cm, Pk/100. <b>6010-1500</b> .....	\$13.75
18.5 cm, Pk/100. <b>6010-1850</b> .....	\$19.95
20.5 cm, Pk/100. <b>6010-2050</b> .....	\$21.50
33 cm, Pk/100. <b>6010-3300</b> .....	\$58.05
Medium flow (60 mL/min), 6 µm retention, <b>Grade 613</b> (Whatman Student; S&S 595, 0903)	
9 cm, Pk/100. <b>6130-0900</b> .....	\$6.65
15 cm, Pk/100. <b>6130-1500</b> .....	\$11.65
20.5 cm, Pk/100. <b>6130-2050</b> .....	\$17.50
33 cm, Pk/100. <b>6130-3300</b> .....	\$38.90
Medium flow (200 mL/min), 10 µm retention, <b>Grade 631</b> (Whatman 4, 4CHR; S&S 604, 2040A/B; Fisher P8)	
20.5 cm, Pk/100. <b>6310-2050</b> .....	\$21.70
33 cm, Pk/100. <b>6310-3300</b> .....	\$58.60
Fast flow (235 mL/min), 25 µm retention, creped, <b>Grade 615</b>	
9 cm, Pk/100. <b>6150-0900</b> .....	\$6.85
12.5 cm, Pk/100. <b>6150-1250</b> .....	\$10.00
15 cm, Pk/100. <b>6150-1500</b> .....	\$12.00
20.5 cm, Pk/100. <b>6150-2050</b> .....	\$18.00
33 cm, Pk/100. <b>6150-3300</b> .....	\$40.00
Very Fast flow (360 mL/min), 35 µm retention, creped, <b>Grade 617</b>	
9 cm, Pk/50. <b>6170-0900</b> .....	\$6.25
12.5 cm, Pk/50. <b>6170-1250</b> .....	\$8.40
15 cm, Pk/50. <b>6170-1500</b> .....	\$9.80
20.5 cm, Pk/50. <b>6170-2050</b> .....	\$14.90
33 cm, Pk/50. <b>6170-3300</b> .....	\$31.90
<b>Qualitative, Wet-Strengthened Grades</b> (for vacuum filtration)	
Medium flow (65 mL/min), 3 µm retention, <b>Grade 909</b> (Whatman 91; S&S 0860, 1574, 3003; Reeve-Angel 201)	
32 cm, Pk/100. <b>9090-3200</b> .....	\$63.25
Fast flow (245 mL/min), 27 µm retention, creped, <b>Grade 962</b> (Whatman 114; S&S 0856, 0905, 3002; Reeve-Angel 202)	
11 cm, Pk/100. <b>9620-1100</b> .....	\$9.40

## Quantitative, Ashless

(Designed for use in analytical and gravimetric analysis)

Slow flow (4 mL/min), 1.5 µm retention, <b>Grade 94</b> (Whatman 42; S&S 589/3)	
9 cm, Pk/100. <b>0940-0900</b> .....	\$27.15
15 cm, Pk/100. <b>0940-1500</b> .....	\$46.70
Medium flow (30 mL/min), 2 µm retention, <b>Grade 74</b> (Whatman 40; S&S 589/2, 589/4)	
9 cm, Pk/100. <b>0740-0900</b> .....	\$27.15
15 cm, Pk/100. <b>0740-1500</b> .....	\$46.70
Fast flow (220 mL/min), 15 µm retention, <b>Grade 55</b> (Whatman 54, 541; S&S 1505)	
9 cm, Pk/100. <b>0550-0900</b> .....	\$30.65
15 cm, Pk/100. <b>0550-1500</b> .....	\$49.25

## Blotting/Chromatography/Electrophoresis Papers

0.19 mm thick, 125 mm/30 min absorption, <b>Grade 601</b> (Whatman 1, 1CHR; S&S 593, 597, 2043A)	
46 x 57 cm, Pk/100. <b>6018-4657</b> .....	\$110.00
0.34 mm thick, 140 mm/30 min absorption, <b>Grade 238</b> (Whatman 3MM, S&S 598, 2316)	
7 x 10 cm, Pk/100. <b>2388-0710</b> .....	\$26.05
20 x 20 cm, Pk/100. <b>2388-2020</b> .....	\$57.25
46 x 57 cm, Pk/100. <b>2388-4657</b> .....	\$209.00
0.42 mm thick, 150 mm/30 min absorption, <b>Grade 237</b> (Whatman 3, 3CHR; S&S 22/2, 591, 2469, 5703)	
46 x 57 cm, Pk/100. <b>2378-4657</b> .....	\$240.00
0.95 mm thick, 180 mm/30 min absorption, <b>Grade 243</b> (Whatman 17CHR; S&S 2668)	
11 x 14 cm, Pk/100. <b>2438-1114</b> .....	\$24.75
20 x 20 cm, Pk/100. <b>2438-2020</b> .....	\$45.95
46 x 57 cm, Pk/100. <b>2438-4657</b> .....	\$352.00
2.48 mm thick, 240 mm/30 min absorption, <b>Grade 320</b> (S&S 2727)	
11 x 14 cm, Pk/50. <b>3208-1114</b> .....	\$39.15
16 x 18 cm, Pk/50. <b>3208-1618</b> .....	\$54.55
20 x 20 cm, Pk/50. <b>3208-2020</b> .....	\$66.55



Smooth tips &  
polypropylene grips

<b>Forceps, Filter, 11.5 cm. 300-025</b> .....	\$4.40
--	--------



<b>Flask Weights</b> , vinyl-coated lead	
C-Shaped lead rings, blue, Heathrow	
Fits 250-1000 mL Flasks, 51 mm ID, 472 g. <b>HS8881B</b> .....	\$27.50
Fits 500-2000 mL Flasks, 54 mm ID, 599 g. <b>HS8881C</b> .....	\$30.35
O-Shaped lead ring, blue, Heathrow	
Fits 1000-4000 mL Flasks, 70 mm ID, 970 g. <b>HS8882D</b> .....	\$30.35
Lead Wire, 300 mm long, 114 g. <b>183240000</b> .....	\$31.70

Continued