

## Stereo Microscopes



420T-430PHF-10

### National

**420-430PHF-10** Stereo, Binocular, 10-40X Zoom ..... \$569.00

- Widefield 10X/20 eyepieces; accept 23 mm reticle
- Binocular, 45° inclined head with dual diopters rotates 360°; 55-75 mm interpupillary adjustment
- Paired achromat objectives provide continuous 10X-40X total magnification; working distance of 79 mm with supplied 1X objective
- Heavy-duty rack and pinion focusing features a slip clutch and tension adjustment
- Dual illuminators offer light from top only (12V/15W halogen incident light with rheostat control and directional beam), bottom only (5W fluorescent transmitted light), or both top and bottom
- Post stand with locking support collar allows viewing specimens up to 93 mm tall
- Frosted glass and reversible black/white plastic stage plates included
- 24.2 cm(D) x 17.1 cm(W) x 43.3 cm(H)
- Lifetime mechanical and optical/3-year electrical warranty

**420T-430PHF-10** Stereo, Trinocular, 10-40X Zoom ..... \$659.00

- Trinocular version of above
- 2X SLR adapter & 0.4X C-mount adapter included



425-445PLL-10

### National

**425-445PLL-10** Stereo, Binocular, 7.5-50X Zoom ..... \$914.00

- Widefield 10X/23 eyepieces; accept 24.95 mm reticle
- Binocular, 45° inclined head with dual diopters rotates 360°; 55-75 mm interpupillary adjustment
- Paired achromat objectives provide continuous 7.5-50X total magnification; working distance of 113 mm with supplied 1X objective
- Heavy-duty rack and pinion focusing features a slip clutch and tension adjustment
- Dual 3W LED illuminators with rheostat control offer cool light from top, bottom, or both top and bottom
- Post stand has a large ergonomic base and locking support collar
- Frosted glass and reversible black/white plastic stage plates included
- Lifetime mechanical and optical/3-year electrical warranty


**425T-445PLL-10** Stereo, Trinocular, 7.5-50X Zoom ..... \$973.00

- Trinocular version of above



Meiji ..... Access [www.meijitechno.com](http://www.meijitechno.com)  
and call us to place an order

## Stands



*Please call regarding availability of items not listed*

### Stereo Boom Stands, Diagnostic Instruments

*All boom stands require a focus mount*

Standard, nylon bearings. **SMS16B** ..... \$599.00

Ball Bearing

Medium-Duty. **SMS6B** ..... \$759.00

Heavy-Duty. **SMS20-\_\*** ..... See note below

Articulating Arm. **SMS25** ..... \$1,409.00

*\*Access [www.spotimaging.com](http://www.spotimaging.com) and call us with microscope model and boom stand version to place an order*



Meiji PBH Stand  
with Adapter Sleeve

### Stereo Dual Illumination Stands, 6V/10W, Meiji

Pole-Type Base, 84 mm focus mount. **PBH** ..... \$563.00

Arm-Type Base, 84 mm focus mount. **ABZH** ..... \$612.00

*Some stereo brands require an adapter sleeve*

Stereo Stand Sleeve for 76 mm dia. bodies ..... \$30.00

Stereo Stage Plates, 94.5 mm

Clear Glass. **MA569** ..... \$45.00

Black & White Plastic. **MA568** ..... \$45.00