


## Homogenizers & Tissue Grinder



Rapidly homogenize, emulsify, suspend and disrupt 0.1 mL to 250 mL samples. Includes a 5 x 50 mm flat-bottom probe & 7 x 50 mm saw tooth probe for use with micro tubes.

### Homogenizers

Handheld. **D1000** ..... \$919.00



Homogenize a variety of sample types, from soft tissue to bone, in 35 seconds or less.

Benchtop, Microtube, BeadBlaster™24. **D2400** ..... \$6,796.00

- Holds (24) 2.0 mL tubes, or (6) 5 mL tubes w/optional tube holder
- Speed adjusts from 4.00 to 7.00 m/s in 0.05 m/s increments
- Stores 50 programs in memory, each with a maximum of 10 cycles
- 2-year warranty

#### Homogenizer Prefilled Tubes

2.0 mL screw-capped tubes are prefilled with your choice of the following beads:

Glass, acid washed

0.1 mm, Pk/50. **D1031-01** ..... \$99.00

0.5 mm, Pk/50. **D1031-05** ..... \$99.00


Zirconium Beads, acid washed/heat treated/nuclease-free

1.5 mm, Pk/50. **D1032-15** ..... \$139.00

3.0 mm, Pk/50. **D1032-30** ..... \$139.00

Stainless Steel Beads, acid washed

2.8 mm, Pk/50. **D1033-28** ..... \$139.00




### Tissue Grinder, Duall® Glass Pestle

1 mL Capacity, **LG10640100** ..... \$78.50

3 mL Capacity, **LG10640102** ..... \$88.35

5 mL Capacity, **LG10640104** ..... \$100.00



Hydrometers ..... Access [www.belart.com](http://www.belart.com)  
and call us to place an order

## Ice Buckets & Pans



Ice Buckets



Ice Pans

Temperature Range: -196 to 100°C

### Ice Buckets with lids, expanded PVC foam, Bel-Art, USA

2.5 L Capacity (Handy)

Black. **188482002** ..... \$90.90

Blue. **188482001** ..... \$90.90

Purple. **188482006** ..... \$90.90

4 L Capacity (N'Icer)

Black. **188484002** ..... \$111.00

Blue. **188484001** ..... \$111.00

### Ice Pans w/o lids, expanded PVC foam, Bel-Art, USA

1 L Capacity (Mini), 15.5 x 15.5 x 8.5 cm(H)

Black. **188481102** ..... \$67.75

Blue. **188481101** ..... \$67.75

Purple. **188481106** ..... \$67.75

4 L Capacity (Midi), 33 x 24.5 x 12.5 cm(H)

Black. **188484102** ..... \$108.00

Blue. **188484101** ..... \$108.00


Purple. **188484106** ..... \$108.00

9 L Capacity (Maxi), 44 x 35 x 13 cm(H)

Black. **188489102** ..... \$127.00

Blue. **188489101** ..... \$127.00

Purple. **188489106** ..... \$127.00



### Tube Racks, Floating, PVC Foam, Bel-Art, USA

Circular. Recommended for ice buckets

44-Hole, 1.5 mL tubes, 22 cm dia. x 1 cm ..... \$38.90

Black/**188480012** Green/**188480014** Red/**188480013**

30-Hole, 14/15 mL tubes, 22 cm dia. x 1.7 cm ..... \$40.80

Black/**188480022** Green/**188480024** Red/**188480023**

Rectangular. Recommended for ice pans

36-Hole, 1.5 mL tubes, 23.5 x 11 x 1 cm ..... \$38.90

Black/**188480112** Green/**188480114** Red/**188480113**

24-Hole, 14/15 mL tubes, 23.5 x 11 x 1.7 cm ..... \$40.80

Black/**188480122** Green/**188480124** Red/**188480123**