




glass w/rubber bulb; 75 mm (3")

**Medicine (Eye) Dropper, Pk/12. 97-4301** ..... \$3.60



**Pipet Fillers**

3-Valve, rubber, black. **207B** ..... \$14.40  
 Wheel-Type, Pipette Pump™  
 Fits Pipets Up To 2 mL, blue. **378970000** ..... \$24.35  
 Fits Pipets Up To 10 mL, green. **378980000** ..... \$25.45  
 Fits Pipets Up To 25 mL, red. **378990000** ..... \$35.05  
 Wheel-Type, Pipette Pump™ II, quick-release ea pk/12  
 Fits Pipets Up To 2 mL, blue. **379111002** ..... \$26.35 299.00  
 Fits Pipets Up To 10 mL, green. **379111010** ..... \$26.35 299.00  
 Fits Pipets Up To 25 mL, red. **379111025** ..... \$35.05 397.00



Dimensions: 29 x 8.8 x 40 cm

**Pipet Rack, magnetized, blue, Heathrow. HS20612M** ..... \$72.45

**Pipets, Mohr (Measuring), Class A, individually certified**


5 mL, (±0.020 mL), Pk/5. **600PM7058H** ..... \$64.80  
 10 mL, (±0.030 mL), Pk/5. **600PM7058I** ..... \$64.80  
 25 mL, (±0.050 mL), Pk/5. **600PM7058J** ..... \$96.60



Dispense fluids accurately and safely with these color-coded and plugged pipets. Negative graduations provide extra volume (not available on 25 mL). RNase-/DNase-/Pyrogen-free, and noncytotoxic. Calibrated "To Deliver" (TD).

**Pipets, Serological, polystyrene, UniPlast™**

Nonsterile  
 5 mL, blue band, pk/25, Cs/375. **1755** ..... \$74.40  
 10 mL, orange band, pk/25, Cs/250. **1775** ..... \$50.85  
 25 mL, red band, pk/10, Cs/100. **1784** ..... \$47.35  
 Sterile Pk/1 (individually paper/plastic wrapped)  
 1 mL, yellow band, Cs/500. **1700** ..... \$74.65  
 2 mL, green band, Cs/500. **1720** ..... \$80.05  
 5 mL, blue band, reverse grad., Cs/250. **1740** ..... \$59.05  
 10 mL, orange band, reverse grad., Cs/200. **1760** ..... \$51.50  
 25 mL, red band, reverse grad., Cs/100. **1780** ..... \$55.65  
 50 mL, black band, reverse grad., Cs/100. **1820** ..... \$139.00



**Pipets, Transfer, LDPE (low-density polyethylene)**

Graduated, Nonsterile, 15 cm long  
 5 mL capacity, 1 mL graduated, 25 drops/mL  
 ★ #1 Bx/500. **PTP-02** ..... \$16.45  
 Cs/10 Bx. **PTP-02** ..... \$153.00  
 7 mL capacity, 3 mL graduated, 14 drops/mL  
 Bx/500. **PTP-01** ..... \$16.45  
 Cs/10 Bx. **PTP-01** ..... \$153.00  
 Graduated, Sterile Pk/1 (individually wrapped), 15 cm long  
 5 mL capacity, 1 mL graduated, Simport  
 Bx/400. **P200-521S** ..... \$39.65  
 Cs/10 Bx. **P200-521S** ..... \$364.00  
 Graduated, Sterile Pk/20, 15 cm long  
 5 mL capacity, 1 mL graduated, Bx/500. **HS206371B** ..... \$33.30  
 Long Stem, Sterile Pk/1 (individually wrapped)  
 22.5 cm (9" long), 6 mL capacity, Cs/400. **139038** ..... \$59.45  
 Long Stem, Nonsterile  
 30 cm (12" long), 23 mL capacity, Bx/100. **139050** ..... \$17.70

Made of borosilicate glass in compliance with ASTM specification E-969, Class A requirements. Individually serialized and certified. Calibrated "To Deliver" (TD).



**Pipets, Volumetric, Class A, individually certified**

2 mL (±0.006 mL), orange band, Pk/5. **600PT7100B** ..... \$65.30  
 5 mL (±0.010 mL), white band, Pk/5. **600PT7100C** ..... \$68.65  
 10 mL (±0.020 mL), red band, Pk/5. **600PT7100D** ..... \$78.40  
 20 mL (±0.030 mL), yellow band, Pk/5. **600PT7100E** ..... \$111.00



These mini, air displacement pipettors provide an economical and accurate method for pipetting fixed volumes. Pipettors feature a two-step plunger and are fully autoclavable.

**Pipettors, Single Channel**

Fixed Volume, pocket-size, Diamond Jr.  
 10 µL, pink. **3322-10** ..... \$18.20  
 20 µL, green. **3322-20** ..... \$18.20  
 50 µL, purple. **3322-50** ..... \$18.20  
 100 µL, orange. **3322-100** ..... \$18.20



Adjustable Volume, Diamond Advance™ (3-year warranty)  
 0.5-10 µL, 0.02 µL increments, red. **3341-10** ..... \$219.00  
 5-50 µL, 0.1 µL increments, yellow. **3341-50** ..... \$219.00  
 10-100 µL, 0.1 µL increments, blue. **3341-100** ..... \$219.00  
 100-1000 µL, 1 µL increments, orange. **3341-1000** ..... \$219.00

Continued