


**Pipettor Holders**

Carousel, Flip & Grip™, Heathrow



- Blue. **HS120026** ..... \$76.20
- Lime. **HS120056** ..... \$76.20



Stands, Heathrow


- 4-Place, clear acrylic, **USA. HS206204** ..... \$39.30
- 6-Place
- Blue, opaque ABS plastic. **HS20613E** ..... \$38.25
- Clear, acrylic, **USA. HS20620** ..... \$50.90  
(82024-544)



**Pipettor Tips, Filtered (Barrier)** ..... Access [www.globescientific.com](http://www.globescientific.com)  
and call us to place an order

**Pipettor Tips, Standard, Sterile, graduated**

- 100-1250 µL, blue, 6 racks of 96, Bx/576. **151153BRS-96** ..... \$38.55



**Reagent Reservoirs**

- 25 mL, white polystyrene
- Sterile Bg/5, Bx/200. **P8025-5S** ..... \$71.30
- 55 mL, v-bottom, Heathrow
- Nonsterile, clear PVC, Pk/100. **HS20521A** ..... \$40.20
- Sterile Bg/1, white polystyrene, Pk/50. **HS20521C** ..... \$69.80



Polypropylene



Stainless Steel

**Pitchers**

- Polypropylene, graduated (mL), tall form, tapered
- 1000 mL, 50 mL graduations, each. **289910000** ..... \$8.80
- 2000 mL, 100 mL graduations, each. **289920000** ..... \$9.95
- 3000 mL, 100 mL graduations, each. **289930000** ..... \$13.80
- Stainless Steel, graduated (mL/oz), tapered, Polar Ware
- 1000 mL/32 oz, each. **T1063** ..... \$73.20

**Rockers, Shakers & Mixers**



BR1000



B3D1020



B3D1320

**Rockers, Fixed Speed, Benchmark (2-year warranty)**

- 2-D (seesaw), 12 x 12" platform w/flat mat
- BenchBlotter™. **BR1000** ..... \$413.00
- Fixed speed of 12 rocks/minute
- Platform angle is ±8°
- Cold room & incubator safe (4-45°C)
- Load capacity is 2.3 kg (5 lbs)
- 3-D (nutating)
- Small, 10.5 x 7.5" platform w/dimpled mat
- MiniMixer™. **B3D1020** ..... \$350.00
- Large, 12 x 12" platform w/dimpled mat
- BioMixer™. **B3D1320** ..... \$422.00
- Fixed speed of 24 rocks/minute
- Platform angle is ±20°
- Cold room & incubator safe (4-45°C)
- Load capacity is 1.8 kg (4 lbs)



BR2000



B3D2300

**Rockers, Variable Speed**

- Benchmark (2-year warranty)
- 2-D (seesaw), 14 x 12" platform w/flat mat
- BenchRocker™ 2D. **BR2000** ..... \$679.00
- Variable speed of 2-30 rocks/minute
- Platform tilt angle adjusts from 0-30°
- Cold room & incubator safe (4-45°C)
- Load capacity is 3.6 kg (8 lbs)
- 3-D (nutating), 12 x 12" platform w/dimpled and flat mats
- BenchRocker™ 3D. **B3D2300** ..... \$722.00
- Stacking Platform. **BR1000-STACK** ..... \$94.00
- 12 x 12" platform w/flat mat has 3" of clearance.
- Stacking Platform Extenders, Pk/4. **BR2000-SP** ..... \$36.10
- Adds 1.25" to platform separation.
- Variable speed of 2-30 rocks/minute
- Platform tilt angle adjusts from 0-30°
- Cold room & incubator safe (4-45°C)
- Load capacity is 2.3 kg (5 lbs)