



Rockers, Variable Speed

Reliable Scientific (5-year warranty), **USA**

2-D (seesaw)

Standard, 11 x 14" platform w/rubber mat

Single Platform. **55S1114** \$576.00

Double Platform, 3/4" apart. **55D1114** \$634.00

Large, 12 x 16" platform w/rubber mat

Single Platform. **55S1216** \$615.00

Double Platform, 3/4" apart. **55D1216** \$672.00

- Variable speed of 40-160 rocks/minute
- Platform tilt angle adjusts from 0-20°
- Cold room & incubator safe (0-65°C)
- Load capacity is 4.5 kg (10 lbs)



Rocker, Tube, 12 x 4.5" grooved mat, Benchmark. **M2100** \$309.00

- Fixed speed of 24 rocks/minute
- Platform tilt angle is ±30°
- Coldroom & incubator safe (4-45°C)
- Accommodates 1.5 to 50 mL tubes up to 2.3 kg (5 lbs)
- 2-year warranty

Orbital Shakers



BT30

Aerate and mix with these low to high capacity orbital shakers.



BT302

Orbital Shakers, Benchmark (2-year warranty)

Low Speed, 14 x 12" platform w/flat rubber mat

Orbi-Blotter™. **BT30** \$710.00

- Speed adjusts from 3-70 rpm
- Orbit diameter measures 19 mm
- Cold room & incubator safe (4-65°C)
- Load capacity is 3 x 2 L, or 4 kg (9 lbs)

Small, 9.5 x 8" platform w/flat rubber mat

Orbi-Shaker™ JR. **BT302** \$809.00

- Speed adjusts from 50-300 rpm in 10 rpm increments
- Orbit diameter measures 19 mm
- Cold room & incubator safe (4-65°C)
- Load capacity is 4 x 1 L, or 4 kg (9 lbs)
- Timer is programmable from 1 min - 99 hrs, 59 min



BT3000

Large, 13 x 12" platform w/flat rubber mat

Orbi-Shaker™. **BT3000** \$1,306.00

- Speed adjusts from 30-300 rpm in 1 rpm increments
- Orbit diameter measures 19 mm
- Cold room & incubator safe (4-65°C)
- Load capacity is 9 x 500 mL, or 5 x 1 L, or 2 x 2 L, or 6.4 kg (14 lbs)
- Timer is programmable from 1 min - 48 hrs/continuous

Optional Magic Clamp™ platform & accessories sold separately - specify configuration

Steel Platform, Magic Clamp™, 14 x 12". **BT3000-MR**... \$218.00

250 mL Flask Magnetic Clamp. **H1000-MR-250** \$29.95

500 mL Flask Magnetic Clamp. **H1000-MR-500** \$30.40

1000 mL Flask Magnetic Clamp. **H1000-MR-1000** \$46.55

Tube Rotator



Tube Rotator, Reliable Scientific, **USA**. Model **24** \$349.00

- Speed is adjustable from 4-60 rpm
- Cold room & incubator safe (0-65°C)
- 5-year warranty

Rotator disks (carousels) sold separately - specify configuration

24-Place Disk, holds 9-15 mm tubes. **24/9-15** \$89.00

12-Place Disk, holds 15-23 mm tubes. **12/15-23** \$89.00

6-Place Disk, holds 23-29 mm tubes. **6/23-29** \$89.00

Vortexers



BV101



BV1000

Vortex Mixers, Benchmark (2-year warranty)

Vornado™, mini vortexer, fixed 2800 rpm. **BV101**_* \$179.00

*Specify Gray, Blue or Purple cup head

BenchMixer™, single cup head, 200-3200 rpm. **BV1000** \$216.00

Flat Head Attachment, 3". **BV1000-FLAT** \$30.90