


Rulers, Meter Sticks & Tape Measure

Stainless Steel



Rulers, 15 cm/6" long, scaled in both metric and English units

Vinyl, flexible

Transparent. 4461	\$0.95
White. 4460	\$0.55
Stainless Steel, 12 mm wide, USA. 605-065	\$5.50

Meter Sticks, lacquered wood

One side is graduated on both edges in centimeters with 1 mm divisions, while the other side is graduated on one edge in inches with 1/8" divisions.

Plain Ends. 605-050	\$4.90
Metal Ends. 605-051	\$6.15

Tape Measure, 152.4 cm/60", reinforced fiberglass. **11036**

One side is graduated on the top edge in inches with 1/16" divisions and on the bottom edge in cm with 1 mm divisions. The opposite side is graduated in inches with 1/16" divisions.



C300-1R



4402269

Sharps Container, red polypropylene

1 L, Biodisposer™, Simport, each. C300-1R	\$3.15
Dimensions: 8.9 x 8.9 x 15.25 cm (3 x 3 x 6" H)	
5.1 L, horizontal entry lid, Prevent™, each. 4402269	\$7.95
Dimensions: 12" x 4.75" x 11"(H)	




Screw-Cap



Snap-Cap

Specimen Containers, polypropylene (PP)

Screw-Cap	
90 mL, sterile, Pk/300. 840175	\$84.20
Snap-Cap, white HDPE bottom/white LDPE lid	
480 mL (16 oz), Cs/100. 271016W	\$57.15
9.9 cm opening x 11.1 cm high	
1000 mL (34 oz), Cs/100. 271032W	\$83.75
11.7 cm opening x 13.7 cm high.	



10"(D) x 12"(W) x 18"(H)

Splash Shield, 3/16" acrylic, angled top, TrippNT. **50707**



Spot Plates, glazed white porcelain, 21 mm diameter wells

6-Well. LPC-161	\$2.35
12-Well, 115 x 90 x 12 mm. 3104187701	\$4.50

Stirring Bars



FDA approved Teflon® (PTFE) coating is a minimum of 1.6 mm (0.062") thick, except for micro bars which have a 1.32 mm (0.050") thick cover.

Stirring Bars, Teflon®, Bel-Art

Centered Sphere, Saturn Spinbar®	
40 mm (1½") Long. 371160002	\$25.35
50 mm (2") Long. 37116003	\$37.90
Micro, Spinbar®	
White	
5 x 2 mm. 371190005	\$5.05
7 x 2 mm. 371190007	\$5.05
10 x 3 mm. 371190010	\$5.55
12.7 x 3 mm. 371190127	\$5.55
15 x 1.5 mm. 371190015	\$4.95
Red	
7 x 2 mm. 371210013	\$5.55
8 x 1.5 mm. 371210001	\$5.55
Octagonal, molded-on pivot ring, Spinbar®	
White	
12.7 x 8 mm (½ x 5/16"). 371101256	\$4.70
15.9 x 8 mm (5/8 x 5/16"). 371100058	\$4.70
25.4 x 9.5 mm (1 x 3/8"). 371100138	\$6.20
38.1 x 9.5 mm (1½ x 3/8"). 371101128	\$7.40
50.8 x 9.5 mm (2 x 3/8"). 371100238	\$9.90
63.5 x 9.5 mm (2½ x 3/8"). 371102128	\$11.40
76.2 x 12.7 mm (3 x ½"). 371100003	\$20.50
Red, 12.7 x 3.2 mm (½ x 1/8"), no pivot ring. 371090034	\$5.75
Stirring Bar Assortment, octagonal, Pk/6. 371600000	\$36.55
Includes two of each size: 25.4 x 8 mm (1 x 5/16"),	
38.1 x 8 mm (1½ x 5/16"), and 50.8 x 8 mm (2 x 5/16").	
Polygonal, molded-on pivot ring, Spinbar®	
12 x 4.5 mm (7/16 x 3/16"). 371220012	\$5.10
Triangular, vanes w/cross-mounted magnets, Spinvane®	
Fits 1 mL Vials. 371340000	\$9.85
Fits 3/5 mL Vials. 371350000	\$9.85

Continued



Ruler, 1 Decimeter, stainless steel, Germany. 605-060 \$10.20
Ruler is 1 dm long, 1 cm wide, and 1 mm thick. Its volume is 1 cm³ and its weight in grams is the specific gravity. One edge of one side is graduated in centimeters with 1 mm divisions.



This beautiful, transparent ornament is a collage of cytological color photomicrographs creatively captured on a 3" glass disk. Suspended from a gold braid.

Suncatcher, "Cytological Crash Course." CCC01 \$22.50



Tube-Wringer®, Heavy-Duty, **USA. 295401** \$24.80
3" aluminum rollers with steel handle