



Accurately determine the correct size of rubber stopper, size 00-12, or cork, size 0-26, for any vessel using this durable plastic tool.
Reverse side has a circular scale for sizing 4-52 mm OD tubing.

Stopper/Tubing Sizer. 133920000 \$13.20



Surgical Drape, nonfenestrated, 18 x 26", Bx/50. 44025-515..... \$25.50
sterile, poly-lined

Tape

Tape



Autoclave, steam

$\frac{3}{4}$ " x 60 Yd Roll. **44073ST048** \$8.50
1" x 60 Yd Roll. **44073ST036** \$10.70

Kraft paper with rubber resin adhesive.



Label, Write-On™, Bel-Art
1" Core, 15 yd (540") roll

$\frac{1}{2}$ " Wide, each \$5.35

Green. **134623005** Yellow. **134622005**
White. **134620005**

$\frac{3}{4}$ " Wide, Pk/4 \$23.20

Green. **134820075** Yellow. **134810075**
White. **134800075**

3" Core, 40 yd roll

$\frac{1}{2}$ " Wide, each \$9.80


Green. **134633005** White. **134630005**

$\frac{3}{4}$ " Wide, each \$13.00

Green. **134633075** Yellow. **134632075**

$\frac{3}{4}$ " Wide, White, Pk/4. **134830075** \$51.95

$\frac{3}{4}$ " Wide, Mixed, Pk/6. **134630600** \$79.20
blue, green, orange, red, white & yellow



Label Tape Dispenser, multi-roll. 134730000..... \$41.35

Test Strips *Call regarding availability*



Fuses

¼ x 1¼"

Fast-Acting, AGC

0.1A/250V. 312-100	\$2.00
0.2A/250V. 312-200	\$2.00
1A/250V. 312-001	\$2.00
2A/250V. 312-002	\$2.00
3A/250V. 312-003	\$2.00
4A/250V. 312-004	\$2.00

Time-Delay (Slow-Blow), MDL

0.7A/250V. 313-700	\$2.00
1A/250V. 313-001	\$2.00
2.5A/250V	\$2.00

5 x 20 mm

Fast-Acting, GMA


0.5A/250V. 235-500	\$2.00
1A/250V. 235-001	\$2.00
3A/250V. 235-003	\$2.00
5A/125V. 235-005	\$2.00

Time-Delay (Slow-Blow), Glass tube

0.5A/250V. 239-500	\$2.00
1A/250V. 239-001	\$2.00
2A/250V. 239-002	\$2.00
3A/250V. 239-003	\$2.00
5A/250V. 239-005	\$2.00

Time-Delay (Slow-Blow), Ceramic

1A/250V. 215-001	\$2.00
-------------------------------	--------



Cone Point



Cup Point

Set Screws, steel

Cone Point, DIN 914

M4 x 0.7 x 5 mm length, 2 mm hex key size. 77702	\$0.50
M4 x 0.7 x 10 mm length, 2 mm hex key size. 77704	\$0.50
M4 x 0.7 x 12 mm length, 2.5 mm hex key size. 77706	\$0.50
M5 x 0.7 x 10 mm length, 2.5 mm hex key size. 77708	\$0.50

Cup Point, DIN 916

M3 x 0.5 x 3 mm length, 1.5 mm hex key size. 77722	\$0.50
M3 x 0.5 x 5 mm length, 1.5 mm hex key size. 77724	\$0.50

Tools



Cutter, diagonal, 5". PLR-766.00 \$10.95




Demagnetize or magnetize parts and tools.

Demagnetizer/Magnetizer, 1½ x 3" opening. DEM-345.00 \$80.95

Gauges



Vernier Caliper



Compasses

Vernier Caliper, digital, 3-button. 600VCD001 \$28.95
Measures in millimeters (0-150 mm with 0.01 mm precision) and inches (0-6"). Provided with a plastic case.

Compasses, replaceable tips, Germany

3" Long. DIV-320.00	\$11.95
6" Long. DIV-326.00	\$16.95
Tips (replacement), Pair. DIV-320.05	\$1.95



Reproducibility enhanced by rapid stop. Provided with padded, plastic case.

Micrometer, 0-25 mm in 0.01 mm div, Germany. 601A \$39.00



Measure inside and outside diameters of glass and plastic tubing, connecting pieces, rods, test tubes, etc. between 1.5 and 16 mm.

Tubing Sizer, metric diameter gauge, PP. 133900000 \$35.10

Continued