



Thermometer/Clock/Humidity Monitor

Traceable®

Model	Temperature Range	Temperature Resolution/Accuracy	Humidity Range	Humidity Res./Accuracy	Unit Size	Replacement Battery	Price
4040	32 to 122°F/ 0 to 50°C	0.1°±1°C	20 to 90%	1% RH/±4%	4 ¼ x 2 ¼ x ½"	AAA	\$45.45

Warranty: 1 year

Timers, Digital, Traceable®

5004 4-Channel




5000 3-Channel



5008 Three-Line



5020 Mini



5042 Fingertip



Model	Channels	Range	Count Direction	Clock	Memory	Replacement Battery/Battery Cover	Price
5004	4	1 sec - 100 hrs	Down/Up	Yes	Yes	*357/**9588	\$29.65
5000	3	1 sec - 100 hrs	Down/Up	Yes	Yes	AAA	\$49.35
5008	3	1 sec - 100 hrs	Down/Up	Yes	Yes	*357	\$34.35
5020	1	1 sec - 100 min	Down/Up	No	No	*357/**8640	\$23.55
5042	1	1 sec - 100 min	Down/Up	No	Yes	*357/**9588	\$15.55

* Battery (replacement): 357 / \$4.00

** Battery Cover (replacement): 9588 / \$5.65

*** Battery Cover (replacement): 8640 / \$5.65

Warranty: 1 year



Cafeteria Style

Trays

Cafeteria Style, polypropylene, 12 x 16¼", NSF, **USA**

Orange. **105CT1216O** \$4.95

Red. **105CT1216R** \$4.95

Black. **105CT1216B** \$4.95



HDPE



Polypropylene

HDPE (high-density polyethylene), **USA. 1053506** \$38.95
66 x 45.7 x 8.9 cm (26 x 18 x 3½")

Contain spills or transport lab products with this rimmed, durable, large white tray. *Not autoclavable.*

Polypropylene, autoclavable

Small, 30 x 20 x 13 cm (12 x 7¾ x 5⅛"). **162600000** \$22.10

Large, 51 x 39 x 13 cm (20 x 15½ x 5⅛"). **162640000** \$50.05