

- Tubes, PCR**, 0.2 mL thin-wall, polypropylene  
Individual, flat caps, Simport
- Blue, Pk/1000. **T325-2VB** ..... \$44.50
  - Natural, Pk/1000. **T325-2VN** ..... \$44.50
- Strips, Amplitube™, Low Profile, individually attached caps
- 8-Tube, Pk/125 (1000 tubes). **T320-2LPN** ..... \$90.25  
Holds 0.1 mL. Tubes have flat caps & natural color.
- Strips, Bio Plas, *Caps are sold separately; available in natural only.*
- 8-Tube, Pk/125 (1000 tubes)..... \$71.60
  - Blue. **5010-4**      Natural. **5010-1**
  - 8-Tube Cap Strips, domed, Pk/300 (2400 caps). **5030-1** .. \$75.15
  - 12-Tube, Pk/100 (1200 tubes)..... \$80.50
  - Blue. **5020-4**      Natural. **5020-1**
  - 12-Tube Cap Strips, domed, Pk/200 (2400 caps). **5040-1** \$75.15



Transport and store biological specimens in these self-standing tubes. External threads minimize contamination.  
Dimensions of tube with cap:  
17 mm diameter x 84 mm high.

- Tubes, Transport**, polypropylene with cap, 10 mL, Simport
- Plain, Cs/1000. **T550-10AT/T550W** ..... \$121.00
  - Graduated, Cs/1000. **T550-10ATPR/T550W** ..... \$152.00
- Tubes have a writing surface and 0.5 mL graduations from 1-9 mL.




- Vial Decapping Pliers**, vinyl grip, **USA**
- For 13 mm Crimp Seals. **D-13**..... \$72.00
  - For 20 mm Crimp Seals. **D-20**..... \$72.00
  - For 30 mm Crimp Seals. **D-30**..... \$72.00

**Tube Storage**



- Baskets**, test tube, aluminum, Heathrow
- Small, 15.2 x 15.2 x 15.2 cm(H). **HS20341A**..... \$24.90
  - Medium, 25.6 x 15.2 x 15.2 cm(H). **HS20341B**..... \$27.95



- Tube Cooler**, benchtop, 0°C, CoolCube™, Benchmark
- Microtube and PCR Plate. **R1000** ..... \$104.00
- Each cube holds (36) 1.5/2.0 mL tubes, or (96) 0.2 mL tubes, or (12) 0.2 mL strips, or (1) 96-well PCR plate and (10) 1.5/2.0 mL tubes. Maintains 0°C for 4 hours.



HS2860A/H

- Tube Storage Boxes** (cryogenic)
- Cardboard, vented
  - 2" Box Only, 5¼ x 5¼ x 2", Pk/12. **HS2860A**..... \$27.05
  - 2" box w/dividers, white
  - 100-Cell (10 x 10), Pk/12. **HS2860A/E** ..... \$47.80
  - ★ 81-Cell (9 x 9), Pk/12. **HS2860A/H** ..... \$47.80
  - 64-Cell (8 x 8), Pk/12. **HS2860A/F** ..... \$47.10
  - 2" Mini box w/dividers, 3 x 3 x 2"
  - 25-Cell (5 x 5), Pk/12. **MB-2/MD-25** ..... \$50.40
  - Holds 13 mm vials and 4.0-5.0 mL cryovials.
  - 3" Box Only, 5¼ x 5¼ x 3", Pk/12. **HS2860B**..... \$29.30
  - 3" box w/dividers
  - 100-Cell (10 x 10), Pk/12. **HS2860B/E** ..... \$49.20
  - 81-Cell (9 x 9), Pk/12. **HS2860B/H** ..... \$49.20
  - 64-Cell (8 x 8), Pk/12. **HS2860B/F** ..... \$48.50
  - Dividers only, fit standard 2 & 3" boxes
  - 100-Cell (10 x 10), Pk/12. **HS2860E** ..... \$25.90
  - Holds 1.5-2.0 mL tubes and 1.2-2.0 mL cryovials.
  - 81-Cell (9 x 9), Pk/12. **HS2860H**..... \$25.90
  - Holds 13 mm vials and 4.0-5.0 mL cryovials.
  - 64-Cell (8 x 8), Pk/12. **HS2860F** ..... \$24.50
  - Holds 14 mm vials.
  - 25-Cell (5 x 5), 25 mm cell, Pk/12. **D-25** ..... \$21.60
  - 16-Cell (4 x 4), 31 mm cell, Pk/12. **D-16** ..... \$21.60
  - 5" Box Only, 5¼ x 5¼ x 5", Pk/10. **HS2860M** ..... \$44.70
  - 5" Box w/dividers
  - 36-Cell (6 x 6), 15 mL tubes, Pk/10. **HS2860M/T** ..... \$66.30
  - 16-Cell (4 x 4), 50 mL tubes, Pk/10. **HS2860M/P** ..... \$63.10



Polycarbonate box stores vials from -196°C to 121°C.

- Arctic Squares®, polycarbonate, Heathrow
- 81-Well (9 x 9), hinged lid, Pk/4. **HS23381B**..... \$43.35
- Holds vials up to 2 mL; mixed. 13.3 x 13.3 x 5.3 cm

Continued